

SHOWER STOOLS

• 42902 BLUE SHOWER STOOL WITH TRAY

Shower stool made of aluminium alloy frame with plastic seat and tray. 360° swivel seat panel and adjustable height 40-52 cm. Ø seat: 33 cm. Weight: 2.2 kg. Load: 136 kg. Multilingual instructions



• 42904 SHOWER STOOL WITH TRAY
Designed for comfort, safety and independence while washing. Adjustable height 39-54 cm, legs made of anodised aluminium. Anti-slip seat surface with a drainage hole is made of premium quality plastic. The plastic tray enhances stool's stability and it is useful for placing toiletries and essentials. Max load: 170 kg. Multilingual instructions: GB, FR, IT, ES, PT, DE, GR.



MADE IN ITALY

SHOWER COMMODE CHAIR

• 42914 SHOWER COMMODE CHAIR - grey/white

Commode chair with double function: it can be used as a toilet and as a shower chair. Made of aluminium and plastic material with PP frame for better load capacity. Removable armrest and backrest, adjustable height in six levels. It can be easily mounted without any tools. Supplied with backrest and seat cushions, pan with lid and handle and foot pads. Load capacity: 150 kg. Weight: 5.1 kg. Seat size: 46x43xh 42 cm. Packing size: 55x47x24 cm.



Height adjustable in 6 levels from 75 to 87.7 cm
Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, PL.

WC - SHOWER COMMODE

• 27757 VALUE WC - SHOWER COMMODE

Powder coated steel frame. Seat cover and bucket made of PP. Removable bucket. Packing size: 62x35x50 cm



VALUE COMMODE WHEELCHAIR

• 27694 VALUE COMMODE WHEELCHAIR

Designed to provide comfort, safety and mobility, this commode wheelchair is an essential aid for individuals with limited mobility. Powder-coated frame made of steel. Equipped with pan and four castors, two with brakes. Weight: 5.8 kg. Maximum load: 100 kg. Packing size: 51x43x17 cm. Multilingual instructions.



FOLDING BEDSIDE RAIL

• 27704 FOLDING BEDSIDE RAIL

It can be attached to the bed to assist and help those people who have difficulty moving, getting up and lying down.

One-click folding design frame: handle can be easily folded on both left and right side.

Easy to use, single hand operation.

- stepless height adjustment

- 3 levels of handle width adjustment

- comfortable PU handle lowers and raises by 180 degrees

- luminous function for night-time use


- made of iron

Size: 64-84x63.5xh 73.5-91.5 cm. Weight: 3.95 kg

Packaging size: 60x46x36 cm

Max load: 125 kg

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, PL, CZ, HU, HR, SI, RO, BG, GR, Arabic.



OVERBED TABLE - 2 SHELVES

• 44832 OVERBED MULTIFUNCTIONAL TABLE with 2 shelves

68 x 40 cm - white
A versatile overbed table, designed to ensure maximum comfort while dining or reading. It features a powder-coated steel frame, adjustable in height thanks to an ergonomic lever with a gas-lift piston. The white plastic-coated MDF surface is divided into two sections, consisting of a fixed top (20 x 40 cm) and a 33° tilting top (47.7 x 40 cm) with double retaining lips to prevent items from falling. Ultra-slim metal base (3 cm thick). Ultra-smooth swivel castors for silent, 360-degree movement. For hospital and domestic use. Suitable for either side of the bed. Delivered disassembled. Size: 68x40xh 70-103.5 cm. Weight: 11 kg. Max. load: 20 kg. Packing size: 80x55x15 cm.

