

ELECTRONIC CHAIR SCALE

- 25020 DE20 ELECTRONIC CHAIR SCALE - 250 kg - Class III - red - standard
- 25021 DE20 ELECTRONIC CHAIR SCALE - 250 kg - Class III - customized*

High precision electronic chair scale with an innovative design. Functions: weight, tare, hold, BMI, auto zero, programmable auto-off. Painted metal tube structure with adjustable ABS display, 4 360° swivelling castors with 4 safety brakes, padded handle for comfortable transport. Adjustable footrests and reclinable armrests up to 90° to facilitate patients access. Dual LCD display with triple simultaneous reading: weight, height and BMI. Compact packaging reduces shipping costs. Supplied with power adaptor, RS232 interface, level, 4 castors with brakes, adjustable armrests and footrests, multilingual manual: GB, IT, ES. Available on gimaitaly.com: FR, PT, DE, NL, SE, NO, FI. **MDR Certified.**

MADE IN ITALY



Class III



CUSTOMIZABLE OPTIONS for 25021*

- customized seating colours: white, grey, black, red
- unit of measure kg/lb
- customized calibration based on gravity factor
- France: customized calibration by area with carnet métrologique

*available in 20 days

TECHNICAL SPECIFICATIONS

Capacity:	250 kg
Sensitivity:	50/100 g
Display:	Double LCD: LCD 20mm weight with 5 digits; LCD 12mm height and BMI
Size - weight:	590x1040xh 960 mm - 22 kg
Seating size:	550x400 mm
Packag. size:	620x810xh 620 mm
Power supply:	12V - 2A power adaptor (included), 6xAAA alkaline batteries or rechargeable battery (optional)



25026



Class III

STANDARD ACCESSORIES

Power adaptor	Adjustable feet	Foldable handrail (25026-7)	Manual: GB, IT, ES.
RS232 interface	Double ramp	4x360° swivel wheels (25027)	Available on gimaitaly.com:
Level	Carrying handle	2 m connection cable (25025)	FR, PT, DE, NL, SE, NO, FI

TECHNICAL SPECIFICATIONS

	25025	25026 - 25027
Capacity - Sensitivity:	300 kg - 100 g	
Display:	Double LCD: LCD 20mm weight with 5 digits; LCD height and BMI 12mm	
Size:	840x1150xh 70mm	900x1150xh 1100mm
Platform size:	800x1150xh 70mm	
Weighing surface:	740x900 mm	
Handrail capacity	-	90 kg
Weight:	24 kg	30.5 kg (25026), 34 kg (25027)
Power supply:	adapter 12V - 2A (included) 6xAAA batteries (not included) or rechargeable battery pack (optional)	

LABORATORY SCALE

- 31030 WUNDER PROFESSIONAL LABORATORY SCALE WJ600 M

Compact, reliable and portable laboratory scale, designed for rapid and accurate precision weighing. ABS structure with a waterproof keyboard and washable stainless-steel weighing plate. Large 20mm backlit LCD display with 5 digits. Functions: net/gross weight, 10x high-resolution mode, subtractive tare, piece counting, weight %, 7 units (g, ct, oz, ozt, dwt, t, gn), external digital calibration, overload alarm, auto-zero and programmable auto-off. Supplied with power supply, rechargeable battery, RS232, level, adjustable feet, windshield included.

MADE IN ITALY



TECHNICAL SPECIFICATIONS

Capacity - sensitivity:	600 g - d=0,01 g	e=0,1 g
Weighing surface:	Ø 115 mm	
Power supply:	rechargeable battery with 40h run time or 12V power cable	
Size - weight:	200x225xh 80 mm - 1.3 kg	
Packing size:	320x320xh 120 mm	

PLATFORM SCALES - CLASS III

- 25025 WUNDER RW02 DIGITAL PLATFORM SCALE - 300 kg
- 25026 WUNDER RW2.0 DIGITAL PLATFORM SCALE with handrail - 300 kg
- 25027 WUNDER RW2.0 MOVE DIGITAL PLATFORM SCALE with handrail and 4 wheels - 300 kg

Large and easily accessible platform with double ramp, ideal for wheelchair patients. Functions: weight, tare, hold, BMI, auto zero, programmable auto-off. **MDR certified.** Dual ABS display with triple simultaneous readout: weight, height, BMI.

- 25005 HEIGHT ROD 60-210 cm for 25026, 25027 - optional



25025



Easy to transport
25027

The 4 wheels support the platform vertically without the need of additional support

MECHANICAL SCALE

- 25011 WUNDER R150A MECHANICAL COLUMN SCALE 150 kg - Class III

Professional mechanical column scale with sturdy and stable structure in strong aluminium coated with epoxy powder. Platform made of non-toxic, anti-slip rubber with removable mat for better sanitation. Zero adjustment function. Weight system with lever mechanism in cast aluminium and zinc plated steel knife edges to guarantee high performance over time. Large analogue dial Ø 180 mm, easy to read.

MDR certified. Supplied with non-slip feet and user manual in GB, IT, ES. Available on gimaitaly.com: FR, PT, DE, NL, SE, NO, FI.

- 25005 HEIGHT ROD 60-210 cm - optional
- 25012 KIT OF 2 TRANSPORT WHEELS - optional



25011



MADE IN ITALY



Class III

TECHNICAL SPECIFICATIONS

Capacity - sensitivity:	150 kg - 1000 g
Platform size:	240-260x300xh 85 mm
Column size:	Ø 60xh 750 mm
Size - weight:	270x400xh 900 mm - 7 kg