


I.V. STANDS




GIMA code	IV STANDS WITH S/S HOOK	Adjustable height	Weight kg	Base Material	Load kg Per hook Total
27814	Stainless steel IV stand	130-220 cm	5.0	mild steel	2 8
27818	Stainless steel AISI 202 IV stand	135-215 cm	3.5	aluminium	3 12
27815	Stainless steel 201 IV stand	110-200 cm	3.6	stainless steel	3 12

Multilingual manual (only 27814-15): GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic

INFUSION STAND

• **27816 IV STAND FOR INFUSION PUMPS**
IV stand made of high quality 201 stainless steel. Lightweight and easy to move, suitable for supporting single or dual-channel infusion pumps. Slip resistant and powder coated steel base with ABS dust cover to prevent dust accumulation, ensuring more hygienic environment. It is equipped with a multifunctional aluminium tray with integrated handle, which can support infusion pump; four stainless steel hooks and five castors Ø 76 mm two of which with brakes. Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic



TECHNICAL SPECIFICATIONS	
Height adjustment:	1500 - 2500 cm
Weight:	6.4 kg
Tray size:	20x17.5 cm
Tray load capacity:	10 kg

MAYO TABLES



• **45828 MAYO TABLE - H-SHAPE**
Mayo instrument table ideal for delivering and operating room. Rectangular frame made of stainless steel, anti-corrosion. H-shaped base equipped with 4 castors, Ø 50 mm. Adjustable height by side knob. Removable tray.



• **45833 MAYO TABLE - 5 spokes - rotating tray**
Frame and removable tray made of stainless steel. Stainless steel base weighted with ballasts, 5 spokes base with castors Ø 50mm. 360° rotating tray at any height using the fixing bushing.



• **45837 HYDRAULIC MAYO TABLE - 5 spokes - with pump and rotating tray**
Frame and removable tray made of stainless steel. Stainless steel base weighted with ballasts, 5 spokes base with castors Ø 50mm. Possibility to rotate the tray 360° at any height by fixing a handwheel to the rotation joint. Adjustable height by means of a hydraulic foot pump.

TECHNICAL SPECIFICATIONS			
	45828	45833	45837
S/S tray size:	54x42 cm	70x45xh 1.7 cm	70x45xh 1.7 cm
Adjustable height:	81-125 cm	90-115 cm	90-120 cm
Weight:	-	16 kg	10 kg
Load:	20 kg	10 kg	23 kg

NEO PLUS CLINICAL TROLLEY

• **45715 NEO PLUS CLINICAL TROLLEY**
Clinical trolley made of aluminium, steel and ABS structure. Upper part is equipped with stainless steel 304 three-sided rail, which prevents items from falling. It is equipped with four castors Ø 100 mm (two with brake), one drawer with silent guide rails and ABS inserts, which can be freely organized, and card holder. Size: 62.5x47.5x92 cm. Max load: table top 80 kg, middle tray 15 kg.

STANDARD ACCESSORIES	
2 ABS waste bins 7 l with lid	2 shelves
1 drawer with dividers	1 fixed tray with inner hollow depression

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, PL, RO, CZ, HU, SE, FI, LT, HR, GR, BG.

