

ROYAL ADJUSTABLE EXAMINATION COUCH - 3 SECTIONS

• **44541 ROYAL ADJUSTABLE EXAMINATION COUCH - blue**
Couch frame is made of painted powder coated steel with strong corrosion resistance: this structure allows to maintain more stability, making patients feel safe.
Couch surface is made of synthertic leather and hand-made cotton mattress, which make it soft and breathable. Easy to clean and disinfect.
It is equipped with manual air levers with two gas springs to adjust the back and leg plate angle, one hand controller to adjust height.
The brake system is composed by four brake pedals connected to four castors Ø 76 mm.
360° spinning castors made of strong, durable and non-slip polyurethane PU.
Size: 201x69xh 48-91 cm. Weight: 68 kg.
Safe working load: 200 kg.
Voltage: AC 220 V, 50 Hz.
Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic



Matching cushion included

Adjustable arm support

Leg section up/down 53°/24° ±5°

Back section up/down 80°/35° ±5°

Height adjustment: 48-91 cm

44541

LOAD 200 kg

2 SECTION EXAMINATION COUCH - STAINLESS STEEL

• **44460 VIGOR STAINLESS STEEL EXAMINATION COUCH**
Examination couch with structure made of stainless steel.
It is equipped with a mechanical lever to adjust backrest up to 65° and 4 anti-slip feet covers. Covering made of flame-retardant leather that can resist to yellowing, mould and to 100,000 wearing times. Padding is made of flame-retardant sponge, 6 cm thickness, high-density foam (density 30 kg/m³).
Size: 190x60xh 68 cm.
Delivered disassembled.
Supplied with couch roll holder (max paper roll height: 60 cm) and multilingual manual (GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic).



Adjustable backrest (0° to 65° ±10%)

Thickness 6 cm

44460

Leg Ø 50 mm

Anti-slip feet cover

Couch roll holder (paper roll not included)

LOAD: 175 kg
LENGTH: 190 cm

2 SECTION EXAMINATION COUCHES WITH DRAWERS

• **44465 GREY PAINTED EXAMINATION COUCH WITH 3 DRAWERS**

• **44467 STAINLESS STEEL EXAMINATION COUCH WITH 3 DRAWERS**

Structure made of enamelled plated steel (44465) or stainless steel (44467) tube Ø 30x30 mm.

Adjustable backrest by rack.

Padding and covering made of washable, flame retardant, plastic material.

Both equipped with three drawers to be placed under the fixed part.

Size: 177x60xh 80 cm

Drawer internal size: 38x45x h 10 cm

Single drawer load capacity: 10 kg

Packing size 1 box 74x21x202 cm + 1 box 77x51x55 cm

Class 1 IM.

Multilingual manual and instructions: GB, FR, IT, ES PT, DE, NL, DK, SE, PL, CZ, HU, SI, HR, RO, BG, GR, Arabic.

Delivered disassembled.

• **27418 COUCH ROLL HOLDER - optional**



MADE IN ITALY

LOAD 150 KG

COUCH PAPER ROLLS - PEDIATRIC

• **27423 MICROEMBOSSED GLUED 2 PLIES PEDIATRIC COUCH**

ROLL - 46m x 50cm - box of 9

Decorated medical rolls for pediatric use.

Produced in pure cellulose, dermatologically tested.

Hypoallergenic inks.

Each box contains 9 rolls with 5 assorted themes.

COUCH PAPER ROLLS	
MATERIAL	Microembossed glued
Colour	white+decoration
Ply	2 ply
Weight	2x17.5 g/m²
Length	46 m
Height	50 cm
N. of sheet	135
Perforation	34 cm
Rolls per box	9

27423

