

NEO PLUS ANESTHESIA TROLLEY

• 45718 NEO PLUS ANESTHESIA TROLLEY

Anesthesia trolley made of ABS and antibacterial materials externally, while the internal structure is composed of aluminium. It is equipped with five big and deep drawers with dividers and four anti-static castors Ø 10 cm (two with brake).

Size: 75x47.5x95 cm.

90x52x138 cm with accessories.

Weight: 44.5 kg with accessories.


Drawers inner size:

2 drawers 43x38x7 cm

2 drawers 43x38x11 cm

1 drawer 43x38x22 cm

Medical Device Class I



Pull-out worktop

Storage basket

Hand sanitizer holder

Removable container

45718

Foot pedal controls lids opening

STANDARD ACCESSORIES

Central lock 2 ABS waste bins with single lid
Pull-out worktop 10 anesthesia boxes
Dividers Hand sanitizer holder, basket
Utility container Sharp container 2l



NEO ANESTHESIA TROLLEY

• 45726 NEO ANESTHESIA TROLLEY

Anaesthesia trolley featuring a steel and aluminium internal structure and an antibacterial ABS outer shell. Easy to clean. It is equipped with 5 drawers with dividers and 4 anti-winding hard-shell castors. Size: 75x47.5x95 cm.

Weight: 35 kg.

Drawers inner size:

2 drawers 43x33.5x6.8 cm

2 drawers 43x33.5x11 cm

1 drawer 43x33.5x22 cm



Dual rows of adjustable storage boxes

Table top covered by transparent glass

Pull-out worktop

Handle on both sides



Oxygen bottle holder

Removable container

Two ABS waste bins with lid

STANDARD ACCESSORIES

Central locking Storage basket
Pull-out worktop Hand sanitizer holder
10 anesthesia boxes Utility container
Sharp container Oxygen bottle holder
2 ABS waste bins with lid 10 disposable seals
Adjustable dividers for drawers



NEO PLUS EMERGENCY TROLLEY

• 45720 NEO PLUS EMERGENCY TROLLEY

Emergency trolley made of ABS and antibacterial materials externally, while the internal structure is composed of steel and aluminium. It is equipped with five big and deep drawers with dividers and 4 castors Ø 10 cm (two with brake).

Size: 75x47.5x95 cm.

Weight: 41 kg with accessories.

Drawers inner size:

2 drawers 43x38x7 cm

2 drawers 43x38x11 cm


1 drawer 43x38x22 cm.

Medical Device Class I

• 45721 NUMBERED

SEALS - box of 100

DOUBLE LOCK
SEAL + KEY



Pull-out worktop

Movable IV pole and defibrillator plate

Table top covered by transparent glass

Handle on both sides

45720

STANDARD ACCESSORIES

IV pole Removable container
Defibrillator plate Utility container
Pull-out worktop CPR board
Dividers Oxygen bottle holder
Power outlet Hand sanitizer holder
Central lock Sharp container 2l
10 numbered seals 2 ABS waste bins with lid



NEO PLUS TREATMENT TROLLEY

• 45716 NEO PLUS TREATMENT TROLLEY

Treatment trolley made of ABS and antibacterial materials externally, while the internal structure is composed of aluminium.

It is equipped with four castors Ø 10 cm (two with brake) and two drawers with ABS inserts, which can be freely organized.

Size: 75x47.5x93 cm. Weight: 23 kg.

Multilingual manual:

GB, FR, IT, ES, PT, DE, NL, PL, RO, CZ, HU, SE, FI, LT, HR, GR, BG.




Table top covered by transparent glass

Handle on both sides

Hand sanitizer holder and storage basket

Inner size: 43x38x11 cm

Internal size: 50x42x6 cm

45716

STANDARD ACCESSORIES

Sharp containers 2 l
2 ABS waste bins with lid
Dividers
Hand sanitizer holder
Storage basket

