

Wi-Fi EARWAX REMOVAL TOOL WITH CAMERA

• 31405 Wi-Fi EARWAX REMOVAL TOOL WITH CAMERA

The earwax removal kit is equipped with 6 LED lights and a 1080 P HD camera to help to see and clean earwax effectively. The camera can be connected to devices via App, to store images and videos.

It comes with 3 kinds of ear scoops with different curvatures, to choose which one fits best:

- metal ear scoop
- silicone ear scoops, to prevent scratching the ear canal
- plastic ear scoops

Spoons and lens can be washed with water or wiped with alcohol.

Suitable for children > 3 years old and adults. Medical device, Class I. Supplied with 1 metal scoop, 4 silicone scoops, 4 plastic scoops, 8 tools for hand-cleaning of ears, charging cable and user manual in GB, FR, IT, ES, PT, DE, NL, PL, SE, HU, RO, GR.

• 31406 ACCESSORY KIT (4 plastic scoops + 4 silicone scoops) for 31405 - box of 10 kits

TECHNICAL SPECIFICATIONS			
Image sensor:	CMOS	Gyroscope:	3-axis
Pixels:	0.3 megapixels	Lens Ø:	3.5 mm
Frame rate:	20FPS	Charging time:	60 min
Focal length:	15-20 mm	Working time:	90 min
DOP (Depth of field):	10-50 mm		
Power supply:	built-in 200mAh rechargeable Li-on battery		
Input power:	DC 5V/0.5 A		
Working frequency:	2.4 Ghz		
Size - weight:	13.55x1.50 cm - 17.8 g		
Protection:	IP67		

CE

- 1 - Ring type for oily ears
- 2 - Rotating type to remove debris of ear canal
- 3 - Small type for children or small ear canal
- 4 - Standard type for adult

31405

Smartphone not included


SCHIMMELBUSCH SYRINGES FOR EAR WAX REMOVAL

BRASS SCHIMMELBUSCH SYRINGES

Reliable and accurate syringes.

Schimmelbusch syringes come complete with 3 cannulas and ear bowl. Individually packed in cardboard box.

Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.



Patented 3-stream action.
Compatible with 25804, 25809
standard Luer Lock fitting (conical shape)

GIMA code	BRASS SYRINGES COMPLETE WITH 3 CANNULAS AND EAR BOWL
25835	Schimmelbusch 50 cc
25836	Schimmelbusch 100 cc
25837	Schimmelbusch 150 cc
25838	Schimmelbusch 200 cc
25839	Spare cannula for Schimmelbusch (conical shape)
25840	Spare cannula for Schimmelbusch (olive shape) small
25841	Spare cannula for Schimmelbusch (olive shape) large
25861	Ear irrigation tip - single use - box of 40



GLASS-METAL EAR SYRINGES FOR EAR WAX REMOVAL

Reusable, glass-metal medical syringe typically used for ear irrigation or fluid aspiration. It features a durable design, is sterilizable and comes with interchangeable cannulas.

Used in ENT settings, it's valued for precision and longevity.

Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.

GIMA code	GLASS-METAL SYRINGES
25843	Janet record 50 cc
25844	Janet record 100 cc
25846	Janet record 200 cc
25839	Spare cannula for Schimmelbusch (conical shape)
25840	Spare cannula for Schimmelbusch (olive shape) small
25841	Spare cannula for Schimmelbusch (olive shape) large
25861	Ear irrigation tip - single use - box of 40

