


HEINE GAMMA® G5, G7, GP SPHYGMOMANOMETERS



The HEINE GAMMA series is designed for life-long accuracy, every model is shock-resistant, latex-free and features a pressure-proof manometer protected by a micro-filter.

- high quality indicator, precisely calibrated, durable and shock resistant
- large Ø 56 mm black-on-white scale for easy, accurate readings
- housing is moulded of shockproof plastic with shock absorbing, anti-slip coating
- large insufflation bulb (50% larger) for fast cuff inflation
- adjustable spoon for both right and left-handed operation
- high-quality, washable single-tube cuff with velcro fastener

HEINE GAMMA® G5

Exclusive air-release valve provides fine and precise air release rate adjustment, quick deflation, and extreme durability.

HEINE GAMMA® G7

Exclusive air-release valve as G5. Upper housing made of aluminium alloy with a shock absorbing ring protecting indicator and glass. Ergonomic design for comfort and ease of use.

GAMMA® GP

Standard sphygmomanometer with push-button valve for fine and precise air-release rate adjustment, ultra-fast air-release, practically wear-free.

GIMA code	HEINE code	HEINE GAMMA SPHYGMOMANOMETERS	Child cuff (13-20 cm)*	Adult cuff-small (20-29 cm)*	Adult cuff (29-41 cm)*	Zipper case
32737	M-000.09.230	HEINE GAMMA G5 - adult			•	•
32741	M-000.09.561	HEINE GAMMA G5 - adult			•	•
32758	M-000.09.231	HEINE GAMMA G5 - pediatric	•		•	•
32679	M-000.09.554	HEINE GAMMA G7 - family kit	•	•	•	•
32739	M-000.09.232	HEINE GAMMA G7 - adult			•	•
32759	M-000.09.242	HEINE GAMMA GP - ad. - push button valve			•	•

* Circumference




WITTY ANEROID SPHYGMOMANOMETERS

Aneroid sphygmomanometer with nylon cuff and ABS plastic case. 2 tube adult size PVC bladder and PVC bulb. Dial diameter 51.9 mm. Dual use air-release valve with spring. Supplied with a black drawstring bag. Packing size: 108x65x168 mm.

32764



WITTY ANEROID SPHYGMOMANOMETERS		
GIMA code	Colour	Languages
32762	black	GB, ES, PT, NL, SE, LT, PL, CZ, HU, HR, RO, GR
32763	black	FR, IT, DE, DK, NO, FI, LV, EE, SK, SI, BG, Arabic
32764	light blue	GB, FR, IT, ES, PT, DE, SE, PL, HU, HR, RO, GR



MADRID SPHYGMOMANOMETER

• 50000 MADRID ANEROID SPHYGMOMANOMETER

Professional, high quality sphygmomanometer with 0 to 300 mmHg scale. Two handling modes for right and left-hand users. Large dial Ø 59.5 mm with PVC shockproof ring. Zinc-alloy handle with nickel plated. Copper male and female connector. Supplied with adult nylon cuff (22-32 cm) with D-ring, 1-tube PVC bladder, ribbed PVC bulb, black nylon bag and multilingual box and manual GB, FR, IT, ES, PT, DE, NL, SE, PL, CZ, HU, HR, RO, BG.



Bright dial



Male and Female connector



Left and Right hand use