

## GIMA WRAP BLOOD PRESSURE MONITOR WITH BLUETOOTH

### • 49870 GIMA WRAP BLOOD PRESSURE MONITOR with Bluetooth

Portable and lightweight blood pressure monitor with Bluetooth connection, suitable for professional and domestic use.

- large and easy-to-read 3.2" LED screen with backlit display
- 3 colour backlit indicator detection
- 2 users, 120 memories each + guest mode, ideal for family use
- IHB detection
- "keep still" indicator for movement error detection
- "cuff worn" detection for correct cuff worn
- built-in rechargeable lithium battery

Supplied with adult cuff (22-42 cm), lithium battery and type-C USB cable, user manual in GB, FR, IT, ES, PT, DE, NL, SE, PL, HU, RO, GR.



#### TECHNICAL SPECIFICATIONS

|                           |                                                                |                |                                            |
|---------------------------|----------------------------------------------------------------|----------------|--------------------------------------------|
| Measuring method:         | Oscillometric                                                  | Display:       | LED screen, 75x35 mm                       |
| Pneumatic pressure range: | 0-295 mmHg                                                     | Cuff size:     | 22-42 cm                                   |
| Blood pressure rate:      | sys 57~255 mmHg; $\pm 3$ mmHg<br>dia 25~195 mmHg; $\pm 3$ mmHg | Size - weight: | 59x123xh 28 mm - 225 g                     |
| Pulse rate:               | 40~199 bpm; $\pm 5\%$                                          | Power supply:  | rechargeable lithium battery 3.7V, 450 mAh |
| Operating mode:           | continuous operation                                           | Charging time: | 2 hours                                    |
|                           |                                                                | Protection:    | internally powered supply, type BF, IP22   |

## BP2A BLOOD PRESSURE MONITOR WITH BLUETOOTH

### • 49865 BP2A BLUETOOTH BLOOD PRESSURE MONITOR

Intelligent automatic upper arm blood pressure monitor designed for daily use. The all-in-one design makes it easy to carry.

**Clinically validated:** combination of Intelligent pressurization technology and the brand-new Lycra high elastic cuff design ensure the most accurate measurement results

**Easy to read and to use:** one click measurement and OLED screen shows real-time monitoring

**Reliability:** it detects body movement and cuff wearing to avoid invalid results

Bluetooth sync allows to monitor, share and save results on ViHealth App, which supports unlimited storage and multi-

user management. Supplied with integrated cuff (22-42 cm), built-in rechargeable battery, charging cable, quick guide and multilingual manual in: GB, FR, IT, ES, DE.



49865

Phone not included

#### TECHNICAL SPECIFICATIONS

|                         |                                                     |
|-------------------------|-----------------------------------------------------|
| Measuring method:       | oscillometric method                                |
| Pressure measur. range: | 0-300 mmHg; $\pm 3$ mmHg or 2% whichever is greater |
| Pulse rate range:       | 40 to 200 bpm; $\pm 2$ bpm                          |
| Storage:                | 50 records                                          |
| Charge input:           | USB Type-C, DC 5V                                   |
| Power supply:           | Li-ion rechargeable battery                         |
| Battery life:           | 500 measurements                                    |
| Charging time:          | 2 hours                                             |
| Bluetooth:              | Bluetooth 4.0 BLE                                   |
| Size - Weight:          | 135x45xh 20 mm - 240 g                              |

## BP3 BLOOD PRESSURE MONITOR + ECG

Download Vihealth App on:



Clear and easy to read  
4.3" colour display

35191



### • 35191 BP3 BLOOD PRESSURE MONITOR + ECG

The BP3 is a digital monitor designed for accurate measurement of blood pressure (BP) and pulse rate (PR) in adults. The 2-in-1 design enables simultaneous monitoring of ECG and blood pressure, providing a comprehensive view of cardiovascular health and helping to prevent complications through multi-dimensional analysis. The ECG recording can be performed using two methods: by placing both hands on the electrodes or by using the right hand and left leg. Multiple user profiles, allow to manage personal health data.

With advanced AI-ECG technology, available on ViHealth App, the device can detect up to 17 types of results with 95.2% accuracy, enabling intelligent screening and early warning of heart abnormalities, including irregular heart rhythms. BP3 supports Bluetooth and WiFi, allowing for real-time synchronized data transmission to compatible health platforms or mobile devices.

Supplied with adult cuff (22-42 cm), air hose, USBC charging cable and manual in GB, FR, IT, ES, PT, DE, PL, RO.

• 35194 CUFF 22-42 cm for 35191 - spare

#### TECHNICAL SPECIFICATIONS

| Blood pressure measurements |                                                    | Device                 |                              |
|-----------------------------|----------------------------------------------------|------------------------|------------------------------|
| Measuring method:           | oscillometric method                               | Battery:               | rechargeable lithium battery |
| Pressure range:             | 0-300 mmHg; $\pm 3$ mmHg                           | Charging time:         | 2 hours                      |
| Pulse rate range:           | 40-200 bpm; $\pm 5\%$                              | Battery life:          | $\geq 200$ times             |
| ECG recording               |                                                    | Input:                 | DC 5V                        |
| Lead type:                  | integrated electrodes                              | Cuff size:             | 22-42 cm                     |
| Lead set:                   | Lead I, Lead II                                    | Wireless connectivity: | Bluetooth 5.2 BLE, WiFi 2.4G |
| Heart rate range:           | 30-250 bpm; $\pm 2$ bpm or 2% whichever is greater | Size - weight:         | 120x121.5xh 49 mm - < 500 g  |
| ECG length:                 | 30s - 5 min                                        | Protection:            | Type BF, IP22                |
| Input impedance:            | $\geq 10$ M $\Omega$ , 10 Hz                       |                        |                              |
| CMRR:                       | $\geq 60$ dB                                       |                        |                              |

