

## X-CHECK PRO AUTOMATIC B.P. MONITOR - BLUETOOTH



### • 49872 X-CHECK PRO AUTOMATIC B.P. MONITOR with Bluetooth

Automatic blood pressure monitor which combines smart technology with medical-grade accuracy to make home monitoring easier and more reliable.

It features a large 4.2" LCD display with a clear 4-colour backlight system. The upgraded algorithm and intelligent inflation technology ensure precise readings, while the irregular heartbeat (IHB) detection helps identify abnormal heart rhythms during each measurement.

A cuff wearing detection system confirms when the cuff is properly positioned with an "OK" indicator, and movement detection ensures that measurements are taken under the right conditions.

- one-button operation
- 2x120 memories with date and time
- average reading of last 3 measurements

It pairs with the free Tuya Smart App, to track and manage health data.

Supplied with adult soft nylon cuff 22-42 cm, pouch and multilingual manual in GB, FR, IT, ES, PT, DE, NL, SE, PL, HU, RO, GR.

### • 49877 ADULT CUFF 22-42 cm - spare

#### TECHNICAL SPECIFICATIONS

Measuring method: oscillometric  
Measuring range: pressure: SYS 57-255 mmHg  
DIA 25-195 mmHg  
pulse: 40-199 bpm  
Display: digital LCD display, 4-colour backlight

Display size: 80x68.9 mm  
Memory: 2x120 memories  
Power supply: 4xAAA batteries or type-C USB cable (both not included)  
Power off: automatic switch off after 60 seconds  
Size - weight: 108x139x62 mm - 316 g (without batteries)

## GIMA AFIB BLOOD PRESSURE MONITOR

Dual measurement modes:  
- Afib mode  
- Normal mode



### • 32876 GIMA AFIB BLOOD PRESSURE MONITOR

This blood pressure monitor integrates patented technology for precise signal analysis, enabling early and accurate atrial fibrillation (AFib) detection. It automatically identifies irregular heartbeats during each measurement.

The device offers one-touch operation with comfortable automatic inflation and deflation, ensuring a smooth and user-friendly measuring experience.

The device features a 2-user memory function with 199 sets of readings per user and guest mode, allowing multiple individuals to store and review their results.

- correct cuff position detection
- large 4.5-inch easy-to-read LED display
- movement detection

Medical Device, Class II. Supplied with adult cuff (22-42 cm), pouch, 4xAAA batteries and multilingual manual GB, FR, IT, ES, PT, DE, SE, PL, HU, RO, HR, GR.

#### TECHNICAL SPECIFICATIONS

Measuring range: Pressure: SYS: 57-255 mmHg  
DIA: 25-195 mmHg  
Pulse: 40-199 mmHg

Cuff size: 22-42 cm  
Power supply: 4xAAA batteries (included)  
Size - weight: 130x112x59 mm - 280 g

## X-CHECK B.P. MONITORS

### • 49874 X-CHECK BLOOD PRESSURE MONITOR GB, FR, ES, DE, NL, SE, LT, PL, CZ, HR, RO, GR

### • 49875 X-CHECK BLOOD PRESSURE MONITOR IT, PT, DK, NO, FI, LV, EE, SK, HU, SI, BG, Arabic

Accurate blood pressure monitor for home use.

- 2 users, 99 memories for each user
- IHB detection and WHO indicator
- average reading of last 3 measurements
- time and date display
- low battery indicator
- it can work with type-C USB cable (not included)

Supplied with cuff (22-42 cm), dustproof plastic bag, 4xAAA batteries and user manual (see above languages)

### • 49877 ADULT CUFF 22-42 cm - spare

#### TECHNICAL SPECIFICATIONS

Measuring method: oscillometric  
Blood pressure rate: sys 57~255 mmHg;  $\pm 3$  mmHg  
dia 25~195 mmHg;  $\pm 3$  mmHg  
Pulse rate: 40~199 bpm;  $\pm 5\%$   
Operating mode: continuous operation  
Display: LCD screen, 58.5x44 mm  
Cuff size: 22-42 cm  
Size-weight: 98x118xh 62.5 mm - 220 g without batteries  
Power supply: 4xAAA batteries DC 6V or DC 5V type C cable  
Automatic switch off: 1 minute  
Protection: internally powered supply, type BF, IP21

