

GIMA EVACUATION STRETCHER

• 56057 GIMA EVACUATION STRETCHER - foldable - black

Gima evacuation stretcher has a special design which allows to transport patients up and down stairs. The stretcher is foldable and made of high strength aluminium alloy and waterproof PVC. Equipped with 2 small rear castors with brakes, 2 front wheels and 2 safety belts made of polypropylene fiber with plastic buckle.

TECHNICAL SPECIFICATIONS

Size:	92x52x105 cm
Folding size:	114x52x26 cm
Packing size:	115x54x27 cm
Max. load:	160 kg
Weight:	12 kg

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.



56057

GIMA STAIR STRETCHERS

• 56050 GIMA STAIR STRETCHER - foldable - blue

• 56053 GIMA STAIR STRETCHER - foldable - yellow/black

Foldable stair stretchers made of high strength aluminium alloy and waterproof PVC. Suitable for patient transport up and down stairs. Both equipped with 4 wheels, 2 with brakes and 2 (56053) or 3 (56050) safety belts made of polypropylene fiber with plastic buckle. Easy to clean and to sterilize. 2 front wheels with brake (56050), 2 rear wheels with brake (56053).

TECHNICAL SPECIFICATIONS

	56050	56053
Size - Weight:	84x52x91 cm - 9.5 kg	84x50x100 cm - 10.6 kg
Folded size:	97x52x17 cm	100x55x27 cm
Max. load:	160 kg	160 kg



56050

56053

GIMA BASKET STRETCHER

• 56100 GIMA BASKET STRETCHER

This basket stretcher is designed for a wide range of emergency situations. It is sturdy and flexible, making it perfect for specialized rescue work. The reliable fittings allow first-aid teams to operate quickly and safely. This stretcher is ideal for outdoor, underground and water rescue. The sling is originally designed for rope-type rescue operations. The unit comes with an adjustable foot-securing mechanism, 4 safety belts (plastic buckle), one set of lifting bridles, mattress and 4 connecting pieces. All materials are fire-proof and highly resistant to wear and corrosion. They are safe for the environment, as they do not release any toxic or polluting substances.

Size: 218x62x16 cm. Weight: 15 kg. Max load: ≤ 280 kg.

Packing size: 220x67x20 cm.

• 56101 LIFTING BRIDLES - spare



Black mattress
8 mm thickness

56100

56101

VACUUM MATTRESS STRETCHERS

• 34102 VACUUM MATTRESS STRETCHER

• 34075 VACUUM MATTRESS STRETCHER PLUS

Vacuum mattresses used to transport injured patients, to prevent compression and mitigate the risk of additional injuries.

High quality and durable TPU material, filled with harmless EPS particulate material and shaped to conform to the human body. They are compatible with X-rays and CT scan.

Healthcare professionals can use the air pump to generate vacuum or negative pressure, facilitating easy and rapid adjustment of softness and hardness based on severity of injury. Lightweight, foldable after deflation and easy to carry.

Both equipped with four handles on both sides, supplied with air pump and carrying bag.

Max. load capacity:

160 kg (34102), 200 kg (34075).

Unfolded size:

204x104x7 cm (34102)

204x100x5 cm (34075)

Weight: 6 kg



34102

34075

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic