

KONIX® BURN GEL & DRESSINGS - STERILE

Konix® Sterile First Aid Gel and Burn Dressing is designed for first aid in all kinds of burn injuries (heat burns, sun burns, electrical burns, scalds, cuts, abrasions, etc.). It stops burning feeling and prevents further burns from progressing on the skin. It quickly cools the burn to prevent further tissue damage and aid the healing process, and offers the patient a rapid and

continuous reduction in sensitivity. It contains Melaleuca (tea tree) oil and Aloe Vera. Class IIb. Suitable for infants, children, adults.

GIMA code	BURN GEL & DRESSINGS - STERILE	Minimum order
57002	Burn dressing 10x10 cm	box of 50
57003	Burn dressing 20x20 cm	box of 10
57004	Burn dressing 30x30 cm	box of 18
57006	Burn dressing 40x60 cm	box of 12
57010	Burn gel - bottle 118 ml	box of 20
57012	Burn gel - sachets 3.5 g	box of 1,000
57013	Burn gel - sachets 3.5 g	box of 20



57012-13 57010 57002 57003 57006

FIRST AID KIT FOR EYES, BURNS, DRESSING

• 34886 FIRST AID KIT FOR EYES, BURNS AND DRESSING

3 in 1 kit that can be used in 3 different emergency situations: eye wash, burns and dressing. It includes 9 products for quick first aid.

Size: 340x80xh 380 mm.
Weight: 1.15 kg (with products), 306 g (empty).

STANDARD ACCESSORIES

1 eye wash solution 500 ml + cap compliant with DIN EN 15154-4	1 scissors for cutting bandages 14 cm
1 PIC 3 case: 8 wipes (3 disinfectants, 2 ammonia, 3 liquid soap)	1 plaster roll 5 m x 2.5 cm
2 bandage 80x100 mm	1 Dry Skin waterproof plasters 40 assorted pieces
	1 blue Soft Next cohesive bandage 3x450 cm
	1 burn gel bottle 50 ml

34886



PROFESSIONAL BURN FIRST AID KIT

• 34885 PROFESSIONAL BURN FIRST AID KIT

First aid kit ideal for all companies, vehicles and workplaces where burns and minor injuries can occur. It contains hydrogel burn sachets that are specifically formulated to soothe and cool minor burns, for effective and prompt response. Case made of shockproof polypropylene with wall support, carrying handle and closure with two rotating clips.

STANDARD ACCESSORIES

4 sterile anti-burn gauzes soaked in BURNSHIELD Hydrogel 10x10 cm	
4 sterile gauze compress 10x10 cm 8 layers	
10 anti-burn gel sachets 3.5 g	2 pairs of sterile gloves
1 DIN lister scissors 14.5 cm	2 sterile TNT drapes for burns 40x60 cm
1 sterile saline solution 100 ml	3 bandage with non-stick pad 100 x 120 mm
1 disposable sterile tweezer 10 cm	1 anti-burn hydrogel bottle 50 ml
1 plaster roll 5m x 2.5 cm	1 skin-coloured self-adhering plaster 6x100 cm

Size: 260x205xh 90 mm.
Weight: 1.2 kg.
MDR certified.



34885

SOFT FIRST AID KIT

• 34140 SOFT FIRST AID KIT

Red nylon bag in shrink wrap with zip that contains first aid kit for minor injuries, useful complement for private cars, leisure activities, DIY activities... MDR certified.

STANDARD ACCESSORIES

1 pair of sterile gloves	1 roll of adhesive plaster 5 m x 1.25 cm
1 elastic bandage 6x4 cm	1 pack of 10 tissues
2 x 5ml disinfectant liquid soap sachets	1 pack of 10 assorted plasters
3 disinfectant wipes	2 sterile gauze sachets 18x40 cm
1 tweezers	1 haemostatic tourniquet
1 pair of scissors	

Size: 170x110x60 mm.
Weight: 200 g



34140

GIMA DISPOSABLE PRESSURE INFUSION BAGS FOR I.V. SET

- 32781 DISPOSABLE PRESSURE INFUSION BAG - 500 ml - box of 50
- 32782 DISPOSABLE PRESSURE INFUSION BAG - 1,000 ml - box of 50
- 32783 DISPOSABLE PRESSURE INFUSION BAG - 3,000 ml - box of 50

Safe, high-quality and reliable devices for intra-arterial infusion treatment, including A-line pressure monitoring. Designed for single-patient use, these bags help prevent and control cross-contamination and hospital-acquired infections.

- pressure relief valve prevents over inflation at 300 mmHg
- large, oval shaped bulb for quick and easy inflation of bladder
- single-handed inflation and deflation design makes it easy to use
- colour-coded gauge ensures accurate pressure monitoring
- three-way stopcock provides precise control of pressure
- dedicated fluid bag hook enables efficient loading and unloading without having to remove the bag from the IV pole.

Suitable for emergency blood transfusions, infusions and various invasive arterial pressure monitoring procedures in emergency department and operating room.

Each unit is leakage tested, latex-free and non-sterile. Individually packed in PE transparent plastic bags.

Multilingual box and manual:
GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic



32781

32782

32783

TECHNICAL SPECIFICATIONS

Pressure range:	0-300 mmHg
Max. pressure:	300 mmHg
Tubing size:	inner Ø: 4 mm outer Ø: 7 mm length: 200 mm