

GIMA STRETCHERS AND SPINAL BOARDS



• **56005 GIMA STRETCHER 2 - foldable in 2 - black**
Foldable emergency stretcher made of high strength aluminium alloy. Lightweight, small-sized, easy to carry and sterilize. Supplied with 2 safety belts made of polypropylene fibre with plastic buckle. Multilingual manual*



• **56010 GIMA STACKABLE STRETCHER - orange**
Stackable stretcher made of high strength aluminium alloy and leather. Durable, waterproof and small sized so that it can be piled up to save space. Supplied with 2 safety belts made of polypropylene fibre with plastic buckle. Multilingual manual*

• **56015 GIMA CAMP BED - foldable - green**
Foldable and lightweight camp bed made of high strength aluminium alloy and 600D polyester. It is mainly used for army and battle field. Supplied with carrying bag. Multilingual manual*



• **56020 GIMA 2 IN 1 FOLDABLE STRETCHER CHAIR - orange**

Emergency stretcher convertible into 3 positions: chair, stretcher and folded. Made of high quality aluminium alloy and vinyl coated nylon. Suitable for transporting patients in narrow spaces, it has two safety belts, two handles and two wheels to facilitate mobility on the ground. Lightweight, waterproof and easy to carry. Multilingual manual*



• **56040 GIMA SPINAL BOARD WITH SAFETY PINS - orange**
Suitable for patient with head, hip, pelvis and spine injuries. It can float on water. Made of HDPE material. Supplied with safety belts. It can be used together with spider straps and head immobilizer. Multilingual manual*

• **56044 GIMA SPINAL BOARD PLUS WITH SAFETY PINS - orange**
Suitable for patient with head, hip, pelvis and spine injuries. It is made of HDPE material with improved design and two bamboo battens inside to strengthen the stretcher. It can be used with head immobilizer. Multilingual manual*

• **56047 GIMA PEDIATRIC SPINAL BOARD - orange**
Suitable for patient with head, hip, pelvis and spine injuries. It can float on water. Made of HDPE material. Supplied with safety belts. It can be used together with spider straps and head immobilizer. Multilingual manual*

TECHNICAL SPECIFICATIONS	56005	56010	56015	56020	56040	56044	56047
Open size cm:	185x50x3	188x51x19	196x75x43	196x50x18	184x45x6	185x46x5	136x35x5
Folded size cm:	92.5x50x7	-	95x20x11	106x50x15	-	-	-
Weight:	6.5 kg	3.5 kg	6.4 kg	9 kg	8 kg	7 kg	4,5 kg
Load:	159 kg	159 kg	159 kg	159 kg	181 kg	227 kg	91 kg

* Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO**, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic. ** Not for 56020

CPR BOARD



• **56528 CPR BOARD - blue**

CPR board provides a raised surface to perform chest compression. The board is rugged, flexible and lightweight. The head shape ensures the patient airway is open. Made of HDPE material.

Size: 60x41.5x8 cm. Weight: 1.55 kg.

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.

NEWS - 70

CARRYING SHEET

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.



• **55995 GIMA FOLDABLE CARRYING SHEET 200x70 cm - yellow**
Foldable carrying sheet ideal for small transfers of patients. Equipped with 10 handles. Made of 600D polyester and PVC coated. Supplied with carrying bag. Size: 200x70 cm (without handles), 203x80 cm (with handles). Weight: 2.3 kg. Max. load: 159 kg. Folded size: 80x35 cm.