


ALUMINIUM AND PLASTIC SCOOPING STRETCHERS



56030

Supplied with 2 safety belts



56035

Multilingual manual (both): GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic

- 56030 GIMA SCOOPING STRETCHER - aluminium - silver
 - 56035 GIMA SCOOPING STRETCHER - HPDE - yellow
- Gima scooping stretchers allow operators to raise patient without having to move the body, minimising the risk of further injury.
- Detachable design for transferring and immobilize fractured and severely injured patients.
- The clutch hinges at both ends of the stretcher allow to split the main frame vertically into two parts.
- mainly used for hospitals, sport centres and ambulances
 - adjustable length according to patient's body
 - bottom end (feet) has a narrow frame structure
 - thermal stability of polymer surface (56035 only) for better comfort.
- Supplied with 2 polypropylene (56030) or 3 polyester safety belts, with plastic (56030) or metal (56035) buckle.

TECHNICAL SPECIFICATIONS

	56030	56035
Open size cm:	168/210x42x6	160/210x44x7
Folded size cm:	120x42x8	120x44x7
Weight:	7.4 kg	8 kg
Max. load:	≤159 kg	≤159 kg

HEAD IMMOBILIZERS

HEAD IMMOBILIZERS

Adult and child models of head Immobilizers, different in shape, material and colour.

All are composed of two head blocks and two straps to prevent any head movement of the injured person during transport.


Both sides have an ear hole.

They can be used with scooping and PE stretchers.



34108


Waterproof plastic coating makes it easy to clean and prevents bacterial growth



34109

CT, MRI, X-ray compatible

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.



34109

GIMA code	Type	Colour	Unfolded size (cm)	Weight (kg)	Material
34108	FERMO 3 - child	blue	26x34x14	0.5	high density waterproof plastic
34109	FERMO 3 - adult	red	40x26x18	1	high density waterproof plastic

IMMOBILISATION BELTS FOR STRETCHERS AND SPINAL BOARDS



33838

33840

33842

33843

33844

32794

34081

33837

33847 set of 4



33849

76 cm
90 cm 90 cm 90 cm

Adjustable spider strap for child and adult total length 173 cm

GIMA code	Type	Colour	Size (cm)		Belt material	Buckle	Compatibility
			Width	Length			
33836	Safety belt*	orange	4.8	180	polyester	metal buckle	PE stretchers
33837	Safety belt*	orange	4.8	180	polypropylene	plastic buckle	PE stretchers
33838	Safety belt*	orange	4.8	128+50***	polyester	metal buckle	PE stretchers
33840	Safety belt B2*	black	5.0	180	polypropylene	plastic buckle	PE stretchers
33842	Safety belt B4*	black	4.7	180	polyester	plastic buckle	PE stretchers
33843	Safety belt B8*	black	4.7	143+25***	polyester	plastic buckle, end with hooks	PE stretchers
33844	Safety belt B9*	black	4.7	180	polyester	metal buckle	PE stretchers
33847	Set of 4 safety belts*	mix	4.5	143+25***	polyester	metal buckle, end with hook	PE stretchers, boards
33849	Adjust. Ragno/Spider*	mix	-	-	polypropylene	-	spinal boards
34081	Set 3 belts	orange	5.0	95+47***	-	quick release plastic buckle	stretchers + boards
32794	Nylon bag - empty	red	-	-	-	-	contains up to 3 belts

*Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic. **Only 33836-8 ***Two parts belt