

SILICONE NURSE WATCHES



NURSE WATCHES

Silicone nurse watches, provided with a safety pin.

Suitable to measure patient pulse, beat and sterilization procedure timing.

Conforming to procedures not allowing nurses to wear wrist watches.

Watch case 35 mm, watch length 85 mm.

Shelf life of battery: 12-15 months from first use.

Manual: GB, FR, IT, ES, PT, DE, GR, PL, Arabic.



SILICONE WATCHES		SILICONE WATCHES	
GIMA code	Round - single colour	GIMA code	Round - single colour with aluminium box
20400	White	20410	White
20401	Blue	20411	Blue
20402	Red		Square - single colour with aluminium box
20403	Purple		
20404	Pink	20455	Blue
20405	Black		
	Round - multicolour		Square - multicolour
20420	Red & white	20460	Violet and sky blue
20421	Patchwork	20461	Hearts
20422	Zebra	20462	Flowers
20423	Bubbles	20463	Leopard
20424	Leaves	20464	Roses
20425	Caterpillar	20465	Foot print
	Silicone straps - spare		Digital**
20418	Kit of 6 silicone straps*	20480	Blue

Spare silicone straps for round watches* *Manual: GB, FR, IT, ES*

*Spare silicone straps for round watches

**Manual: GB, FR, IT, ES

20480



NURSE STAINLESS STEEL SCISSORS AND FORCEPS



Nurse bandage scissors - 14 cm



Straight artery forceps - 16 cm



Straight scissors blunt/sharp - 14 cm

*These powder coated instruments are not autoclavable, can be sterilized in water at maximum 110°C

Model	BANDAGE* 14 cm	SCISSORS 18 cm	STRAIGHT* SCISSORS - 14 cm	ARTERY* FORCEPS - 16 cm
Standard	34129**	34130	26725***	26717
Standard with belt clip	20500	-	20520	-
Green ring	20501	-	20521	20541
Lilac ring	20502	-	20522	20542
Foot print	20509	-	20529	20549
Leopard	20511	-	20531	20551
Bubbles	20512	20492	20532	20552
Hearts	20513	-	20533	20553
Caterpillar	20514	20494	20534	-

**15 cm

***14.5 cm

STAINLESS STEEL UTILITY AND BANDAGE SCISSORS



Plastic coloured finger-contour handles for maximum comfort.

All purpose, first aid scissors, easily cut bandages, leather and plaster casts. Autoclavable at 136°C. Instructions: GB, FR, IT, ES, PT, DE.

UTILITY SCISSORS - box of 10				
12 cm - 4.7"	16.5 cm - 6.5"	19 cm - 7.5"	Colour	
20560	20570	20580	Black	
-	20571	20581	Blue	
-	20572	20582	Red	
-	20573	20583	Purple	
-	20574	20584	Pink	
20590	Steel chain for scissors and forceps			




20590



BI-LAYER NON-WOVEN FABRIC STERILE SURGICAL DRAPES

DRAPES



All surgical drapes are made of light-blue bi-layer water-proof non-woven fabric 54 g/m². Individually packed in PE bag, sterile (Ethylene oxide), latex free


GIMA code	DRAPES	Minimum order
23600	Drape 50x50	box of 100
23601	Drape 50x75	box of 100
23602	Drape 75x90	box of 75
23604	Drape 100x100	box of 50
23606	Drape 120x150	box of 50
23608	Drape 150x200	box of 25

DRAPES WITH ADHESIVE SIDE




GIMA code	DRAPES WITH ADHESIVE SIDE	Minimum order
23610	Drape 50x50	box of 100
23611	Drape 50x75	box of 100
23612	Drape 75x90	box of 75
23614	Drape 100x100	box of 50
23617	Drape 180x200	box of 25
23618	Drape 150x240	box of 25

DRAPES WITH CENTRAL ADHESIVE HOLE



Drape with hole Ø 10 cm




Drape with hole 6x9 cm

GIMA code	DRAPES WITH CENTRAL ADHESIVE HOLE	Minimum order
23620	Drape 50x75 - hole 6x9 cm	box of 75
23622	Drape 75x90 - hole 6x9 cm	box of 75
23623	Drape 75x90 - hole Ø 10 cm	box of 75
23624	Drape 100x100 hole Ø 10 cm	box of 50
23625	Drape 100x150 - hole 9x18 cm	box of 25
23626	Drape 120x150 - hole Ø 10 cm	box of 25


NON-WOVEN FABRIC STERILE OPHTHALMIC SURGERY DRAPES

Hole 6x9 cm with incision film.
Single bag to receive fluids and adaptable bar.



- 23634 OPHTHALMIC SURGERY DRAPE 80x100 cm - sterile - box of 40

Hole 10x10 cm with incision film. Double bag to receive fluids with adaptable bar



- 23635 OPHTHALMIC SURGERY DRAPE 160x160 cm - sterile - box of 15

Hole Ø11 cm with incision film.
Double bag 35x73 cm to receive fluids with adaptable bar and adaptable bar for nasal septum. Wrapped in crepe paper.



- 23636 OPHTHALMIC SURGERY DRAPE 160x260 cm - sterile - box of 10


SURGERY COTTON DRAPES



- 26495 DRAPE 90x150 cm
 - 26496 DRAPE 150x150 cm
 - 26497 DRAPE 250x150 cm
- Made of 100% high quality cotton, green colour.
Made in Europe.

INCISION FILM - STERILE

Its purpose is to define the anatomic part of the patient where the surgery will be done.
It is applied at the point where the wound is, to protect it from possible contamination.
Incision film made of adhesive polyurethane.



GIMA code	SIZE AxB	SIZE AxC	Minimum order
23628	45x50 cm	45x40 cm	box of 10
23630	30x35 cm	30x25 cm	box of 50
23631	45x50 cm	45x40 cm	box of 50
23632	60x65 cm	60x55 cm	box of 50

POLYESTER OR MICROFIBRE SURGERY DRAPES


100% TEXTURIZED POLYESTER DRAPES

Studied for replacing the classic cotton drape in operating rooms.

Advantages:

- it does not have particles release of dust or cotton fibers that could cause infections during surgery
- various usability on the field:
 - as a sheet under the patient
 - as a towel
 - to tamponage, considering its high absorbency
 - as a pre or post surgery cover, thanks to the capability of keeping a good body temperature

Polyester texturized fabric



- 41000 POLYESTER DRAPE 90x150 cm - light blue
- 41001 POLYESTER DRAPE 150x150 cm - light blue
- 41002 POLYESTER DRAPE 250x150 cm - light blue

MICROFIBRE DRAPES

(99% polyester + 1% carbon fibre)

The microfibre line is conceived for operating room.

The characteristics and abilities are various:

- the fabric does not absorb liquids, maintaining the operating field dry;
- the fabric does not have particles release of dust or fibre that could cause infections on the patient.

GIMA fabric is composed of polyester and carbon fibres. The addition of carbon fibres serves as protection from potential electric discharges.


- 41009 MICROFIBRE DRAPE 90x150 cm - green
- 41010 MICROFIBRE DRAPE 150x150 cm - green
- 41011 MICROFIBRE DRAPE 250x150 cm - green

Sizes of drapes may vary ±5% due to fabric processing during production phase

Microfibre fabric




STERILE SURGERY SETS



Non woven bi-layer fabric universal set - sterile

• 23640 UNIVERSAL SET - sterile - box of 10


- 1 table cover 150x200 cm with central absorbent band
- 1 Mayo stand cover 80x145 cm
- 1 standard surgical drape, 150x240 cm with adhesive tape
- 1 standard surgical drape, 180x200 cm with adhesive tape
- 2 standard surgical drapes, 75x90 cm with adhesive tape
- 4 towels 35x40 cm
- 1 adhesive strip 10x50 cm
- 1 tube holder 2.5x20 cm



Non woven bi-layer fabric orthopaedy set for knee arthroscopy - sterile

• 23642 SET FOR KNEE ARTHROSCOPY - sterile - box of 8


- 1 table cover 150x200 cm with absorbent central band
- 1 Mayo stand cover 80x145 cm
- 1 standard surgical drape, 150x180 cm
- 1 orthopaedy drape for arthroscopy at knee 230x320 cm with elastic hole and bag for liquid collection
- 1 stockinette 22x75 cm
- 3 adhesive strips 10x50 cm
- 4 towels 35x40 cm



Non woven SMS fabric obstetrics set for caesarean section - sterile

• 23644 SET FOR CESAREAN SECTION - sterile - box of 8

- 1 table cover 150x200 cm with central absorbent band
- 1 Mayo stand cover 80x145 cm
- 1 drape for caesarean section 260x200x310 cm with incision film and bag for liquid collection
- 1 drape "baby blanket" 75x90 cm




• 23650 Table cover 150x200 cm - sterile - box of 20




• 23654 Adhesive pocket for tools 35x42 cm - box of 100



• 23651 Mayo stand cover 80x145 cm - sterile - box of 40




• 23655 Camera drape 13x250 cm - box of 100



Non woven bi-layer fabric gynaecology set

• 23641 GYNAECOLOGY SET - sterile - box of 8


- 1 table cover 150x200 cm with central absorbent band
- 1 Mayo stand cover 80x145 cm
- 1 gynaecology drape for lithotomy position 290x240 cm with adhesive hole, built-in
- leggings, bag for liquid collection with filter and drain valve
- 1 drape 75x90 cm
- 1 incision film 35x50 cm
- 3 adhesive strips 10x50 cm
- 4 towels 35x40 cm



Non woven SMS fabric orthopaedy set for hand surgery - sterile

• 23643 SET FOR HAND SURGERY - sterile - box of 8


- 1 table cover 100x150 cm with central absorbent band
- 1 standard surgical drape, 100x100 cm
- 2 towels 35x40 cm
- 1 T-shape drape for hand surgery 370x190x290 cm with elastic hole 5 cm and cable fixing Velcro strips



Non woven SMS fabric urology set for TUR - sterile

• 23645 UROLOGY SET FOR TUR - sterile - box of 8

- 1 table cover 100x150 cm with central absorbent band
- 1 urology drape for TUR, 200x280x220 cm with perineal and subpubic holes,
- bag for liquid collection with filter, adhesive tape and straps
- 1 camera drape 13x250 cm
- 2 towels 35x40 cm



• 23652 Drape, 200x240 cm with adhesive U shape cut - box of 15



• 23656 Bag to cover instruments 75x90 cm - box of 100



• 23653 Pair of leggings 75x130 cm for urology/gynaec. - box of 20



• 23657 Adhesive strip 10x50 cm - box of 200

TUNGSTEN CARBIDE INSTRUMENTS - GOLD LINE

GOLD LINE - TUNGSTEN CARBIDE INSTRUMENTS

A top quality theatre use Tungsten Carbide Instruments made of stainless steel AISI 410 (AISI 420 for scissors) with very hard tips made in Germany to help surgeons perform precision surgeries with best results. Highly reliable for hardness and corrosion resistance. Instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.

TUNGSTEN CARBIDE NEEDLE HOLDERS - AISI 410



26544 Webster needle holder - 12.5 cm



26546 Crile Wood needle holder - 15 cm

26547 Crile Wood needle holder - 18 cm

26545 Crile Wood needle holder - 20 cm



26539 Mayo Hegar needle holder - 14 cm

26530 Mayo Hegar needle holder - 16 cm

26531 Mayo Hegar needle holder - 18 cm

26532 Mayo Hegar needle holder - 20 cm

26535 Mayo Hegar needle holder - 24 cm

26533 Mayo Hegar needle holder - 26 cm

26534 Mayo Hegar needle holder - 30 cm



26540 Olsen Hegar needle holder - 14 cm

26541 Olsen Hegar needle holder - 16.5 cm


26542 Olsen Hegar needle holder - 19 cm



26536 Mathieu needle holder - 14 cm

26537 Mathieu needle holder - 17 cm

26538 Mathieu needle holder - 20 cm



26548 Gillies needle holder - 16 cm



26549 Derf needle holder - 12 cm

TC MICRO NEEDLE HOLDER - AISI 410



26550 Castroviejo Needle Holder straight - 14 cm - plain tip

26552 Castroviejo Needle Holder curved - 14 cm - plain tip

26553 Castroviejo Needle Holder straight - 14 cm - rough tip

26555 Castroviejo Needle Holder curved - 14 cm - rough tip

26551 Castroviejo Needle Holder straight - 18 cm - rough tip

26556 Castroviejo Needle Holder straight - 21 cm - rough tip

26557 Castroviejo Needle Holder curved - 21 cm - rough tip



26554 Castroviejo - 13 cm



Smooth

38600 Barraquer Micro needle holder - 13 cm - smooth tips - with lock



38610 Arruga needle holder - straight - 16 cm plain tips

38611 Arruga needle holder - curved - 16 cm plain tips

TUNGSTEN CARBIDE FORCEPS - AISI 410



26560 Adson - 12 cm

26561 Adson 1x2 teeth - 12 cm



26562 Micro Adson forceps - 12 cm - 1x2 teeth

26563 Micro Adson forceps - 15 cm - 1x2 teeth



26558 Gillies Dissecting forceps - 15 cm - 1x2 teeth



26565 Potts Smith Dissecting forceps - 15 cm

26566 Potts Smith Dissecting forceps - 18 cm

26567 Potts Smith Dissecting forceps - 20 cm

26568 Potts Smith Dissecting forceps - 23 cm

26569 Potts Smith Dissecting forceps - 25 cm

LONG LASTING JOINT GERMAN TIPS



38552 Cushing forceps - 18 cm



38554 Cushing Taylor forceps - 18 cm - rough tips



38560 Gerald forceps - 17 cm - rough tips



38565 Mc Indoe forceps - 15 cm - rough tips

TC WIRE CUTTERS - AISI 410



26570 Wire Cutter - 14 cm - for soft wires 0-1 mm

26571 Wire Cutter - 18 cm - for hard wires cutting up to 1.6 mm

26572 Wire Cutter - 23 cm - for hard wires cutting up to 2.0 mm



26575 Wire holding - 18 cm

TC ORGAN HOLDING FORCEPS - AISI 410



26564 Duval - 20 cm

TUNGSTEN CARBIDE INSTRUMENTS - GOLD LINE

TUNGSTEN CARBIDE SCISSORS - AISI 420

DOUBLE JOINT GERMAN TIPS

- 26500 Straight scissors blunt/sharp - 14.5 cm
 26501 Straight scissors blunt/blunt - 14.5 cm
 26502 Straight scissors sharp/sharp - 14.5 cm
 26503 Curved scissors blunt/sharp - 14.5 cm
 26504 Curved scissors blunt/blunt - 14.5 cm
 26505 Curved scissors sharp/sharp - 14.5 cm

- 26526 Iris scissors - straight - 11.5 cm
 26527 Iris scissors - curved - 11.5 cm

- 26528 Goldman Fox scissors - straight - 13 cm
 26529 Goldman Fox scissors - curved - 13 cm

- 26506 Universal scissors - 12 cm

- 26510 Metzenbaum scissors - straight - 14 cm
 26511 Metzenbaum scissors - straight - 18 cm
 26512 Metzenbaum scissors - straight - 20 cm
 26514 Metzenbaum scissors - curved - 14 cm
 26515 Metzenbaum scissors - curved - 18 cm
 26516 Metzenbaum scissors - curved - 20 cm

- 26520 Mayo scissors - straight - 14.5 cm
 26521 Mayo scissors - straight - 18 cm
 26524 Mayo scissors - straight - 23 cm
 26522 Mayo scissors - curved - 14.5 cm
 26523 Mayo scissors - curved - 18 cm
 26525 Mayo scissors - curved - 23 cm

- 38530 Lister scissors - 14 cm
 38532 Lister scissors - 18 cm

- 38510 Mayo Stille scissors - straight - blunt - 15 cm
 38511 Mayo Stille scissors - curved - blunt - 15 cm
 38514 Mayo Stille scissors - straight - blunt - 17 cm
 38515 Mayo Stille scissors - curved - blunt - 17 cm

- 38520 Kelly scissors - straight - sharp - 16 cm
 38521 Kelly scissors - curved - sharp - 16 cm

- 26507 Debakey Potts Smith scissors - 19 cm - angle 25°
 26508 Debakey Potts Smith scissors - 19 cm - angle 45°
 26509 Debakey Potts Smith scissors - 19 cm - angle 60°
 26517 Debakey Potts Smith scissors - 23 cm - angle 25°
 26518 Debakey Potts Smith scissors - 23 cm - angle 45°
 26519 Debakey Potts Smith scissors - 23 cm - angle 60°

KEVORKIAN - S/S



MINI-TISCHLER - S/S



TISCHLER - S/S




SUPER CUT SCISSORS



- 26488 METZENBAUM SCISSORS SUPER CUT - curved - 15 cm
 • 26489 METZENBAUM SCISSORS SUPER CUT - curved - 18 cm

SUPER CUT SCISSORS
 Super Cut scissors have a specially designed razor-sharp upper blade edge that cuts effortlessly through tissue. The lower blade has micro-fine serrations to hold tissue and prevent slippage.

Tungsten carbide honed insert combined with razor-sharp edge



- 26485 MAYO SCISSORS SUPER CUT WITH T.C. - curved - 14 cm
 • 26486 MAYO SCISSORS SUPER CUT WITH T.C. - curved - 17 cm
 • 26487 MAYO SCISSORS SUPER CUT WITH T.C. - curved - 20 cm

STAINLESS STEEL NEEDLE HOLDERS

CRILE WOOD NEEDLE HOLDER



- 26852 Crile Wood needle holder - 15 cm
26878 Crile Wood needle holder - 18 cm
26879 Crile Wood needle holder - 20 cm

MAYO HEGAR NEEDLE HOLDER



- 26818 Mayo Hegar needle holder - 14 cm
26732 Mayo Hegar needle holder - 16 cm
26819 Mayo Hegar needle holder - 18 cm
26733 Mayo Hegar needle holder - 20 cm

MATHIEU NEEDLE HOLDER



- 26741 Mathieu needle holder - 14 cm
26730 Mathieu needle holder - 16 cm
26731 Mathieu needle holder - 20 cm

STAINLESS STEEL SURGICAL SCISSORS

BLUNT/SHARP SCISSORS



- 26754 Straight scissors blunt/sharp - 11.5 cm
26725 Straight scissors blunt/sharp - 14.5 cm
26729 Straight scissors blunt/sharp - 16 cm
26853 Straight scissors blunt/sharp - 18 cm
26726 Straight scissors blunt/sharp - 20 cm
26755 Curved scissors blunt/sharp - 11.5 cm
26727 Curved scissors blunt/sharp - 14.5 cm
26756 Curved scissors blunt/sharp - 16 cm
26854 Curved scissors blunt/sharp - 18 cm
26728 Curved scissors blunt/sharp - 20 cm

BLUNT/BLUNT SCISSORS



- 26855 Straight scissors blunt/blunt - 11.5 cm
26740 Straight scissors blunt/blunt - 14.5 cm
26856 Straight scissors blunt/blunt - 16 cm
26857 Straight scissors blunt/blunt - 18 cm
26887 Straight scissors blunt/blunt - 20 cm
26888 Curved scissors blunt/blunt - 11.5 cm
26742 Curved scissors blunt/blunt - 14.5 cm
26889 Curved scissors blunt/blunt - 16 cm
26890 Curved scissors blunt/blunt - 18 cm
26891 Curved scissors blunt/blunt - 20 cm

SHARP/SHARP SCISSORS



- 26892 Straight scissors sharp/sharp - 11.5 cm
26744 Straight scissors sharp/sharp - 14.5 cm
26893 Straight scissors sharp/sharp - 16 cm
26894 Straight scissors sharp/sharp - 18 cm
26895 Straight scissors sharp/sharp - 20 cm
26896 Curved scissors sharp/sharp - 11.5 cm
26746 Curved scissors sharp/sharp - 14.5 cm
26897 Curved scissors sharp/sharp - 16 cm
26898 Curved scissors sharp/sharp - 18 cm
26899 Curved scissors sharp/sharp - 20 cm

LISTER BANDAGE SCISSORS



- 34125 Lister scissors - 11 cm
34129 Lister scissors - 16 cm
34130 Lister scissors - 18 cm

MAYO-STILLE SURGICAL SCISSORS




- 26846 Mayo-stille scissors - straight - 14.5 cm
26847 Mayo-stille scissors - straight - 18 cm
26848 Mayo-stille scissors - straight - 20 cm
26849 Mayo-stille scissors - curved - 14.5 cm
26850 Mayo-stille scissors - curved - 18 cm
26851 Mayo-stille scissors - curved - 20 cm

METZENBAUM SURGICAL SCISSORS



- 26840 Metzenbaum scissors - straight - 14.5 cm
26841 Metzenbaum scissors - straight - 18 cm
26843 Metzenbaum scissors - curved - 14.5 cm
26844 Metzenbaum scissors - curved - 18 cm

STITCH REMOVER FORCEPS/SCISSORS



- 26737 Spencer scissors - 13 cm for suture


- 26722 Michel stitch remover forceps - 12 cm

IRIS SCISSORS



- 26822 Iris scissors - straight - 11.5 cm
26824 Iris scissors - curved - 11.5 cm

MICRO RIBBON SCISSORS



- 26956 Ribbon scissors - straight - 9.5 cm
26957 Ribbon scissors - curved - 9.5 cm
26958 Ribbon scissors - straight - 12 cm
26959 Ribbon scissors - curved - 12 cm


DELIVERY SCISSORS



- 26827 Umbilical scissors Usa model - 10.5 cm


- 26874 Brown Sadler episiotomy scissors - 14 cm

DERMATOLOGY PROBES






- 26825 Katsch - 14 cm
26826 Gillies hook - 13 cm - small
26821 Unna round - 14 cm
26820 Unna Oval - 14 cm
26823 Round forceps - 9 cm

PROBES, SCALPELS



- 26698 Scalpel straight blade - 13 cm
26703 Scalpel straight blade - 17.5 cm
26699 Scalpel pot-bellied blade - 13 cm
26704 Scalpel pot-bellied blade - 17.5 cm
26701 Probe D-end - 15 cm - Ø 2 mm
26702 Butterfly probe - 14.5 cm

STAINLESS STEEL SURGICAL FORCEPS

TOWEL FORCEPS	HAEMOSTATIC FORCEPS	DRESSING SPONGE HOLDING FORCEPS
 26743 Backhaus forceps - 9 cm 26736 Backhaus forceps - 11 cm 26747 Backhaus forceps - 13.5 cm  26739 Backhaus Forceps for TNT - 13 cm	 26714 Mosquito forceps - straight - 12.5 cm 26872 Mosquito forceps - straight - 14 cm 26715 Mosquito forceps - straight - 16 cm 26716 Mosquito forceps - curved - 12.5 cm 26873 Mosquito forceps - curved - 14 cm	 26745 Pean forceps - 14 cm 26720 Pean forceps - 16 cm 26710 Pean forceps - 20 cm
DISSECTING FORCEPS		
 26690 Anatomy forceps - 12 cm 26705 Anatomy forceps - 14 cm 26691 Anatomy forceps - 16 cm 26706 Anatomy forceps - 18 cm 26752 Anatomy forceps - 20 cm  26692 Surgery forceps - 12 cm - 1x2 teeth 26707 Surgery forceps - 14 cm - 1x2 teeth 26693 Surgery forceps - 16 cm - 1x2 teeth 26708 Surgery forceps - 18 cm - 1x2 teeth 26753 Surgery forceps - 20 cm - 1x2 teeth  <i>Adson without teeth</i> 26750 Adson forceps - 12 cm 26751 Adson forceps - 12 cm - 1x2 teeth 26694 Micro Adson forceps - 12 cm 26695 Micro Adson forceps - 12 cm - 1x2 teeth  26724 Forceps thin end - 12 cm  26977 English Toe plain dissecting forceps - 14 cm	 26757 Kelly forceps - straight - 14 cm 26717 Kelly forceps - straight - 16 cm 26758 Kelly forceps - curved - 14 cm 26718 Kelly forceps - curved - 16 cm  <i>Straight</i> <i>Curved</i> 26828 Kocher forceps - straight - 14 cm - 1x2 26711 Kocher forceps - straight - 16 cm - 1x2 26829 Kocher forceps - straight - 18 cm - 1x2 26712 Kocher forceps - straight - 20 cm - 1x2 26831 Kocher forceps - curved - 14 cm - 1x2 26832 Kocher forceps - curved - 16 cm - 1x2 26833 Kocher forceps - curved - 18 cm - 1x2 26834 Kocher forceps - curved - 20 cm - 1x2  26749 Klemmer forceps - straight - 14 cm 26709 Klemmer forceps - straight - 16 cm 26835 Klemmer forceps - straight - 18 cm 26713 Klemmer forceps - straight - 20 cm 26836 Klemmer forceps - curved - 14 cm 26837 Klemmer forceps - curved - 16 cm 26838 Klemmer forceps - curved - 18 cm 26839 Klemmer forceps - curved - 20 cm  26978 Mixer forceps - 23 cm	 26800 Foerster forceps - 20 cm 26801 Foerster forceps - 25 cm  26805 Foerster "small ring" forceps - 25 cm  26845 Duval 14 cm x 10 mm  <i>Curved</i> <i>Straight</i> 26734 Crile forceps - straight - 16 cm 26735 Crile forceps - curved - 16 cm  26721 Allis forceps - 15 cm - 4x5 teeth 26759 Allis forceps - 19 cm - 5x6 teeth  26986 Cheron forceps - 25 cm
ANESTHESIA FORCEPS	OB/GYN FORCEPS	DRESSING SWAB FORCEPS
 34005 Collin Forceps - 20 cm  34001 Magill Forceps - 15 cm - infants 34002 Magill Forceps - 20 cm - children 34003 Magill Forceps - 25 cm - adults	 26979 Museux Vulsellum forceps 8 mm - straight - 24 cm	 26987 Rampléy forceps - 18 cm 26988 Rampléy forceps - 24 cm  26811 Bozeman dressing forceps straight - 26 cm 26812 Bozeman dressing forceps curved - 26 cm
	BONE CUTTING FORCEPS	
	 26983 Luer forceps - 15 cm - sharp edge  26984 Ruskin Liston forceps - 18 cm	

S/S INSTRUMENT KITS



26768

GIMA code	INSTRUMENTS KIT IN ALUMINIUM BOX OR POLYESTER BAG	26701	26702	26703	26704	26705	26707	26709	26722	26723	26724	26725	26730	26732
26760	Trousse standard - polyester bag 9 instr.	•	•	•	•	•	•	•	•	•	•	•	•	•
26761	Trousse classica - polyester bag 10 instr.	•	•	•	•	•	•	•	•	•	•	•	•	•
26762	Trousse suprema - polyester bag 11 instr.	•	•	•	•	•	•	•	•	•	•	•	•	•
26768	Trousse suprema - aluminium box 11 instr.	•	•	•	•	•	•	•	•	•	•	•	•	•

PROBE	BUTTERFLY BLADE	STRAIGHT BLADE	POT-BELLIED BLADE	ANATOMY FORCEPS	SURGERY FORCEPS	KLEMMER FORCEPS	STITCH REMOVER	HOOK SHARP	FORCEPS THIN END	SCISSORS BLUNT/SHARP	MATHIEU NEEDLE HOLDER	MAYO-HEGAR HOLDER
26701	26702	26703	26704	26705	26707	26709	26722	26723	26724	26725	26730	26732

S/S DRESSING, SUTURE AND DELIVERY SURGERY SETS



• 26963 DRESSING SET - 8 PIECES



• 26965 SUTURE SET - 12 PIECES



• 26966 MINOR DELIVERY SET - 25 PIECES

SET CONTENTS	26963	26965	26966
Kidney Dish 20 cm (8") - medium	•	•	•
Kidney Dish 25 cm (10") - large	•	•	•
Gallipot 150 ml - 6 oz - medium	•	•	•
Mayo Hegar needle holder - 18 cm	•	•	•
Brown Sadler episiotomy scissors - 14 cm	•	•	•
Kocher artery forceps - straight - 16 cm	•	•	•
Kocher artery forceps straight - 18 cm	•	•	•
Kocher artery forceps curved - 16 cm	•	•	•
Dressing Forceps (Thumb) - 16 cm	•	•	•
English TOE dissecting forceps - plain - 14 cm	•	•	•
English TOE dissecting forceps - plain - 18 cm	•	•	•
Dissecting forceps toothed 1x2 - 14 cm	•	•	•
Suture scissors - Heaths Ligature - 15 cm	•	•	•
Mayo scissors - straight - 15 cm	•	•	•
Operating scissors sharp/blunt - straight - 15 cm	•	•	• 2 pcs.
Operating scissors sharp/blunt - curved - 15 cm	•	•	• 2 pcs.
Ramsey sponge holding forceps - str. - 18 cm	•	•	•
Ramsey sponge holding forceps - str. - 24 cm	•	•	• 2 pcs.
Museux Vulsellum forceps 8 mm - str. - 24 cm	•	•	•
Cusco vaginal specula - small	•	•	•
Cusco vaginal specula - medium	•	•	•
Cusco vaginal specula - large	•	•	•
Sims vaginal specula - small	•	•	•
Sims vaginal specula - medium	•	•	•
Sims vaginal specula - large	•	•	•
Umbilical scissors USA Patent - 10.5 cm	•	•	•
Instrument tray with lid 25x30x6.2 cm	•	•	•
Instrument tray with lid 30x20x5 cm	•	•	•
Instrument tray with lid 45x30x6.2 cm	•	•	•

GYNAECOLOGY INSTRUMENTS AND KIT



26807 Novak suction cannula - 23 cm -

Ø 4 mm

26799 Adaptor for Novak cannula

26804 ★ Sims hystrometer - malleable - 33 cm - silver

26816 Martin hystrometer - 30 cm - malleable - brass handle, copper wire

26803 ★ Bosch/Sims scissors - 23 cm

26802 ★ Schroeder forceps - 25 cm (Pozzi-Palmer)

26801 ★ Foerster forceps - 25 cm

• 26810 GYNAECOLOGY KIT
Includes 4 instruments marked with ★ in a nylon bag.



26817

26817 Set of chrome plated Hegar Dilators
Include 14 dilators Ø 4-17 mm

60x25 mm
70x30 mm70x30 mm
70x35 mm85x35 mm
90x40 mm

• 26884 SET OF 3 SIMS RETRACTORS




26810

Rigid cannula with sharp tip

26813 Sims uterine curette - sharp - 8 mm - 26 cm

26814 Sims uterine curette - sharp - 9 mm - 26 cm


26815 Sims uterine curette - sharp - 14 mm - 26 cm



26806 Schumacher biopsy forceps - 24 cm

26808 Schubert biopsy forceps - 21 cm

26809 Kogan endospeculum - 24 cm



26811 Bozeman dressing forceps straight - 26 cm

26812 Bozeman dressing forceps curved - 26 cm

ENT STAINLESS STEEL INSTRUMENTS AND KIT



Tongue depressors

- 25620 Mayo tongue depressor - 17 cm
 25621 Bruenings tongue depr - fenestrated - 19 cm ●
 25623 Tobold tongue depressor



Nose/ear speculum

- 26774 Voltolini nose speculum - adjustable - Size 2 ●
 26775 Voltolini nose speculum - adjustable - Size 0
 26770 Set of 3 ear speculum Hartmann - Ø 4, 5, 6 mm ●
 26784 Thudicum nose clip




Nose Ear forceps and scissors

- 26773 Lucae ear forceps - 14 cm ●
 26783 Troeltsch ear forceps - 12 cm
 26779 Hartmann ear forceps - 14 cm
 26781 Heymann scissors - 18 cm



Suction cannulas

- 26795 Frazier nose/ear suction cannula Ø 2 mm (18 cm)
 26796 Frazier nose/ear suction cannula Ø 3 mm (18 cm)
 26793 Frazier nose suction cannula Ø 4 mm (18 cm)
 26794 Frazier nose suction cannula Ø 5 mm (18 cm)
 26797* Rosen ear suction cannula Ø 1.5 mm (15.5 cm)
 26798* Rosen ear suction cannula Ø 3 mm (16.1 cm)
 *with adaptor with LUER tip and thumb valve
 26799 Adaptor with Leur tip for 26797-8 - spare



Hooks-Loops

- 26771 Lucae Hook - 14 cm (5'4") ●
 26769 Lucae Hook - fine - 14 cm (5'4")
 26790 Billeau Loop - 16 cm
 26785 Swab holder - 16 cm

● 26788 MICRO EAR CUP-SHAPED POLYPUS FORCEPS

A = 80 mm, B = 4 mm



● 26789 MICRO EAR SCISSORS BELLUCCI

A = 80 mm, B = 4 mm
 C = 0,8 mm tip thickness




● 26791 MICRO EAR PUNCH - introduction instrument

A = 80 mm, B = 4 mm



Killian/Hartmann nose forceps

- 26792 Killian/Hartmann nose forceps - 14 cm (5'4")/15 mm
 26778 Killian/Hartmann nose forceps - 14 cm (5'4")/22 mm
 26772 Killian/Hartmann nose forceps - 14 cm (5'4")/35mm ●
 26776 Killian/Hartmann nose forceps - 14 cm (5'4")/50mm
 26764 Killian/Hartmann nose forceps - 14 cm (5'4")/75mm
 26777 F. O. Killian nose forceps - 14 cm/75mm fiber optic
 26765 Tieck Halle nose speculum - 13.5 cm



Hartmann ear forceps

GIMA code	mm			Picture
	A	B	C	

26782	90	11	0.6	
-------	----	----	-----	--

Hartmann micro forceps

26786	80	6	0.3	
-------	----	---	-----	--

26787	80	4	0.6	
-------	----	---	-----	--




- 26780 ENT KIT Includes all instruments marked with ● inside an aluminium box (Code 26665) 21x10xh 5 cm.

LARYNGEAL MIRRORS AND SPIRIT LAMP



31614



31612

- 31616 ADJUSTABLE MIRROR HANDLE - 12.5 cm
 ● 31615 ADJUSTABLE MIRROR HANDLE - 8 cm

31615+31606

31616+31605

GIMA code	LARYNGEAL MIRROR - STRAIGHT - 13 cm	Minimum order
31600	Mirror number 000 diameter 10 mm	10 pcs. mixed sizes
31601	Mirror number 00 diameter 12 mm	
31602	Mirror number 0 diameter 14 mm	
31603	Mirror number 1 diameter 16 mm	
31604	Mirror number 2 diameter 18 mm	
31605	Mirror number 3 diameter 20 mm	
31606	Mirror number 4 diameter 22 mm	
31607	Mirror number 5 diameter 24 mm	
31608	Mirror number 6 diameter 26 mm	

● 31612 SET 9 LARYNGEAL MIRRORS 000-6

Set of 9 laryngeal mirrors from 000 to 6, all mirrors have a different diameter: from 10 mm to 26 mm.
 Instructions: GB, FR, IT, ES, PT.

minimum order: 1 set of 9 mirrors

● 31614 SPIRIT LAMP - 100 ml

Lamp is made of AISI 304 stainless steel.
 Brass cap to close flame area to avoid discoloration.

GOUGE BLADES AND HANDLES

CARBON STEEL GOUGE BLADES - sterile

Sterile, sharp, surgical gouge blades are used by podologists for hand, feet and skin care. Each blade is sterilized by gamma



GIMA code	CARBON STEEL GOUGE BLADES - STERILE	code	
22501	Blades 1 - box of 50	22506	Blades 6 - box of 50
22502	Blades 2 - box of 50	22508	Blades 8 - box of 50
22503	Blades 3 - box of 50	22510	Blades 10 - box of 50
22504	Blades 4 - box of 50	22512	Blades 12 - box of 50
22505	Blades 5 - box of 50	22515	Blades 15 - box of 50

radiation and protected by VCI paper in a peel - open foil sachet. Gouge blades fit into handle no. 3 and special gouge handle. Wide range of sizes.



GIMA code	STAINLESS STEEL HANDLES
22520	Special gouge handle - stainless steel
26913	Gouge/scalpel handle N°3 - stainless steel - see page 61


*Compatible with blades from 10 to 15

PODIATRY-CHIROPODY STAINLESS STEEL INSTRUMENTS

SCISSORS




TWEEZERS



PEDICURE FILES



CUTTERS/NIPPERS



PODIATRY STAINLESS STEEL KITS



• 39800 HEART SHAPE KIT - blue

Cuticle nipper 10 cm, Nail nipper 10 cm, Nail scissors 9.5 cm, Nail file 10 cm, Pusher double ended 12 cm, Tweezers 10 cm

• 39805 STEEL FRAME KIT

Nail nipper 10 cm, Nail Scissor 9.5 cm, Cuticle scissors 9.5 cm, Cuticle pusher double ended 12 cm, Cuticle pusher and cleaner 13 cm, Tweezers 9 cm, Nail file with white plastic handle 9 cm



• 39810 PROFESSIONAL KIT

Nail nipper Moon Shape 14 cm, Side cutter 13 cm, Cuticle Nipper 10 cm, Dressing Scissors 14 cm, B.P. handle No. 4 13.5 cm, Splinter Tweezers 12 cm, Loop 14 cm, Knife & Cleaner double ended 16 cm, Nail file 12 cm, Nail file with white plastic handle 13 cm, Spatula double ended 15 cm



DENTAL STAINLESS STEEL INSTRUMENTS

EXTRACTING-FORCEPS - FOR JAWS



60040
set of 10 perforated
impression trays
(5 lower + 5 upper)



60045
set of 10 solid impression
trays (5 lower + 5 upper)
Above impression trays sets
have below sizes

UPPER SIZES mm		LOWER SIZES mm	
U1	U2	L1	L2
62	60	63	61
53	54	57	60
52	52	53	53
U3	U4	L3	L4
73	69	78	74
68	68	71	71



fig. 1
60000
upper
17.5 cm



fig. 2
60001
upper
17.5 cm



fig. 7
60002
upper
17.5 cm




fig. 17
60003
upper
17.5 cm



fig. 18
60004
upper
17.5 cm



fig. 29
60005
upper
17.5 cm



fig. 51
60006
upper
18.5 cm




fig. 67A
60007
upper
18.5 cm




fig. 16
60015
lower
18.5 cm



fig. 22
60016
lower
15 cm



fig. 74N
60017
lower
15.1 cm



fig. 79
60018
lower
17.5 cm



fig. 151
60019
lower
17 cm

EXPLORER AND PROBES



60101 Explorer - fig. 17/23



60100 Explorer probe



60105 WHO probe

GRACEY AND Mc CALL CURETTES



60120 Mc Call curettes - fig. 13/14



60125 Gracey curette - fig. 1/2 ant.



60127 Gracey curette - fig. 7/8 medium



60130 Gracey curette - fig. 11/12 medium



60131 Gracey curette - fig. 13/14 post.

SCALER



60140 Scaler - fig. 1




60141 Mc Call scaler - fig. 18



60142 Mitchel scaler



60144 Langer scaler - fig. 1/2




60145 Langer scaler - fig. 3/4



60148 Hygienist scaler - fig. H6/7



60150 Scaler - fig. 2045



60155 Periodontia instr. - fig. ST 2/3



60156 Periodontia instr.


CARVER



60170 Cement spatula



60171 Wiland modelling carver



60172 Interproximal modelling carver



60173 Beale modelling carver - mod. 7A



60174 Fahrenstock - modelling carver - 13 cm

FORCEPS




60190 London college forceps - 15 cm

CARTRIDGE SYRINGE



60195 Cartridge Syringe - 2.2 cc
with thread in millimeter with aspiration.
Compatible with Gima and Terumo
hypodermic needles page 69-70



3 mm

60200



4 mm

60201



60202



60203



60204 **60205**

Winter
root
elevator




60210

Amalgam
carrier



Molt
Mouth gag
60220 - 11 cm
60221 - 14 cm



60225

Bone Cutter
- 15 cm



60226

Beyer Forceps
- 18 cm

■ PEHA® DISPOSABLE STAINLESS STEEL INSTRUMENTS - STERILE

Single-use instruments made of matt brushed, high-quality stainless steel. Made in EU, they offer similar characteristics as multi-use instruments in terms of intended performance, haptics and design.

Sterile, individually packed and ready-to-use for a variety of surgical interventions. They allow precise handling during procedures. Peha® color lock® safety marking allows quick and reliable differentiation from multi-use instruments.

The self-adhesive traceability label on the packaging facilitates easy, fast and less error prone data entry into the patient records.

Excellent choice for health care facilities with no processing unit.

Recyclable and disposable in existing clinical waste containers after use.

Multilingual box: GB, FR, IT, ES, PT, DE, NL, PL, GR, HU, CZ, SI, RU, Arabic.

HARTMANN

CE 0123



STERILE SCISSORS



38701

- **38700** Operating scissors sharp/sharp straight - 13 cm - 991083 - box of 25
- **38701** Operating scissors sharp/blunt straight - 14.5 cm - 991081 - box of 25
- **38702** Operating scissors blunt/blunt straight - 14.5 cm - 991082 - box of 25



- **38705** Metzenbaum scissors blunt/blunt curved - 14.5 cm - 991084 - box of 25



- **38712** Braun-Stadler episiotomy scissors 14.5 cm - 991080 - box of 20



38708

- **38708** Iris scissors straight - 11.5 cm - 991086 - box of 25
- **38709** Iris scissors curved - 11.5 cm - 991087 - box of 25

STERILE FORCEPS



- **38746** Standard tissue forceps straight - 14 cm - 991063 - box of 25



- **38744** Adson tissue forceps straight - 12 cm - 991062 - box of 25



- **38748** Standard anatomic forceps straight - 14 cm - 991064 - box of 25




- **38740** Adson anatomic forceps straight - 12 cm - 991060 - box of 25

STERILE HAEMOSTATIC FORCEPS




- **38742** DeBakey forceps straight - 15.5 cm - 991061 - box of 25



- **38733** Halsted Mosquito forceps anatomic curved - 12.5 cm - 991041 - box of 25

STERILE NEEDLE HOLDER



38720


- **38720** Mayo-Hegar needle holder 12 cm - 991031 - box of 25
- **38721** Mayo-Hegar needle holder 14 cm - 991030 - box of 25
- **38722** Mayo-Hegar needle holder 16 cm - 991034 - box of 25

STERILE CURETTE



- **38726** Curette sharp spoon double end - 16.5 cm - 991010 - box of 25

STERILE HAEMOSTATIC FORCEPS



- **38730** Pean forceps anatomic straight - 14 cm - 991040 - box of 25



- **38738** Micro-Mosquito forceps anatomic straight - 12.5 cm - 991044 - box of 25



- **38735** Kocher tissue forceps straight - 14 cm - 991042 - box of 25

■ STAINLESS STEEL RETRACTORS



26865

- **26858** Weitlaner-Loktite retractor - sharp 13 cm - 3x4 prongs
- **26859** Weitlaner-Loktite retractor - blunt 13 cm - 3x4 prongs
- **26865** Weitlaner Retractor - sharp 20 cm - 3x4 prongs
- **26866** Weitlaner Retractor - blunt 20 cm - 3x4 prongs



- **26867** Farabeuf retractors - set of 2 - 12 cm


- **26868** Frazier Hook - sharp - 13 cm

- **26869** Frazier Hook - blunt - 13 cm

- **26723** Hook - 15 cm

- **26860** Masing Retractor - 14 cm

- **26989** Senn-Muller - 15.5 cm



- **26877** Roux Retractor set of 5 (3+2)

- **26862** Volkman- 3 teeth - sharp - 22 cm

- **26863** Volkman- 3 teeth - blunt - 22 cm



- **26861** Langenbeck Retractor - 22 cm

DISPOSABLE STAINLESS STEEL INSTRUMENTS - STERILE

A complete range of disposable instruments made of sturdy stainless steel alternative to plastic disposable instruments. Each instrument has a blue coloured part to avoid mix with reusable instruments and is individually packed in a sterile pouch.

Sterilized by Gamma radiation.
Sold in box of 10 or 25 pouches.
Instructions: GB, FR, IT, ES, DE, PT, PL, NL, SE, HU, GR, Arabic.

BLUE HANDLES
means disposable
for instant identification



STERILE SCISSORS - STAINLESS STEEL

- **38800** Surgical scissors sharp/sharp - straight - 13 cm - box of 25
- **38801** Surgical scissors sharp/blunt - straight - 13 cm - box of 25
- **38809** Metzenbaum scissors - straight - 16 cm - box of 25
- **38810** Metzenbaum scissors - curved - 16 cm - box of 25
- **38812** Metzenbaum scissors - curved - 18 cm - box of 25
- **38830** Iris scissors sharp/sharp - straight - 11.5 cm - box of 25
- **38831** Iris scissors sharp/sharp - curved - 11.5 cm - box of 25
- **38820** Mayo scissors - straight - 15 cm - box of 25
- **38822** Mayo scissors - straight - 17 cm - box of 25
- **38840** Lister bandage scissors - 14 cm
- **38841** Lister bandage scissors - 16 cm
- **38842** Lister bandage scissors - 18 cm - box of 25
- **38835** Umbilical scissors - straight - 10.5 cm - US model - box of 25

STERILE FORCEPS - STAINLESS STEEL

- **38850** Dressing tissue forceps - straight - 13 cm - box of 25
- **38852** Dressing tissue forceps - straight - 16 cm - box of 25
- **38851** Dressing forceps - straight - 13 cm 1x2 teeth - box of 25
- **38860** Pean forceps - curved - 14 cm
- **38861** Pean forceps - curved - 16 cm
- **38862** Pean forceps - curved - 18 cm - box of 25
- **38865** Spencer Wells forceps - straight - 13 cm - box of 25
- **38875** Haemostatic Halstead-Mosquito artery forceps - straight - 12.5 cm - box of 25
- **38876** Haemostatic Halstead-Mosquito artery forceps - curved - 12.5 cm - box of 25
- **38870** Adson tissue forceps - straight - 12 cm - box of 25
- **38871** Adson tissue forceps - straight - 12 cm 1x2 teeth - box of 25
- **38880** McIndoe dissecting forceps - straight - 15 cm - box of 25
- **38885** Splinter plain forceps - 12.5 cm - box of 25
- **38887** Foerster polypus forceps - 25 cm - box of 10
- **38890** Magill forceps - 20 cm for children - box of 10
- **38891** Magill forceps - 25 cm for adult - box of 10


STERILE NEEDLE HOLDERS - STAINLESS STEEL

- **38895** Micro Hartmann ear forceps - 15 cm - 3 mm - box of 10
- **38897** Micro Bellucci ear scissors - 15 cm - 3 mm - box of 10
- **38920** Mayo-Hegar needle holder - 14 cm - box of 25
- **38921** Mayo-Hegar needle holder - 16 cm - box of 25
- **38930** Crile Wood needle holder - 15 cm - box of 25

DISPOSABLE STAINLESS STEEL PROCEDURE PACKS - STERILE

- **38950 SUTURE REMOVAL PACK - box of 25**
Pack includes: Iris scissors - straight - 11.5 cm
Adson tissue forceps - toothed - 12 cm
- **38951 SUTURE PROCEDURE PACK - box of 25**
Pack includes: Spencer ligature scissors - straight - 9 cm
Adson tissue forceps - toothed - 12 cm

FOLEY CATHETERS - LATEX




GIMA code	LATEX STERILE 100% SILICONE COATED FOLEY CATHETERS - SINGLE USE	Minimum order
Two-way catheters - balloon 5-10 ml - length 40 cm - soft valve		
22054	2-way Foley catheter ch/fr 12 - balloon 5-10 ml	box of 10
22055	2-way Foley catheter ch/fr 14 - balloon 5-10 ml	box of 10
22056	2-way Foley catheter ch/fr 16 - balloon 5-10 ml	box of 10
22057	2-way Foley catheter ch/fr 18 - balloon 5-10 ml	box of 10
22058	2-way Foley catheter ch/fr 20 - balloon 5-10 ml	box of 10
22059	2-way Foley catheter ch/fr 22 - balloon 5-10 ml	box of 10

Ball 5-10 ml means that the ball could be from 5 ml to 10 ml

Two-ways catheters - balloon 30 ml - length 40 cm - soft valve		
22064	2-way Foley catheter ch/fr 12 - balloon 30 ml	box of 50
22065	2-way Foley catheter ch/fr 14 - balloon 30 ml	box of 50
22066	2-way Foley catheter ch/fr 16 - balloon 30 ml	box of 50
22067	2-way Foley catheter ch/fr 18 - balloon 30 ml	box of 50
22068	2-way Foley catheter ch/fr 20 - balloon 30 ml	box of 50
22069	2-way Foley catheter ch/fr 22 - balloon 30 ml	box of 50
22070	2-way Foley catheter ch/fr 24 - balloon 30 ml	box of 50
Three-ways catheters - balloon 30 ml - length 40 cm - soft valve		
22077	3-way Foley catheter ch/fr 18 - balloon 30 ml	box of 50
22078	3-way Foley catheter ch/fr 20 - balloon 30 ml	box of 50
22079	3-way Foley catheter ch/fr 22 - balloon 30 ml	box of 50
22080	3-way Foley catheter ch/fr 24 - balloon 30 ml	box of 50

RECTAL CATHETERS



GIMA code	PVC STERILE RECTAL CATHETERS length: 38 cm - with conical adapter	Minimum order
22169	Rectal catheters ch/fr 22 - purple	box of 100
22170	Rectal catheters ch/fr 24 - dark blue	box of 100
22171	Rectal catheters ch/fr 26 - white	box of 100
22172	Rectal catheters ch/fr 28 - green	box of 100
22173	Rectal catheters ch/fr 30 - grey	box of 100

CATHETER LUBRICANT GEL - STERILE

• **22104 CATHETER LUBRICANT GEL 12 ml - sterile - box of 25**
Sterile, clear*, water soluble gel used before applying catheter, cystoscope or other medical instrument to urethra. It covers urethra generating a slick layer between urethral mucosa and catheter or medical instrument. Antiseptic effect protecting patient from potential infection risk, it contains Lidocaine 2%, chlorhexidine 0.05%.

*After Gamma sterilization it may change the colour over time from transparent to yellowish or pinkish.



Z-STRAP 3 IN ONE - FOR FIXING URINE BAGS, TUBES OR CATHETERS

• **28727 Z-STRAP FIXING SYSTEM FOR URINE BAG, TUBE, CATHETER - box of 2**

Comfortable and discreet all in one universal strap for fixing a urine bag, a urine bag tube or a catheter. Special designed Z-strap is made of top-quality medical garment preventing irritations; silicone beading on internal side keeps the strap in place preventing any inconvenience.


- **one size fits all (70 cm long)**, can be cut to desired size
- **wide band (5 cm width)**, distributes pressure and weight better, especially for diabetics
- **easy application**, for nurse and user, fast wear and removal
- **skin friendly**, no use of adhesive / no need to shave hair
- **fixing straps for drainage tubes** to prevent slipping
- **prevents injuries** caused by sudden movement

FOLEY CATHETERS - SILICONE




GIMA code	100% SILICONE STERILE FOLEY CATHETERS - SINGLE USE	Minimum order
Two-way catheters - balloon 30 ml - length 40 cm		
22096	2-way Foley catheter ch/fr 16 - balloon 30 ml	box of 10
22097	2-way Foley catheter ch/fr 18 - balloon 30 ml	box of 10
22098	2-way Foley catheter ch/fr 20 - balloon 30 ml	box of 10
22099	2-way Foley catheter ch/fr 22 - balloon 30 ml	box of 10
22100	2-way Foley catheter ch/fr 24 - balloon 30 ml	box of 10

NELATON CATHETERS



GIMA code	PVC DEHP FREE* STERILE NELATON CATHETERS two lateral eyes with smooth edges Single use	Minimum order
MALE - length 38 cm		
22003	Male Nelaton catheter ch/fr 10 - black	box of 100
22004	Male Nelaton catheter ch/fr 12 - white	box of 100
22005	Male Nelaton catheter ch/fr 14 - dark green	box of 100
22006	Male Nelaton catheter ch/fr 16 - orange	box of 100
22007	Male Nelaton catheter ch/fr 18 - red	box of 100
FEMALE - length 18 cm		
22022	Female Nelaton catheter ch/fr 8 - light blue	box of 100
22023	Female Nelaton catheter ch/fr 10 - black	box of 100
22024	Female Nelaton catheter ch/fr 12 - white	box of 100
22025	Female Nelaton catheter ch/fr 14 - dark green	box of 100
22026	Female Nelaton catheter ch/fr 16 - orange	box of 100

SUCTION CATHETERS





Length		PVC DEHP FREE* STERILE SUCTION CATHETERS with conical adapter	Minimum order
50 cm	90 cm		
22110	22120	Suction catheter ch/fr 6 - light green	box of 100
22112	22122	Suction catheter ch/fr 8 - light blue	box of 100
22113	22123	Suction catheter ch/fr 10 - black	box of 100
22114	22124	Suction catheter ch/fr 12 - white	box of 100
22115	22125	Suction catheter ch/fr 14 - dark green	box of 100
22116	-	Suction catheter ch/fr 16 - orange	box of 100
22117	-	Suction catheter ch/fr 18 - red	box of 100
22118	-	Suction catheter ch/fr 20 - yellow	box of 100

*DEHP is an additive used in the plastic industry for PVC which could be harmful to human reproductive system. Latex is totally DEHP free

- fits all size of urine bags / leg bags

- reusable can be machine washed, does not shrink

- prevents pressure wounds, wide strap less pressure per cm²



28727

URINE BAGS - WITH/WITHOUT NON RETURN VALVE



GIMA code	STANDARD URINE BAGS	Minimum order
28690	Urine bag 2000 cc - tube 90 cm - polybag	box of 30*
28688	Urine bag 2000 cc - tube 90 cm - bulk	box of 250
28692	Urine bag 2000 cc - tube 130 cm - polybag	box of 30*
28689	Urine bag 2000 cc - tube 130 cm - bulk	box of 200

*Polybag of 30 bags with tipps (Italian reimbursement label) and 1 hanger

URINE BAGS WITH NON RETURN VALVE		
28691	Urine bag 2000 cc - tube 90 cm - polybag	box of 30*
28694	Urine bag 2000 cc - tube 90 cm - bulk	box of 250
28693	Urine bag 2000 cc - tube 130 cm - polybag	box of 30*
28695	Urine bag 2000 cc - tube 130 cm - bulk	box of 200

*Polybag of 30 bags with tipps (Italian reimbursement label) and 1 hanger

URINE LEG BAGS WITH NON RETURN VALVE		
28707	Urine leg bag 750 cc - tube 10 cm - polybag complete with push pull outlet device	box of 30**
28708	Urine leg bag 750 cc - tube 35 cm - polybag complete with push pull outlet device	box of 30**


**Polybag of 30 bags with 1 set of button straps

Instructions in GB, FR, IT, ES, DE, PT, GR, PL. Non sterile.
Tube Ø: internal 5.05 mm, external 6.35 mm. Latex and DEHP free.



GIMA code	HANGER HOLDER FOR URINE BAGS	Minimum order
28698	Hanger holder with reusable cover for urine bag	box of 10
28699	Hanger holder for urine bag - blue - bulk	box of 200

MALE EXTERNAL LATEX CATHETERS AND CATHETER FIXATIONS



External male catheters made from natural medical latex. The graded cap design eliminates the possibility of twisting of the drainage system and blockage of the urine flow. No leaking or ballooning is possible.

Multilingual polybag and instructions (GB, FR, IT, ES, PT, DE, NL, PL, CZ, RO, HU, SE, FI, GR, BG) excluding 28716

GIMA code	MALE EXTERNAL LATEX CATHETERS	Size	Minimum order
28713	Male external catheter Ø 20 mm	XS	box of 30
28709	Male external catheter Ø 25 mm	S	box of 30
28710	Male external catheter Ø 30 mm	L	box of 30
28711	Male external catheter Ø 35 mm	XL	box of 30
28712	Male external catheter Ø 40 mm	XXL	box of 30

KIT: MALE EXTERNAL CATHETER + ADHESIVE STRIP IN A POUCH		
28724	Kit: catheter Ø 30 mm + hydrocolloid strip	L box of 30
28725	Kit: catheter Ø 35 mm + hydrocolloid strip	XL box of 30

MALE EXTERNAL CATHETER FIXATIONS		
28715	Biadhesive foam strip 15x2x0.3 cm	box of 30
28716	Biadhesive hydrocolloid strip 12x1.8x0.1 cm	box of 30
28719	Safe & Secure reusable foam strap	bag of 5

STERILE PAEDIATRIC URINE BAG FOR BOY AND GIRL



• 28685 PAEDIATRIC URINE BAG 100 ml - sterile - box of 100

Paediatric urine bag, individually wrapped in polybag. They can be used for boy and girl. Graduated 10 to 100 ml, sterile.



STERILE MEDICATION, SUTURE, CATHETERIZATION KITS




• 26930 DRESSING KIT 1 - sterile



• 26929 DRESSING KIT 2 - sterile



• 26933 CATHETERIZATION KIT - sterile



• 26936 SUTURE REMOVAL KIT 1 - sterile



• 26945 SUTURE REMOVAL KIT 2 - sterile



• 26932 SUTURE REMOVAL KIT 3 - sterile

KIT CONTENT	Dressing		Suture		Suture removal			Catheters
	26930	26929	26931	26949	26932	26936	26945	26933
Field cloth	•		•		•			•
Cushing surgery forceps			•					
Needle Holder			•					
Kocher forceps			•	•				
Anatomy forceps	•		•	•	•			
Plastic forceps			•	•		•	•	
Haemostatic flow clamp								
Ring shaped forceps	•	•	•		•			•
Styptic clamp (Kelly forceps)			•					
Haemostatic clamp			•	•				
Straight scissors			•					
Adson anatomy forceps 12 cm		•		•				
Light green plastic scissors		•		•		•		
Littauer scissors					•			
Gauze swabs 7.5x7.5 cm 8 ply	•	•	•		•			
Gauze swabs 10x10 cm 8 ply			•	•	•			
4 or 5 Gauze plugs	•		•	•	•			
3 Gauze plugs		•						•
Gauze swabs 10x20 cm - 16 plies		•						
Gauze swab 18x40 cm						•		
Latex gloves - medium								
Vinyl gloves			•		•			•
Disposable big bag			•		•			•
Medication plaster 15x12 cm					•			
WNW (TNT)+PE drape 32x50 cm		•						
WNW (TNT)+PE drape 50x60 cm				•				
WNW (TNT) compresses								•
WNW (TNT)+PE drape 50x60 cm with hole			•	•				•
Water syringe								•
Single dose lubricant								•
2 disinfectant wipes							•	•
2 adhesive WNW (Woven Non Woven) strips				•				
Waste bag		•						
Scalpel with handle							•	
2 cotton applicators							•	
PE bag minigrip				•				



• 26931 SUTURE KIT 1 - sterile



• 26949 SUTURE KIT 2 - sterile

Labellinstructions in GB, IT for 26930-1-2-3
Labellinstructions in GB, FR, IT, ES for 26929, 26936,
26945, 26949

GENERAL SURGERY PROCEDURE PACK



26985

• 26985 GENERAL SURGERY PACK - sterile

Ready to use pack including all necessary for a general surgery procedure.

Peel open PE + tyvek envelope.

Sterile EO.

KIT CONTENT

1 table cover 160x200 cm	4 hand towels
2 surgical gowns XL (1 standard, 1 reinforced)	1 bowl 120 ml
2 adhesive drapes 75x90 cm	1 kidney dish
1 foot drape (ADH side) 150x180 cm	1 scalpel n°15
1 head drape (ADH side) 150x240 cm	5 gauze swabs 10x20 cm
2 syringes 10 ml	4 gauze balls 30 mm
2 needles 21Gx1.5"	2 foam sticks 15 cm
1 scalpel n°11	1 basic tray

SURGICAL SKIN MARKER - STERILE

• 33181 P3 SKIN MARKER® - fine tip - sterile - box of 10

The P3® Surgical Skin Marker features a versatile tip for precise fine line marking preferred in plastic surgery or when held at an angle, broad line marking for applications in general, neuro and vascular surgery, as well as radiology and radiation therapy. The marker's high visibility Gentian violet ink is non-toxic, non tattooing and surgical prep-resistant. Markings can be removed from the skin with isopropyl (rubbing) alcohol. The package is easy to open, HIBC bar coded and includes a separate six inch, flexible ruler calibrated in both inches and centimeters. Sterile EO. Individually packed.

CE 0169



TWO TIPS IN ONE:




For heavier lines, hold the pen at an angle and use the shoulder of the tips (2 mm)



Hold the pen vertically while applying a fine line (0.2 mm)

MADE IN THE UK

SURGICAL SKIN MARKERS - STERILE



Multilingual box and instructions: GB, FR, IT, ES, PT, DE, NL, DK, PL, CZ, FI, GR.

• 33182 SURGICAL SKIN MARKER - single tip - sterile - box of 100

Single tip marker ideal for plastic surgery. Hypoallergenic and violet non-toxic ink. Marking can be removed by alcohol and liquid soap. Intended for single use. Includes a separate six inch, flexible ruler calibrated in both inches and centimeters. Individually packed.

TECHNICAL SPECIFICATIONS		
Tip size:	1.0 mm	Ink colour: violet
Diameter:	9.5 mm	Washable: no
Length:	167 mm	



• 33183 SURGICAL SKIN MARKER - double tip - sterile - box of 100

Double tip marker ideal for plastic surgery. Hypoallergenic and violet non-toxic ink. Marking can be removed by alcohol and liquid soap. Intended for single use. Includes a separate six inch, flexible ruler calibrated in both inches and centimeters. Individually packed.

TECHNICAL SPECIFICATIONS		
Tip size:	fine tip: 0.5mm	Length: 140 mm
	thick tip: 1.0mm	Ink colour: violet
Diameter:	11 mm	Washable: no

SKIN MARKER

• 33179 SKIN MARKER - dual tips gentian violet - box of 10

Double tip skin marker with polypropylene cap and body. Polyester fiber tips.

Ultra fine tip utility marker to write on adhesive labels and infusion bags.

Multilingual box and manual: GB, FR, IT, DE, ES, PT, GR, PL, SE, RO, HU, Arabic.



Regular tip surgical skin marker for use on patient skin. Washable with water and soap. Size 138xØ 12.5 mm.



DERMATOGRAPH MARKERS

- 33175 DERMATOGRAPH MARKER - blue - box of 6
- 33176 DERMATOGRAPH MARKER - white - box of 6
- 33177 DERMATOGRAPH MARKER - red - box of 6
- 33178 DERMATOGRAPH MARKER - mixed colours - box of 6*

*2 pencils each colour

Special dermatograph pencil to mark the skin for identification of areas to be treated.

Made of cedar wood with wax lead. Size 175xØ 7 mm. Suitable for general surgery, dermatology, orthopaedy, radiology, ...

Pencils are sold non sharpened in box of 6 with 1 pencil sharpener included.

Multilingual box and manual:

GB, FR, IT, DE, ES, PT, GR, PL, SE, RO, HU, Arabic.



33175

33176



33177


33178



SUTURE TRAINING PAD

• 23020 SUTURE TRAINING PAD with wounds with mesh

Silicone suture training pad with 14 wounds. 3 layers: skin, subcutis, muscle (extra layer mesh under the skin). Soft true to life skin texture. Multilingual label and manual: GB, FR, IT, ES, PT.



23020

Size 17.2x12x1.2 cm

3M STERI-STRIP™

Adhesive skin closure strips, sterile, breathable and comfortable to wear. See page 410 for: 3M reinforced or elastic Steri-strip, suture strips



SUTURE TRAINING KIT

• 23022 SUTURE TRAINING KIT

Suture training kit includes: silicone pad (23020), 6 tools and 4 blades made of stainless steel, 8 sutures and a black PVC leather case.

Intended for teaching and training only, not for human or vet use.

Multilingual label and manual:

GB, FR, IT, ES, PT, DE, SE, DK, GR, RO.

23022



Kit size 19.5x22.5x3.5 cm

STANDARD ACCESSORIES

1 Case

1 Suture pad

6 Instruments:

1 Needle holder - 12.5 cm

1 Adson forcep - 12.5 cm

1 Iris scissor - 12.5 cm

1 Scalpel handle n.3 - 12.5 cm

1 Mosquito forcep - 12.5 cm

1 Dressing forcep - 12.5 cm

4 Blades:

2 Scalpel blades n.10

2 Scalpel blades n.11

8 Sutures:

2 Nylon sutures 3/0

2 Nylon sutures 4/0

2 Silk sutures 3/0

2 Silk sutures 4/0

SUTURES COMPARISON HEALTHIUM - ETHICON - B.BRAUN

GENERIC NAME OF SUTURES	Healthium	ETHICON	B. BRAUN SHARING EXPERTISE
CATGUT PLAIN & CHROMIC	-	GUT SUTURE PLAIN & CHROMIC	-
POLYGLYCOLIC ACID	TRUGLYDE	-	SAFIL
FAST ABSORBING POLYGLYCOLIC ACID	TRUGLYDE FAST	-	SAFIL QUICK
POLYGLACTIN 910	TRUSYNTH	VICRYL	NOVONSYN
FAST ABSORBING POLYGLACTIN 910	TRUSYNTH FAST	VICRYL RAPID	-
ANTIBACTERIAL COATED POLYGLACTIN 910	-	VICRYL PLUS	-
MONOFILAMENT POLYDIOXANONE	PD SYNTH	PDS OR PDS II	MONOPLUS
MONOFILAMENT POLIGLECAPRONE 25	MONOGLYDE	MONOCRYL	MONOSYN
MONOFILAMENT POLYAMIDE NYLON	TRULON	ETHILON / NUROLON	DAFILON
BRAIDED & COATED SILK	TRUSILK	MERSILK / PERMA-HAND SILK	SILKAM
MONOFILAMENT POLYPROPYLENE	TRULENE	PROLENE	PREMILENE/OPTILENE
BRAIDED & COATED POLYESTER	TRUBOND	ETHIBOND/MERSILENE	PREMICRON/DAGROFIL
STEEL SUTURES	TRUSTEEL	ETHISTEEL	STEELEX
BARB SUTURES	TRUBARB	STRATAFIX	-

SURGICAL NON ABSORBABLE SUTURES - STERILE

MDR Certified

TRUBOND - non absorbable surgical suture, silicone coated and braided polyester for permanent support and strength

Equivalent Ethicon brand: ETHIBOND/MERSILENE


Technical details:

- braided polyester, coated with silicone for a smooth knot tie down
- high tensile strength for reliability
- low tissue reactivity - undergoes no significant change in body
- excellent handling and long lasting
- sterilized by ethylene oxide
- shelf life of 5 years

Applications:

- used for permanent and removable fixtures
- prosthetic implantations
- "Face lifts"
- cardiovascular surgery
- ideal for valve repair/ replacement

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUSTEEL - non absorbable monofilament steel surgical suture
Equivalent Ethicon brand: ETHISTEEL

Technical details:

- monofilament stainless steel LVM 316 grade wire
- optimal compatibility with body tissue
- low or no tissue reactivity
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

Applications:

- orthopedic surgeries
- sternum closure

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUBOND - box of 12

GIMA code	Gauge	Metric	Colour	Cirde	Needle profile	Needle length	Suture length
22770	2/0	3	Green	1/2	Taper cutting double arm	17 mm	90 cm
22771	3/0	2	Green	1/2	Round body	25 mm	90 cm
22772	3/0	2	Green	3/8	Cutting	16 mm	38 cm
22773	4/0	1.5	Green	3/8	Round body	16 mm	76 cm
22774	5/0	1	Green	3/8	Round body double arm	20 mm	76 cm
22775	0	3.5	Green	1/2	Reverse cutting	50 mm	100 cm
22776	1	4	Green	1/2	Reverse cutting	50 mm	100 cm
22777	2	5	Green	1/2	Taper cutting heavy	45 mm	100 cm
22778	5	7	Green	1/2	Trocar point heavy	55 mm	100 cm

TRUSTEEL - box of 48

GIMA code	Gauge	Metric	Colour	Cirde	Needle profile	Needle length	Suture length
22790	1	4	Natural	1/2	Round body blunt point	40 mm	45 cm
22791	2	5	Natural	1/2	Taper cutting	40 mm	45 cm
22792	4	6	Natural	1/2	Conventional cutting	48 mm	45 cm
22793	5	7	Natural	1/2	Taper cutting	48 mm	45 cm
22794	6	8	Natural	1/2	Round body blunt point	44 mm	45 cm
22795*	6	8	Natural	1/2	Conventional cutting	48 mm	75 cm
22796	7	9	Natural	1/2	Conventional cutting	48 mm	45 cm

*This code only is sold in box of 24

SURGICAL NON ABSORBABLE SUTURES - STERILE MDR Certified

TRULON - polyamide nylon non absorbable monofilament suture

Equivalent Ethicon brand: ETHILON / NUROLON

Technical details

- monofilament polyamide
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

Applications

- for cuticular and subcuticular closure
- orthopedic surgery
- ophthalmic
- general closure
- micro surgery

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRULON - box of 12							
GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22800	2/0	3.0	Blue	3/8	Reverse Cutting	26 mm	45 cm
22801	2/0	3.0	Black	3/8	Reverse Cutting	45 mm	70 cm
22803	3/0	2.0	Blue	3/8	Conventional Cutting	20 mm	45 cm
22804	3/0	2.0	Blue	3/8	Reverse Cutting	24 mm	45 cm
22805	3/0	2.0	Blue	3/8	Conventional Cutting	25 mm	45 cm
22807	3/0	2.0	Black	3/8	Reverse Cutting	19 mm	45 cm
22808	3/0	2.0	Black	3/8	Reverse Cutting	24 mm	45 cm
22810	4/0	1.5	Blue	3/8	Reverse Cutting	19 mm	45 cm
22811	4/0	1.5	Black	3/8	Conventional Cutting	10 mm	38 cm
22812	4/0	1.5	Black	3/8	Reverse Cutting	19 mm	45 cm
22814	5/0	1.0	Blue	3/8	Reverse Cutting	16 mm	70 cm
22815	5/0	1.0	Blue	3/8	Reverse Cutting	19 mm	70 cm
22816	5/0	1.0	Black	3/8	Reverse Cutting	19 mm	45 cm
22817	5/0	1.0	Black	3/8	Reverse Cutting	19 mm	75 cm
22820	6/0	0.7	Blue	3/8	Reverse Cutting	12 mm	45 cm
22821	6/0	0.7	Blue	3/8	Reverse Cutting	16 mm	45 cm
22822	6/0	0.7	Black	3/8	Reverse Cutting	16 mm	45 cm

TRULENE - monofilament non absorbable suture - polypropylene

Equivalent Ethicon brand: PROLENE

Technical details

- monofilament polypropylene
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

Applications

- hernia repair
- permanent fixtures
- general surgery
- micro surgery
- orthopedic surgery
- cardiovascular surgery
- plastic surgery
- can also be used in ophthalmic and neurological procedures

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUSILK - braided black silk non absorbable multifilament suture

Equivalent Ethicon brand: MERSILK / PERMA-HAND SILK

Technical details

- made from 100% protein fiber spun
- coating of wax/silicone
- sterilized by ethylene oxide
- shelf life of 5 years

Applications

- general surgery
- ophthalmic surgery
- cardiovascular surgery
- gastrointestinal surgery
- orthopedic surgery

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUSILK - box of 12							
GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22830	0	3.5	Black	1/2	Round Body	30 mm	76 cm
22831	2/0	2.0	Black	1/2	Taper Cutting	19 mm	75 cm
22832	2/0	3.0	Black	3/8	Reverse Cutting	26 mm	45 cm
22833	2/0	3.0	Black	1/2	Reverse Cutting	22 mm	45 cm
22834	2/0	3.0	Black	1/2	Taper Cutting	25 mm	75 cm
22835	2/0	3.0	Black	1/2	Round Body	30 mm	76 cm
22836	3/0	2.0	Black	1/2	Conventional Cutting	22 mm	45 cm
22837	3/0	2.0	Black	1/2	Round Body	22 mm	45 cm
22838	3/0	2.0	Black	3/8	Reverse Cutting	19 mm	70 cm
22839	3/0	2.0	Black	3/8	Reverse Cutting	24 mm	45 cm
22840	3/0	2.0	Black	1/2	Taper Cutting	19 mm	75 cm
22842	3/0	2.0	Black	1/2	Taper Cutting	25 mm	75 cm
22844	4/0	1.5	Black	1/2	Reverse Cutting	17 mm	45 cm
22845	4/0	1.5	Black	1/2	Taper Cutting	19 mm	75 cm
22846	4/0	1.5	Black	1/2	Taper Cutting	25 mm	75 cm
22847	4/0	1.5	Black	Straight	Round Body	19 mm	45 cm
22848	4/0	1.5	Black	3/8	Reverse Cutting	13 mm	45 cm
22849	4/0	1.5	Black	3/8	Reverse Cutting	19 mm	45 cm
22851	5/0	1.0	Black	3/8	Reverse Cutting	19 mm	70 cm
22852	6/0	0.7	Black	3/8	Reverse Cutting	16 mm	45 cm

TRULENE - box of 12							
GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22870	1	4.0	Blue	1/2	Round Body	30 mm	70 cm
22871	0	3.5	Blue	1/2	Round Body	30 mm	70 cm
22872	2/0	3.0	Blue	3/8	Reverse Cutting	24 mm	70 cm
22874	3/0	2.0	Blue	3/8	Reverse Cutting	24 mm	70 cm
22875	3/0	2.0	Blue	Straight	Conventional Cutting	60 mm	70 cm
22877	4/0	1.5	Blue	3/8	Reverse Cutting	19 mm	45 cm
22879	4/0	1.5	Blue	3/8	Reverse Cutting	19 mm	75 cm
22881	5/0	1.0	Blue	3/8	Reverse Cutting	16 mm	70 cm
22882	5/0	1.0	Blue	3/8	Reverse Cutting	19 mm	70 cm
22884	6/0	0.7	Blue	3/8	Reverse Cutting	13 mm	45 cm
22885	6/0	0.7	Blue	3/8	Conventional Cutting	16 mm	45 cm

SURGICAL ABSORBABLE SUTURES - STERILE

TRUBARB - absorbable surgical needled-suture (synthetic) knotless tissue-closure device

Equivalent Ethicon brand: STRATAFIX

Technical details

- monofilament polydioxanone
- secure and effective triangular end-stopper
- single angle cut uni-directional helical barb design
- world class 300 series QNPL needle
- sterilized by ethylene oxide
- shelf life of 3 years

Applications

- gynaecological and urology procedures
- general surgery procedures
- oncology procedures
- plastic and cosmetic surgeries

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUBARB - box of 6							
GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22890	0	3.5	Violet	1/2	Taper point	37 mm	30 cm
22891	0	3.5	Violet	1/2	Taper point	37 mm	45 cm
22892	2/0	3	Violet	1/2	Taper point	37 mm	30 cm
22893	2/0	3	Violet	1/2	Taper point	37 mm	45 cm
22894	2/0	3	Violet	1/2	Taper point	26 mm	30 cm
22895	3/0	2	Violet	1/2	Taper point	26 mm	45 cm
22896	3/0	2	Violet	3/8	Reverse cutting	19 mm	45 cm

SURGICAL ABSORBABLE SUTURES - STERILE

MDR Certified

MONOGLYDE - monofilament absorbable suture poliglecaprone 25

Equivalent Ethicon brand: MONOCRYL

Technical details

- monofilament
- poliglecaprone 25
- synthetic absorbable suture
- sterilized by ethylene oxide
- shelf life of 3 years
- wound support: 21-28 days
- mass absorption: 90-120 days

Applications

- used for subcuticular repair and general soft-tissue approximation
- bowel anastomoses
- closure of the peritoneum
- closure of uterus and vaginal cuff
- not recommended for use in cardiovascular or neurological tissues, microsurgery or ophthalmic surgery
- since this is an absorbable suture, it should not be used for extended approximation of tissue (fascia)

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



MONOGLYDE - box of 12

GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22901	0	3.5	Undyed	1/2	Round body	30 mm	70 cm
22902	2/0	3.0	Undyed	1/2	Round body	30 mm	70 cm
22903	3/0	2.0	Undyed	3/8	Slim blade cutting	25 mm	70 cm
22904	5/0	1.0	Undyed	3/8	Round body	19 mm	45 cm

TRUSYNTH FAST - braided & coated polyglactin 910 fast

Equivalent Ethicon brand: VICRYL RAPID

Technical details

- synthetic absorbable surgical suture
- braided and coated polyglactin 910

- sterilized by ethylene oxide
- shelf life of 5 years
- absorption period: 42 days

Applications

- general soft tissue approximation

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUSYNTH FAST - box of 12

GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22910	2/0	3.0	Undyed	3/8	Reverse cutting	24 mm	70 cm
22911	3/0	2.0	Undyed	1/2	Round body	17 mm	45 cm
22912	3/0	2.0	Undyed	3/8	Reverse cutting	19 mm	45 cm
22913	3/0	2.0	Undyed	1/2	Round body	17 mm	70 cm
22914	3/0	2.0	Undyed	1/2	Round body	22 mm	70 cm
22915	3/0	2.0	Undyed	3/8	Reverse cutting	14 mm	75 cm
22916	3/0	2.0	Undyed	3/8	Reverse cutting	24 mm	70 cm
22917	4/0	1.5	Undyed	3/8	Reverse cutting	19 mm	45 cm
22918	4/0	1.5	Undyed	3/8	Reverse cutting	19 mm	70 cm
22919	5/0	1.0	Undyed	3/8	Reverse cutting	16 mm	70 cm
22920	6/0	0.7	Undyed	3/8	Reverse cutting	11 mm	45 cm

TRUGLYDE FAST - braided and coated polyglycolic acid fast undyed

Technical details

- braided and coated
- coated with calcium stearate and polycaprolactone
- sterilized by ethylene oxide
- shelf life of 5 years
- wound support: 7-10 days
- mass absorption: 40-42 days

Applications

- for episiotomy
- for circumcision
- closure of skin

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUGLYDE FAST - box of 12

GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22980	1	4.0	Undyed	1/2	Round body	40 mm	90 cm
22981	2/0	3.0	Undyed	1/2	Round body	35 mm	75 cm
22983	3/0	2.0	Undyed	3/8	Reverse cutting	24 mm	70 cm
22984	3/0	2.0	Undyed	3/8	Reverse cutting	20 mm	75 cm
22986	4/0	1.5	Undyed	3/8	Reverse cutting	16 mm	70 cm
22988	4/0	1.5	Undyed	3/8	Reverse cutting	20 mm	75 cm
22989	5/0	1.0	Undyed	1/2	Round body	17 mm	75 cm

PDSYNTH - monofilament absorbable suture polydioxanone

Equivalent Ethicon brand: PDS or PDS II

Technical details

- monofilament
- polydioxanone
- synthetic absorbable suture
- sterilized by ethylene oxide
- shelf life of 5 years
- wound support: 35-42 days
- mass absorption: 180-210 days

Applications

- for rectus sheath closure
- fascia closure and gastrointestinal tract
- pediatric cardiovascular surgery
- geriatric surgeries
- orthopedic surgery
- general surgery

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



PDSYNTH - box of 12

GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22940	1	4.0	Violet	1/2	Round body	40 mm	90 cm
22941	0	3.5	Violet	1/2	Round body	30 mm	70 cm
22942	0	3.5	Violet	1/2	Round body	40 mm	70 cm

TRUSYNTH - braided & coated polyglactin 910

Equivalent Ethicon brand: VICRYL

Technical details

- braided and coated polyglactin 910
- coating of calcium stearate and polyglactin 910
- sterilized by ethylene oxide
- shelf life of 5 years
- wound support: 28-30 days
- mass absorption: 56-70 days

Applications

- ophthalmic procedures
- mucosal closures
- skin closure
- dental surgery
- orthopedic surgery
- general soft tissue approximation

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUSYNTH - box of 12

GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22950	1	4.0	Violet	1/2	Round body	40 mm	90 cm
22951	0	3.5	Violet	1/2	Round body	40 mm	90 cm
22952	2/0	3.0	Violet	1/2	Taper cutting	19 mm	75 cm
22953	2/0	3.0	Violet	1/2	Round body	30 mm	90 cm
22954	3/0	2.0	Violet	1/2	Taper cutting	19 mm	75 cm
22955	4/0	1.5	Violet	1/2	Taper cutting	19 mm	75 cm

TRUGLYDE - braided & coated polyglycolic acid

Technical details

- braided and coated polyglycolic acid
- coating of calcium stearate and polyglactin 310
- sterilized by ethylene oxide
- shelf life of 5 years
- wound support: 28-30 days
- mass absorption: 70-90 days

Applications

- for soft tissue approximation during general, gynecological and orthopedic surgeries

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUGLYDE - box of 12

GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length
22960	1	4.0	Violet	1/2	Round body	40 mm	90 cm
22961	0	3.5	Violet	1/2	Round body	40 mm	90 cm
22962	2/0	3.0	Violet	3/8	Reverse cutting	24 mm	70 cm
22963	2/0	3.0	Violet	1/2	Round body	30 mm	90 cm
22964	3/0	2.0	Violet	3/8	Reverse cutting	19 mm	45 cm
22965	3/0	2.0	Violet	3/8	Reverse cutting	19 mm	70 cm
22966	3/0	2.0	Undyed	3/8	Reverse cutting	19 mm	75 cm
22967	3/0	2.0	Undyed	1/2	Round body	22 mm	75 cm
22968	4/0	1.5	Violet	3/8	Reverse cutting	19 mm	45 cm
22969	4/0	1.5	Violet	3/8	Reverse cutting	19 mm	70 cm
22970	4/0	1.5	Undyed	3/8	Reverse cutting	16 mm	45 cm
22971	4/0	1.5	Undyed	1/2	Round body	17 mm	75 cm
22972	5/0	1.0	Undyed	1/2	Round body	17 mm	75 cm
22973	5/0	1.0	Undyed	3/8	Reverse cutting	13 mm	45 cm

ETHICON SUTURES - ABSORBABLE AND NON ABSORBABLE




NON ABSORBABLE STERILE SURGICAL SUTURES

GIMA **ETHICON** **Gauge** **Metric** **Colour** **Circle** **Needle length** **Needle** **Suture length** **Minimum order**

Non absorbable surgical sutures are indicated for soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures

ETHILON - monofilament suture



22296	1669H	3/0	2.0	Black	3/8	19 mm	FS-2	45 cm	box of 36
22297	EH7145H	4/0	1.5	Black	3/8	19 mm	FS-2	45 cm	box of 36
22298	1162H	4/0	1.5	Black	3/8	16 mm	FS-3	45 cm	box of 36
22299	2808G	8/0	0.4	Black	3/8	6.5 mm	BV130-5	13 cm	box of 12
22300	660H	6/0	0.7	Black	3/8	16 mm	FS-3	45 cm	box of 36
22301	661H	5/0	1.0	Black	3/8	19 mm	FS-2	45 cm	box of 36
22302	662SLH	4/0	1.5	Black	3/8	19 mm	FS-2S	45 cm	box of 36
22303	663H	3/0	2.0	Black	3/8	24 mm	FS-1	45 cm	box of 36
22304	EH7506H	3/0	2.0	Black	3/8	19 mm	FS-2	45 cm	box of 36
22305	FH5171H	5/0	1.0	Black	3/8	19 mm	FS-2	75 cm	box of 36
22306	668H	5/0	1.0	Black	3/8	13 mm	C-2	45 cm	box of 36
22307	698H	5/0	1.0	Black	3/8	13 mm	P-3	45 cm	box of 36
22308	669H	3/0	2.0	Black	3/8	24 mm	FS-1	75 cm	box of 36

PERMA-HAND - silk sutures - braided



22310	18501G	4/0	1.5	Black	1/2	17 mm	J-1	45 cm	box of 12
22311	632G	3/0	2.0	Black	1/2	22 mm	X-1	45 cm	box of 12
22312	737G	2/0	3.0	Black	1/2	22 mm	X-1	45 cm	box of 12
22313	18507G	4/0	1.5	Black	Straight	19 mm	ST-4	45 cm	box of 12
22314	641G	4/0	1.5	Black	3/8	13 mm	P-3	45 cm	box of 12
22315	18509G	3/0	2.0	Black	1/2	22 mm	SH-1 Plus	45 cm	box of 12
22316	EH7149G	4/0	1.5	Black	3/8	19 mm	FS-2	45 cm	box of 12
22317	K872	3/0	2.0	Black	1/2	17 mm	RB-1	75 cm	box of 36
22318	683H	4/0	1.5	Black	3/8	19 mm	FS-2S	45 cm	box of 36
22319	684H	3/0	2.0	Black	3/8	24 mm	FS-1	45 cm	box of 36

PROLENE - blue pigmented monofilament sutures



22320	EH7150H	4/0	1.5	Blue	3/8	19 mm	FS-2	75 cm	box of 36
22321	EH7151H	4/0	1.5	Blue	3/8	19 mm	FS-2	45 cm	box of 36
22322	W8873T	6/0	0.7	Blue	3/8	16 mm	PC3Prime	45 cm	box of 24
22323	8661H	5/0	1.0	Blue	3/8	19 mm	FS-2	45 cm	box of 36
22324	8665H	3/0	2.0	Blue	3/8	19 mm	FS-3	45 cm	box of 36
22325	8682H	4/0	1.5	Blue	3/8	19 mm	PS-2	45 cm	box of 36
22326	8683H	4/0	1.5	Blue	3/8	19 mm	FS-2	45 cm	box of 36
22327	8687H	3/0	2.0	Blue	3/8	19 mm	PS-2	45 cm	box of 36
22328	8698H	5/0	1.0	Blue	3/8	13 mm	P-3	45 cm	box of 36
22329	8699H	4/0	1.5	Blue	3/8	13 mm	P-3	45 cm	box of 36

ABSORBABLE STERILE SURGICAL SUTURES

GIMA **ETHICON** **Gauge** **Metric** **Colour** **Circle** **Needle length** **Needle** **Suture length** **Minimum order**

MONOCRYL - complete absorption in 91-119 days



Monocryl is indicated for general soft tissue approximation and/or ligation, but not for cardio, neuro, micro and ophthalmic surgery

22350	Y493H	5/0	1.0	Undyed	3/8	13 mm	P-3 prime	45 cm	box of 36
22351	MCP3205G	4/0	1.5	Undyed	3/8	16 mm	PS-3	70 cm	box of 12

VICRYL PLUS BRAIDED - complete absorption in 56-70 days



Vicryl Plus is indicated for use in general tissue approximation and/or ligation except for ophthalmic, neuro, cardiovascular

22360	VCP392ZH	4/0	1.5	Violet	3/8	19 mm	FS-2	45 cm	box of 36
22361	VCP305H	3/0	2.0	Violet	1/2	17 mm	RB-1 Plus	70 cm	box of 36
22362	VCP452H	3/0	2.0	Violet	3/8	24 mm	FS-1	70 cm	box of 36
22363	VCP453H	2/0	3.0	Violet	3/8	24 mm	FS-1	70 cm	box of 36
22364	VCP311H	3/0	2.0	Violet	3/8	22 mm	SH-1	70 cm	box of 36
22365	VCP317H	2/0	3.0	Violet	1/2	26 mm	SH	70 cm	box of 36
22366	VCP494H	4/0	1.5	Undyed	3/8	13 mm	P-3	45 cm	box of 36

VICRYL RAPID BRAIDED - complete absorption in 42 days



Vicryl Rapid Braided indicated only for superficial soft tissue approximation of the skin and mucosa, where only short term (7-10 days) wound support is required.


22370	V76H	3/0	2.0	Undyed	1/2	17 mm	J-1	45 cm	box of 36
22371	V2920G	4/0	1.5	Undyed	3/8	19 mm	FS-2	45 cm	box of 12
22372	V2930G	3/0	2.0	Undyed	3/8	19 mm	FS-2	45 cm	box of 12
22373	V2150G	3/0	2.0	Undyed	1/2	17 mm	RB-1 Plus	70 cm	box of 12
22374	V2190G	3/0	2.0	Undyed	1/2	22 mm	SH-1 Plus	70 cm	box of 12
22375	VR2298	3/0	2.0	Undyed	3/8	19 mm	FS-2	75 cm	box of 36
22376	V2910H	5/0	1.0	Undyed	3/8C	19 mm	FS-2	45 cm	box of 36

VICRYL



22380	V453H	2/0	3.0	Violet	3/8	24 mm	FS-1	70 cm	box of 36
22381	V452H	3/0	2.0	Violet	3/8	24 mm	FS-1	70 cm	box of 36
22382	V304H	4/0	1.5	Violet	1/2	17 mm	RB-1	70 cm	box of 36
22383	V305H	3/0	2.0	Violet	1/2	17 mm	RB-1	70 cm	box of 36
22384	V311H	3/0	2.0	Violet	1/2	22 mm	SH-1	70 cm	box of 36
22385	18531G	4/0	1.5	Violet	Compaure	25 mm	Davis black	45 cm	box of 12
22386	V994H	4/0	1.5	Violet	1/2C	17 mm	V5	70 cm	box of 36
22387	V422H	3/0	2.0	Undyed	3/8C	24 mm	FS-1	70 cm	box of 36

PDS II




22390	Z340H	0	3.5	Violet	1/2	36 mm	CT-1	70 cm	box of 36
22391	Z317H	2/0	3.0	Violet	1/2	26 mm	SH	70 cm	box of 36
22392	Z292ZE	4/0	1.5	Undyed	3/8	19 mm	FS-2	45 cm	box of 24
22393	Z494E	4/0	1.5	Undyed	3/8	13 mm	P-3	45 cm	box of 24
22394	Z496E	4/0	1.5	Undyed	3/8	19 mm	PS-2	45 cm	box of 24

DISPOSABLE SKIN STAPLER



Multilingual manual and box:
GB, FR, IT, ES, PT, DE.



- **25891 DISPOSABLE STERILE SKIN STAPLER** - box of 5
Easy to use device with 35 staples made of High Quality 317 L Stainless Steel for medical use.
High sewing up and stapling speed help shorten the time of sewing and avoid cross infection of contagious diseases. Reduced tissue reaction, thinner scar after operation and faster wound healing. Painless and convenient staples removal after the incision has healed up.
For skin closure in operating room, delivery, emergency, outpatient surgery, clinics and physician's office.

3M™ PRECISE™ SKIN STAPLERS


- **25896 3M PRECISE™ SKIN STAPLER** 5 staples - sterile - box of 12
- **25897 3M PRECISE™ SKIN STAPLER** 15 staples - sterile - box of 12

The disposable Precise skin staplers offer numerous benefits over sutures


- significantly quicker, simpler to use and safer than sutures
- better cosmetic effect
- six times faster than suturing
- no risk of accidental stick injury

- less painful than sutures
- easy to remove.

It comes pre-loaded with 5 or 15 staples (6.5x4.5 mm). Suitable for simple, recent, open, unbruised wounds. Staples are particularly suitable for scalp wounds as they do not catch on surrounding hair and are more easily identifiable at the time of removal.



SKIN STAPLE REMOVERS



25893

- **25893 DISPOSABLE SKIN STAPLE REMOVER** 11 cm - sterile - box of 50
Suitable for use with most brands of skin staples. It combines functionality with cost effectiveness.

26722

26737

- **26737 S/S SPENCER SCISSORS** - 13 cm
- **26722 S/S MICHEL STITCH REMOVER FORCEPS** - 12 cm



- **25892 DISPOSABLE REMOVER FORCEPS** - sterile - box of 20 pcs.
Practical disposable forceps for removing regular and wider size staples. Plastic handle and stainless steel tips.



CE 0086



25895

- **25895 SR-3 DISPOSABLE SKIN STAPLE REMOVER** - sterile - box of 10 pcs.
The scissor-like action of the 3M SR-3 provides ergonomic grip and removes staples in the same direction they were implanted, making removal simple. Plastic handle and stainless steel tips.

BONEWAX - STERILE

- **23000 TRUWAX STERILE SURGICAL BONEWAX** 2.5 g - box of 12

Equivalent Ethicon brand: **BONE WAX**

Truwax is a sterile bone wax used to control bleeding from bone surfaces during surgical procedures by providing hemostasis from bone injuries and cuts. It works like a mechanical barrier, through a tamponade effect, it blocks pores and channels of cut or damaged bone. Made from sterile mixture of beeswax paraffin wax and isopropyl palmitate. Materials ensure excellent malleability. Multilingual manual and box:

GB, FR, IT, ES, PT, DE, GR.

MDR Certified

TECHNICAL SPECIFICATIONS

Composition:	Bees wax I.P., White Hard Paraffin Wax, I.P., Iso Propyl Palmitate, U.S.P.
Colour - Odour:	opaque - waxy odour
Sterilization:	gamma radiation
Shelf life:	5 years
Packing:	2.5 g per unit pack



23000

CLINICEL OXIDIZED REGENERATED CELLULOSE

CLINICEL OXIDIZED REGENERATED CELLULOSE

Equivalent Ethicon brand: **SURGICEL**

Features

- absorbable hemostat
- made from oxidized regenerated cellulose (plant based product)
- malleable and flexible

Applications

- Neurosurgery: operative intervention into the cerebrum, ventricular system, spinal surgery
- Cardiac Surgery: valve replacement surgery, sterna medullary bleeding, approaches to the skull base, cerebrovascular surgery
- Gynecological surgeries: uterine wall perforation
- ENT: posterior epistaxis
- General surgeries: iatrogenic perforations

TECHNICAL SPECIFICATIONS

Effective Hemostasis:	in 3-5 minutes
Sterilization:	gamma radiation
Shelf life:	3 years



23005

23008

Multilingual manual and box:
GB, FR, IT, ES, PT, DE, GR.

MDR Certified

GIMA code	CLINICEL STANDARD	box of 6	cm	inches
23003	Clinicel standard	5.1x7.6	2x3	
23004	Clinicel standard	5.1x35.5	2x14	
23005	Clinicel standard	10x20	4x8	
GIMA code	CLINICEL FIBRIL TYPE	box of 6	cm	inches
23007	Clinicel fibril type	2.5x5.1	1x2	
23008	Clinicel fibril type	5.1x10	2x4	

SURGICAL PREP USE RAZORS - STAINLESS STEEL BLADE



27081



27082

GIMA code	NON STERILE DISPOSABLE RAZORS stainless steel blades	Blade length	Minimum order
27080	Single blade - 98.6x41 mm	35 mm	box of 100
27081	Double blade - 98.6x41 mm	35 mm	box of 100
27082	Single blade* with comb - 103x47 mm	38 mm	box of 50

*Single blade on both sides

NON STERILE DISPOSABLE RAZORS

Non sterile razors with stainless steel blade.

All razors are supplied in a multilingual dispenser box:

GB, FR, IT, ES, PT, DE, GR, PL, SE, RO, HR, FI, SI, LV, CZ, HU, BG, Arabic.

Single blade

Double blade

Single blade on both sides



27080



27081



27082

3M™ SURGICAL CLIPPER WITH PIVOTING HEAD*



Ergonomic handle with pivoting head* and charger

3M

27054

• 27054 3M™ SURGICAL CLIPPER STARTER KIT (CLIPPER + CHARGER) - 9667L-E

The 3M surgical clipper with its unique pivoting head* is a safe, reliable and easy to use device providing a comfortable clipping experience for patients. *When used with 3M clipper blades 9660 (27055)

SURGICAL CLIPPER FEATURES WITH DIFFERENT SURGICAL CLIPPER BLADES

	27055	27056
Pivoting head	✓	✓
Ergonomically designed handle that's easy to hold	✓	✓
Easy-to-attach blade	✓	✓
Cuts through all types of head and body hair, wet or dry, in a single pass	✓	✓
Quickly and effectively removes hair in trauma and neuro cases	✓	✓
Away-from-the-skin cutting blade minimizes potential for nicks or cuts, which helps reduce the risk of SSIs (surgical site infections)	✓	✓
Li-ion battery designed to hold a charge longer and battery indicator light	✓	✓
Approximately 160 minutes of run time on a four-hour recharge	✓	✓
Submersible for easy cleaning (meets IPX7 rating of submersible in water no deeper than 1 m for up to 30 min)	✓	✓

• 27055 3M™ SURGICAL CLIPPER BLADES - 9660 - pivoting - box of 50

• 27056 3M™ SURGICAL CLIPPER BLADES - 9690 - non pivoting - box of 20

Single use clipper blade individually sealed



27055

Clipper blade for head and heavy body hair



27056



SURGICAL NON ABSORBABLE MESH - TRANSPARENT -STERILE

Multilingual manual and box:
GB, FR, IT, ES, PT, DE, GR.

TRULENE (Polypropylene) MESH - Non-Absorbable

Equivalent Ethicon brand: PROLENE MESH

Features

- non-adhering to tissue and biologically inert
- gives strong support for indefinite periods
- excellent tensile strength and pliability
- easy to draw through tissues
- inert material
- does not harbor micro-organisms
- least thrombogenic

Applications

- hernia repair, plastic, general and gastrointestinal surgery

Technical details

- monofilament polypropylene non-absorbable surgical sutures
- the product does not require any coating is sold untreated
- sterilized by ethylene oxide
- Shelf life of 5 years.

MDR Certified

TRULENE MESH - non absorbable		
GIMA code	Size	Minimum order
23012	6 x 11 cm	box of 6
23013	7.6 x 15 cm	box of 6
23014	10 x 15 cm	box of 3
23015	15 x 15 cm	box of 3
23016	30 x 30 cm	box of 1

GIMA DISPOSABLE STAINLESS STEEL OR CARBON STEEL BLADES


"GIMA" DISPOSABLE BLADES - sterile - dispenser box of 100

Surgical blades are available in high quality carbon steel or stainless steel. Manufactured in accordance with requirements of British Standards BS 2982:1992, ISO 7740:1985, ISO 7153-1:1991, EN 10088-3: 2005.

Each blade is hallmarked and individually packed in aluminium foil. Carbon Steel blades are VPI paper protected.



Blade size is clearly identified on the reverse side of the foil peel-pack and on product unit boxes.

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, HU, SK, SE, FI, DK, RO, BG, GR, HR, SI, Arabic. **MDR Certified.**



DIAMANTINE™ EDGE HIGH PERFORMANCE BLADES

High quality, super sharp blades with zirconium nitride coating fused on the cutting edge to satisfy the most demanding surgeons requirements. Durable cutting edge, that stays sharp incision after incision. The diamond like material creates a very hard edge that will not dull even under the most demanding cutting needs. Bio-compatible and highly corrosion resistant. Compatible with autoclavable scalpel handles N. 3 and N. 4. Stainless steel, sterile. **MDR Certified.**



Carbon steel	Stainless steel	Diamantine stainless steel	GIMA STERILIZED BLADES	Minimum order
27000	27020	27010	Sterilized blades N. 10	box of 100
27001	27021	27011	Sterilized blades N. 11	box of 100
27002	27022	-	Sterilized blades N. 12	box of 100
27003	27023	27013	Sterilized blades N. 15	box of 100
27009	-	-	Sterilized blades N. 15c	box of 100
27004	27024	27014	Sterilized blades N. 20	box of 100
27005	27025	27015	Sterilized blades N. 21	box of 100
27006	27026	27016	Sterilized blades N. 22	box of 100
27007	27027	27017	Sterilized blades N. 23	box of 100
27008	27028	27018	Sterilized blades N. 24	box of 100

26913




GIMA code	SCALPEL HANDLES - STAINLESS STEEL	
26913	Scalpel Handle N° 3 for blades 10-15	
26914	Scalpel Handle N° 4 for blades 20-25	

GIMA DISPOSABLE STAINLESS STEEL SCALPELS

"GIMA" DISPOSABLE SCALPELS - sterile - box of 10

S/S blade fitted into an ABS plastic handle. Individually packed in sterilized envelope. Safety blade guard. **MDR Certified.**



27060-27068

STANDARD


Protection guard	Standard	Safety	Diamantine high performance	GIMA STERILIZED DISPOSABLE STAINLESS STEEL SCALPELS (box of 10 pcs.)
27040	27060	27070	22460	Disposable scalpels - blade 10
27041	27061	27071	22461	Disposable scalpels - blade 11
27042	27062	-	-	Disposable scalpels - blade 12
27043	27063	27073	22463	Disposable scalpels - blade 15
27049	-	-	-	Disposable scalpels - blade 15c
27044	27064	27074	22464	Disposable scalpels - blade 20
27045	27065	27075	22465	Disposable scalpels - blade 21
27046	27066	27076	22466	Disposable scalpels - blade 22
27047	27067	27077	22467	Disposable scalpels - blade 23
27048	27068	27078	22468	Disposable scalpels - blade 24

• 27058 STITCH CUTTER - sterile - box of 10



27058 Disposable stainless steel stitch cutter

"GIMA" SAFETY PROTECTIVE SHIELD SCALPELS - sterile
Protective scalpels with a clear view shield that protects against injuries before use, between steps of procedure and after disposal. Latex and phthalates free. **MDR Certified**



27070-27078


Multilingual box and instructions: GB, FR, IT, ES, PT, DE, NL, SI, SK, DK, LT, LV, EE, RU, RO, PL, BG, CZ, HU, SE, FI, GR, HR, Arabic. * Only for standard scalpels

DIAMANTINE HIGH PERFORMANCE SCALPELS - sterile
High-Performance Edge blade for extreme and durable sharpness

- sensory grip handle enhances cutting action and control of the line of incision
- each handle has a graduated 6 cm ruler to help determine incision length (indicative only)
- protective sheath covers blade to help prevent accidental cuts
- latex and phthalates free

MDR Certified

DIAMANTINE HIGH PERFORMANCE SCALPELS



22460-22468

Stainless steel Diamantine Edge high performance coated blade. Hardness over 2,500 Vickers

Size is clearly marked on blade and handle for easy identification

SWANN-MORTON DISPOSABLE BLADES, SCALPELS AND STITCH CUTTERS

SWANN-MORTON DISPOSABLE BLADES AND SCALPELS

Swann-Morton blades and scalpels are renowned their exceptional quality, precision, consistency and reliability. Swann-Morton manufactures in the UK using high quality

STANDARD SCALPELS MDR Certified **MADE IN THE UK**

Front



Stainless steel blade

Graduated handle

22431

Back

Plastic handle scalpels with stainless steel blades, individually packed in sterilized envelope in box of 10 pcs. Multilingual box: GB, FR, IT, ES, DE, DK, NL, GR, PT, PL, FI, RU, NO, SE.

GIMA code	SWANN-MORTON SCALPELS WITH S/S BLADE	Minimum order
22430	Disposable scalpels - blade 10 - sterile	box of 10
22431	Disposable scalpels - blade 11 - sterile	box of 10
22432	Disposable scalpels - blade 12 - sterile	box of 10
22433	Disposable scalpels - blade 15 - sterile	box of 10
22435	Disposable scalpels - blade 20 - sterile	box of 10
22436	Disposable scalpels - blade 21 - sterile	box of 10
22437	Disposable scalpels - blade 22 - sterile	box of 10
22438	Disposable scalpels - blade 23 - sterile	box of 10
22439	Disposable scalpels - blade 24 - sterile	box of 10

SAFETY SCALPELS WITH RETRACTABLE BLADE



22441

GIMA code	SWANN-MORTON SCALPELS WITH RETRACTABLE STAINLESS STEEL BLADE	Minimum order
22440	Disposable scalpels - blade 10 - sterile	box of 25
22441	Disposable scalpels - blade 11P - sterile	box of 25
22443	Disposable scalpels - blade 15 - sterile	box of 25
22448	Disposable scalpels - blade 23 - sterile	box of 25



22450

GIMA code	STITCH CUTTER WITH RETRACTABLE S/S BLADE	Minimum order
22450	Stitch cutter with retractable blade - sterile	box of 25

steel that gives exceptional cutting performance. They operate stringent quality control procedures which ensure unrivalled blade consistency and reliability. This is why Swann-Morton is the preferred choice of many medical professionals worldwide.



Carbon steel blades individually packed in metal foil and box of 100

SWANN-MORTON CARBON STEEL BLADES - STERILE - box of 100

GIMA code	SWANN-MORTON CARBON STEEL BLADES - STERILE - box of 100
22400	Blade 10
22401	Blade 11
22402	Blade 12
22403	Blade 15
22405	Blade 20
22406	Blade 21
22407	Blade 22
22408	Blade 23
22409	Blade 24

Multilingual box: GB, FR, IT, ES, DE, DK, NL, GR, PT, PL, FI, RU, NO, SE

GIMA code	CARBON STEEL STITCH CUTTER BLADE - STERILE	Minimum order
22420	Carbon steel stitch cutter blade	box of 100



22424

Handle without CE



22426

Handle with CE and UKCA

With CE	Without CE	SWANN-MORTON S/S HANDLES	Compatibility with blade No.
22425	22423*	Graduated S/S n° 3 handle	10, 11, 12, 15
22426	22424*	Graduated S/S n° 4 handle	20, 21, 22, 23, 24

*Box of 10

PARAGON BLADES AND SCALPELS

"PARAGON" DISPOSABLE S/S BLADES - sterile

High quality blades manufactured with best British stainless steel guarantee excellent sharpening and perfect cut control. Individually packed in metal foil. Gamma rays sterilized.



"PARAGON" DISPOSABLE SCALPELS WITH S/S BLADE - sterile

Grip on the plastic handle allows stability for long straight incisions. Individually packed in sterilized envelope, reinforced with blade protection, in cardboard box of 10 units.

Blades box of 100	Scalpels box of 10	STERILE, S/S BLADE
26901	26920	Blade 10
26902	26921	Blade 11
26903	26922	Blade 12
26904	26923	Blade 15T
26905	26924	Blade 20
26906	26925	Blade 21
26907	26926	Blade 22
26908	26927	Blade 23
26909	26928	Blade 24

Multilingual box: GB, FR, IT, ES, DE, DK, NL, GR, PT, PL, FI, RU, NO, SE

CE 2797

26920-26928

BLADE REMOVER

• 22455 SURGICAL BLADE REMOVER - sterile - box of 50


The Swann-Morton Blade Remover allows for safe and easy removal of all sizes scalpel blade from both No.3 and No.4 handles assisting in the reduction of accidental sharps injuries.



Locate the blade inside remover

Close remover

Withdraw the handle and dispose remover with blade



SAFETY HANDLE

• 26915 SAFETY S/S HANDLE N 3 for blades 10-15

• 26916 SAFETY S/S HANDLE N 4 for blades 20-25

Innovative handle allows hand-free removal of blade in compliance with OSHA and CDC.

MDR Certified

SHARP INJURY PREVENTION WITH SAFETY HANDLE AND BLADES



Safety handle N°3 for safety blades N°10-15

GIMA SINGLE-USE STERILE DERMAL CURETTES



Each curette is individually sealed in a peel-apart see-through package



Sterile, ready-to-use, disposable curettes eliminate need for re-sharpening and sterilization.

Seamless, stainless steel sharp cutting edge.

Easy dynamic gripping handle for sure grip.

Size embossed on each curette for easy identification.


Wide range of sizes.

Made with no natural rubber or latex.

Clinical Areas: dermatology, surgical, podiatry, oncology, veterinary.


Multilingual box and instructions:

GB, FR, IT, ES, PT, DE, PL, GR, RU, FI, SE, NL, RO, BG, HR, HU, NO, Arabic.



GIMA code	DERMAL CURETTES - STERILE	Minimum order
22572	Dermal curettes Ø 2 mm	box of 10
22573	Dermal curettes Ø 3 mm	box of 10
22574	Dermal curettes Ø 4 mm	box of 10
22575	Dermal curettes Ø 5 mm	box of 10
22577	Dermal curettes Ø 7 mm	box of 10

GIMA SINGLE-USE STERILE DERMAL CURETTES WITH SCRAPER



New dermal curette recommended by dermatologists, for skin and wound care.

Each curette is individually sealed in a peel-apart see-through package.

Clinical Areas: dermatology, surgical, podiatry, oncology, veterinary.

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, GR, RU, FI, SE, NL, RO, BG, HR, HU, NO, Arabic.

GIMA code	DERMAL CURETTES WITH SCRAPER - STERILE	Minimum order
22562	Dermal curettes with scraper Ø 2 mm	box of 10
22563	Dermal curettes with scraper Ø 3 mm	box of 10
22564	Dermal curettes with scraper Ø 4 mm	box of 10
22565	Dermal curettes with scraper Ø 5 mm	box of 10
22567	Dermal curettes with scraper Ø 7 mm	box of 10

GIMA SINGLE-USE, STERILE BIOPSY PUNCHES

Sterile, ready-to-use, disposable punches eliminate need for sterilization. Seamless, stainless steel sharp cutting edge, ribbed handle for sure grip and control.

Size embossed on each punch for easy identification.

Wide range of sizes. Each punch is individually sealed in a peel-apart see-through package.

Clinical areas: dermatology, gynaecology, podiatry.

Made with no natural rubber or latex.

Multilingual box and instructions:

GB, FR, IT, ES, PT, DE, PL, GR, RU, FI, SE, NL, RO, BG, HR, HU, NO, Arabic.

22590



GIMA code	BIOPSY PUNCHES - STERILE	Blade length mm	Minimum order
22580	Biopsy Punches Ø 1 mm	8	box of 10
22581	Biopsy Punches Ø 1.5 mm	8	box of 10
22582	Biopsy Punches Ø 2 mm	8	box of 10
22583	Biopsy Punches Ø 2.5 mm	8	box of 10
22584	Biopsy Punches Ø 3 mm	8	box of 10
22585	Biopsy Punches Ø 3.5 mm	8	box of 10
22586	Biopsy Punches Ø 4 mm	8	box of 10
22587	Biopsy Punches Ø 5 mm	8	box of 10
22588	Biopsy Punches Ø 6 mm	8	box of 10
22589	Biopsy Punches Ø 7 mm	8	box of 10
22590	Biopsy Punches Ø 8 mm	8	box of 10



STIEFEL DERMAL CURETTES AND BIOPSY PUNCHES

Ready for use, sterile, individually packed. Extremely sharp stainless steel half edge for soft tissue. Ergonomic handle for a better grip.



Size engraved on each curette for easy identification.

GIMA code	STIEFEL code	DERMAL CURETTES - STERILE	Minimum order
22642	STIEF-CUR4	Dermal curettes Ø 4 mm	box of 10
22644	STIEF-CUR7	Dermal curettes Ø 7 mm	box of 10

Disposable sterile instrument to perform an accurate skin biopsy. Stainless steel seamless sharp edge. Ribbed handle for more precise grip.



Size engraved on each punch for easy identification.

GIMA code	STIEFEL code	BIOPSY PUNCHES - STERILE	Minimum order
22650	STIEF-BP2	Biopsy punches Ø 2 mm	box of 10
22651	STIEF-BP3	Biopsy punches Ø 3 mm	box of 10
22652	STIEF-BP35	Biopsy punches Ø 3.5 mm	box of 10
22653	STIEF-BP4	Biopsy punches Ø 4 mm	box of 10
22654	STIEF-BP5	Biopsy punches Ø 5 mm	box of 10
22655	STIEF-BP6	Biopsy punches Ø 6 mm	box of 10
22656	STIEF-BP8	Biopsy punches Ø 8 mm	box of 10

SINGLE-USE STERILE DERMAL CURETTES



Stainless steel sharp cutting edge for soft or weak tissue. Size embossed on each curette for easy identification. Sterile ready-to-use. Use once and discard eliminates need for re-sharpening and sterilization. Wide range of sizes.

Used in the following minor surgical skin procedures:

Curettage of basal cell carcinoma, verruca and any other benign neoplasms of the skin.

Clinical Areas: dermatology, surgical, podiatry, veterinary.

MADE IN JAPAN


GIMA code	DERMAL CURETTES	Minimum order
22602	Dermal Curettes Ø 2 mm	box of 20
22603	Dermal Curettes Ø 3 mm	box of 20
22604	Dermal Curettes Ø 4 mm	box of 20
22605	Dermal Curettes Ø 5 mm	box of 20
22606	Dermal Curettes Ø 7 mm	box of 20



SINGLE-USE STERILE BIOPSY PUNCHES - LATEX FREE

Stainless steel sharp cutting edge, ribbed handle for sure grip. Size embossed on each punch for easy identification. Suggested for dermatology, gynaecology, podiatry. Each punch is individually sealed in a peel-apart see-through package.


Standard	With plunger	BIOPSY PUNCHES	Blade length mm	Minimum order
22610	22630	Biopsy Punches Ø 1 mm	7	box of 20
22611	22631	Biopsy Punches Ø 1.5 mm	7	box of 20
22612	22632	Biopsy Punches Ø 2 mm	7	box of 20
22613	-	Biopsy Punches Ø 2.5 mm	7	box of 20
22614	22634	Biopsy Punches Ø 3 mm	7	box of 20
22615	-	Biopsy Punches Ø 3.5 mm	7	box of 20
22616	22636	Biopsy Punches Ø 4 mm	7	box of 20
22617	-	Biopsy Punches Ø 5 mm	7	box of 20
22618	-	Biopsy Punches Ø 6 mm	7	box of 20
22619	-	Biopsy Punches Ø 8 mm	7	box of 20



Internal Plunger System
Ejects lodged skin specimen inside the metal lumen of the punch



Punches are individually packed in a transparent, sturdy, blister in a 20 punches multilingual dispenser (GB, FR, DE, ES, PT, JP)



URINE CONTAINERS



Graduated urine containers made of polypropylene with red screw cap and frosted label.

One multilingual instruction in each carton: GB, FR, IT, ES, PT.

GIMA code	URINE CONTAINERS - BULK	Minimum order
25958	Urine containers 120 ml	box of 300
URINE CONTAINERS - individually wrapped		
25960	Urine containers 60 ml - cleanroom ISO8	box of 500
25961	Urine containers 60 ml - sterile	box of 500
25962	Urine containers 120 ml - cleanroom ISO8	box of 250
25963	Urine containers 120 ml - sterile	box of 250

PLASTIC BED PAN - URINAL



• 25979 PLASTIC BED PAN - box of 20

Made of polypropylene and individually wrapped.

• 25978 MEN PLASTIC URINAL 1,000 ml - box of 60

Made of polyethylene and individually wrapped.

Bed pan and plastic urinal have multilingual polybag: GB, FR, IT, ES, PT, PL, GR, SE.

BIOLOGICAL COLLECTORS



Coprotainer® graduated faeces containers made of polypropylene with red screw cap and frosted label, individually wrapped.

One multilingual instruction in each carton: GB, FR, IT, ES, PT.

GIMA code	BIOLOGICAL CONTAINERS	Minimum order
25965	Faeces container 30 ml - cleanroom ISO8	box of 500
25966	Faeces container 30 ml - sterile	box of 500
25967	Faeces container 60 ml - cleanroom ISO8	box of 500
25968	Faeces container 60 ml - sterile	box of 500

24 HOURS URINE CONTAINERS



24 h graduated urine containers made of polyethylene.

One multilingual instruction in each carton: GB, FR, IT, ES, PT.

GIMA code	24 h URINE CONTAINERS	Minimum order
25980	Bottle 2500 ml	box of 30
25981	Container 2000 ml with ergonomic handle	box of 54
25982	Tank 2500 ml	box of 30
25984	Tank 3000 ml with aspiration system for vacuum tubes and ergonomic handle	box of 30

CONTAINERS IN SINGLE BOX - PHARMACY LINE



Graduated, sterile urine and faeces containers
Individually wrapped sterile containers with paper label (25987) or frosted label (25969, 25997) and screw cap in single box.

Graduated 24 h urine tank

Polyethylene urine tank for 24 hour urine collection, in single box

Sterile urine test tube and vacuum tube

Individually wrapped urine test tube with paper label in single box. Graduated 0.5, 1, 2.5, 5, 10 ml.

All products have multilingual box: GB, FR, IT, ES.

GIMA code	CONTAINERS	Material	Minimum order
25987	Urine container 120 ml - red cap	Polypropylene	box of 100
25969	Urine container 120 ml - yellow cap	Polypropylene	box of 100
25975	Urine vacuum tube 10 ml	Polypropylene	box of 100
25990	Urine test tube 12 ml	Polystyrene	box of 100
25995	Urine tank 24 h 2500 ml	Polyethylene	box of 27
25997	Faeces container 60 ml	Polypropylene	box of 100

PLASTIC HOLLOWARE - AUTOCLAVABLE 135°C



PLASTIC HOLLOWARE - autoclavable hundreds of times

For the protection and delivery of sterile instruments to operating rooms. For the administration of procedures, holding cleansing solutions and the administration of medicines and fluids.

4 good reasons to choose plastic line:

- 50% cheaper than stainless steel
 - autoclavable at 135°C with no limits
 - attractive modern design
 - lighter than equivalent stainless steel
- Made of polypropylene material, easy to clean, resistant to most healthcare chemicals.



- 26630 GALLIPOT/LOTION BOWL 60xh 40 mm - 50 ml - graduated
- 26631 GALLIPOT/LOTION BOWL 80xh 60 mm - 200 ml - graduated
- 26632 GALLIPOT/LOTION BOWL 100xh 70 mm - 300 ml - graduated
- 26633 GALLIPOT/LOTION BOWL 150xh 70 mm - 500 ml - graduated
- 26635 KIDNEY DISH 200xh 45 mm - 500 ml - graduated
- 26636 KIDNEY DISH 250xh 55 mm - 750 ml - graduated
- 26638 INSTRUMENT TRAY 200x150x51 mm
- 26639 INSTRUMENT TRAY 300x250x52 mm
- 26640 COMPARTMENT TRAY 270x180x41 mm

VALUE PLASTIC HOLLOWARE - AUTOCLAVABLE 134°C

Economic line of reusable plasticware products.

Made of polypropylene, easy to clean and resistant to most

healthcare chemicals. Autoclavable at 121°C (15 min), 126°C (10 min) or 134°C (3 min).



- 37701 GALLIPOT 63x42 mm - 90 ml
- 37702 GALLIPOT 90x60 mm - 240 ml



- 37705 KIDNEY TRAY 6" 155x75 mm
- 37706 KIDNEY TRAY 8" 205x100 mm
- 37707 KIDNEY TRAY 10" 260x125 mm
- 37708 KIDNEY TRAY 12" 306x140 mm



- 37710 MEASURING JUG 500 ml
 - 37711 MEASURING JUG 1,000 ml
 - 37712 MEASURING JUG 2,000 ml
- Raised graduation for easy reading.



- 37715 LABORATORY TRAY 375x300xh 75 mm



- 37717 COMPARTMENT TRAY 266x175xh 42 mm



Provided with cover which fits into groove of the tray to avoid minimum spillage of material.

- 37718 INSTRUMENT TRAY 220x150xh 70 mm

RECYCLED CELLULOSE PAPER HOLLOWARE - WATERPROOF

Single use urine and liquid containers, made of recycled cellulose paper, paper binders and adjuvants. Safe, clean, atoxic and latex free. Biocompatibility certificate in accordance with ISO 10993-10:2010. Containers can be subjected to ETO standard sterilization. Multilingual instruction and label: GB, FR, IT, ES, PT, DE, PL, RO, GR, SE.

Material:
Intended use:
Resistance:
Disposal:
Medical device:

TECHNICAL SPECIFICATIONS

moulded paper pulp
containment and transport of urine, blood, vomit
36 hours to water, 4 hours to organic liquids
biodegradable
MDR 2017/745 - Class I non sterile

MEDICAL USE



- 26645 DISPOSABLE KIDNEY DISH 0.65 l - box of 300 pcs
- Size: 24.7x12.2x4.4 cm



- 37680 DISPOSABLE BASIN 3l - box of 100 pcs
- Size: 29.9x25.4x10.2 cm



- 37681 DISPOSABLE GRADUATED JUG 1.5 l - box of 150 pcs
- Size: 23.9x16.7x11.6 cm



- 37682 DISPOSABLE MALE PORTABLE URINAL 0.9 l - box of 100 pcs
- Size: 26x11.5x11 cm

S/S KIDNEY DISHES - AUTOCLAVABLE 121°C

Kidney dishes have
curved edges

26642

26644

26641

26643

ECONOMIC S/S KIDNEY DISHES

S/S Kidney dishes, autoclavable 121°C, with 0.5 mm thickness

GIMA code	ECONOMIC S/S KIDNEY DISHES	Capacity
26641	Tray 6" - 162x77x31 mm	190 ml
26642	Tray 8" - 207x98x39 mm	400 ml
26643	Tray 10" - 247x122x43 mm	600 ml
26644	Tray 12" - 309x149x59 mm	1500 ml

BUCKETS WITH LID - STAINLESS STEEL

• 26574 COMMODE BUCKET WITH LID 5 l

- Ø 307xh 167 mm

• 26576 BUCKET WITH LID 12 l

- Ø 290xh 300 mm

• 26577 BUCKET WITH LID 15 l

- Ø 322xh 307 mm

Stainless Steel AISI 304 buckets
autoclavable at 121°C.

26574

26576

26577

S/S AISI 304 - 18/8 KIDNEY DISHES - THICKNESS 0.6 MM - AUTOCLAVABLE 121°C

Flat edges

26610

Flat edges

26611

26612

Curved edges

GIMA code	SHALLOW KIDNEY DISHES	Capacity
26610	Tray 8" - 207x128xh 33 mm	480 ml
26611	Tray 10" - 254x141xh 33 mm	750 ml
26612	Tray 12" - 280x141xh 40 mm	1600 ml

26580

26581

26582

Curved edges

GIMA code	DEEP KIDNEY DISHES	Capacity
26580	Tray 6" - 162x77x31 mm	190 ml
26581	Tray 8" - 207x98x39 mm	400 ml
26582	Tray 10" - 247x122x43 mm	600 ml
26583	Tray 12" - 309x149x59 mm	1500 ml

26584

26585

26589

Curved edges

GIMA code	DEEP WITH COVER* KIDNEY DISHES	Capacity
26584	Tray 7" - 200x95x35 mm	300 ml
26585	Tray 10" - 247x122x43 mm	600 ml
26589	Tray 12" - 309x149x59 mm	1500 ml

*Thickness of cover 0.5 mm

S/S AISI 304 - 18/8 WASH BASINS AND GALLIPOTS - AUTOCLAVABLE 121°C

26586

26587

GIMA code	WASH BASIN	Capacity
26586	Ø 310x h 72 mm - shallow	3 l
26588	Ø 405x h 95 mm - shallow	7 l
26587	Ø 318x h 84 mm - deep	4 l

26621

26622

GIMA code	ROUND GALLIPOTS	Capacity
26595	Ø 107xh 56 mm	325 ml
26596	Ø 158xh 75 mm	940 ml
26621	Ø 208xh 84 mm	1940 ml
26622	Ø 258xh 92 mm	3200 ml

26618

26619

26616

26617

GIMA code	GALLIPOTS WITH LIP	Capacity
26616	Ø 56xh 32 mm	50 ml
26617	Ø 88xh 46 mm	170 ml
26618	Ø 128xh 60 mm	450 ml
26619	Ø 158xh 75 mm	940 ml

STAINLESS STEEL HOLLOWARE - AUTOCLAVABLE 121°C - AISI 304



STAINLESS STEEL HOLLOWARE
Stainless steel AISI 304 holloware with minimum thickness 0.6 mm.

GIMA code	INSTRUMENT TRAYS - ANGLED EDGES	
	Without cover	
26603	306x196x h 50 mm	
26604	355x254x h 50 mm	
26597	380x304x h 50 mm	
26598	440x320x h 64 mm	
GIMA code	With cover	
	Perforated with cover	
26607	306x196x h 50 mm	
26608	355x254x h 50 mm	
26609	440x320x h 64 mm	
26646	223x126x h 45 mm	
26647	355x254x h 50 mm	

Supplied with cover



GIMA code	INSTRUMENT TRAYS ROUND EDGES	
	Without cover	
26599	223x126x h 45 mm	
26602	264x172x h 47 mm	
GIMA code	With cover	
26605	223x126x h 45 mm	
26606	264x172x h 47 mm	



GIMA code	DENTAL TRAY	
26600	208x109x h 15 mm	
GIMA code	INSTRUMENT TRAYS	
26601	210x160x h 25 mm	
26629	313x218x h 31 mm	



GIMA code	MAYO TRAYS - SHALLOW	
	Non perforated	
26613	254x165x h 18 mm	
26614	350x252x h 16 mm	
26615	432x295x h 19 mm	
GIMA code	Perforated	
26651	350x252x h 16 mm	



- 26573 ROUND BED PAN WITH COVER - straight handle - Ø 300xh 85 mm



- DRESSING JAR WITH LID
- 26623 DRESSING JAR - 0.5 l - Ø 102x65 mm
- 26624 DRESSING JAR - 1 l - Ø 103x126 mm
- 26625 DRESSING JAR - 2 l - Ø 127x161 mm



- 26626 FORCEPS JAR - Ø 55xh 140 mm
- 26627 FORCEPS JAR - Ø 55xh 180 mm
- 26628 THERMOMETER JAR - Ø 33xh 80 mm
- 26620 MEDICINE CUP - 60 ml - Ø 53xh 53 mm



- 26590 BED PAN FRACTURE TYPE - with ring - 418x292x85 mm



- 26594 BED PAN FRACTURE TYPE - with tube handle - 418x292x85 mm



- 26593 BED PAN PERFECTION TYPE - 354x289x97 mm



- 26578 MALE URINAL - 1.2 l



- 26579 FEMALE URINAL - 1.2 l




- 26591 MALE VERTICAL URINAL - Ø 114x220 mm - 1.5 l

HYPODERMIC SYRINGES IN 3 PARTS WITH AND WITHOUT NEEDLE

High-quality range of medical devices ensuring an optimal safety. The range of Hypodermic syringe INJ/LIGHT have been manufactured in compliance with the most important international standards in conformity with European Directive 93/42/EEC and with UNI EN ISO 9001/2015 - UNI EN ISO 13485/2016.

The INJ/LIGHT line is completely non-toxic and non-pyrogenic, made of transparent polypropylene with graduated black colour scale and individually packed in blister pack. Needles in AISI 304 stainless steel triple bevelling. Box in GB, FR, IT, ES, DE, PT, GR, RO, RU, YU, SA.

SYRINGE WITHOUT NEEDLE

GIMA code	Capacity	Minimum order
LUER LOCK CONE		
23806	1 ml	box of 100
23807	5 ml	box of 100
23808	10 ml	box of 100
23809	20 ml	box of 50
23810	30 ml	box of 50
23811	60 ml	box of 25
		
CENTRIC LUER CONE		
23813	3 ml	box of 100
23814	5 ml	box of 100
23815	10 ml	box of 100
23812	20 ml	box of 50
		
ECCENTRIC LUER CONE		
23816	10 ml	box of 100
23817	20 ml	box of 50
23819	60 ml	box of 25
		
CATHETER LUER CONE		
23820	60 ml	box of 20
		
23821	100 ml	box of 25
		

Eccentric Luer Central Luer Luer Lock



SYRINGE WITH NEEDLE

GIMA code	Capacity	Assembled needle	Minimum order
CENTRIC LUER CONE			
23822	1 ml - tuberculin	26G - 0.45x13 mm	box of 100
23824	2.5 ml	21G - 0.8x38 mm	box of 100
23825	2.5 ml	22G - 0.7x32 mm	box of 100
23827	5 ml	21G - 0.8x38 mm	box of 100
23828	5 ml	22G - 0.7x32 mm	box of 100
23830	10 ml	21G - 0.8x38 mm	box of 100
23831	10 ml	22G - 0.7x32 mm	box of 100
23832	20 ml	21G - 0.8x38 mm	box of 50



SYRINGE WITH SAFETY NEEDLE

GIMA code	Capacity	Assembled needle	Minimum order
CENTRIC LUER CONE			
23835	2.5 ml	22G - 0.7x32 mm	box of 100
23836	5 ml	22G - 0.7x32 mm	box of 100
23837	10 ml	21G - 0.8x38 mm	box of 50
23838	10 ml	22G - 0.7x32 mm	box of 50



**STOP ACCIDENTAL
NEEDLE STICKS**

ECCENTRIC LUER CONE WITH NEEDLE

GIMA code	Capacity	Assembled needle	Minimum order
23833	10 ml	21G - 0.8x38 mm	box of 100
23834	20 ml	21G - 0.8x38 mm	box of 50



STERILE EO

PHT

LATEX

CE

1370

Medical Device class IIa

HYPODERMIC NEEDLES AND SCALP-VEIN SETS - STERILE

A line of hypodermic needles and scalp vein sets ensuring minimal trauma during penetration thanks to the triple bevelling. Each element of the needle is treated for painless injection. Needles are in stainless steel AISI 304 and totally non-toxic and non-pyrogenic.

Needles have Luer-slip cone, can be mounted on syringes with Luer-slip or Luer-lock. The scalp vein set needles are provided with soft flexible grained wings. The tube is made of 30 cm PVC with luer-lock. Individually packaged in blister. Box and instructions: GB, FR, IT, DE, ES, PT, GR, RU (Scalp-vein sets), GB, IT, ES, PT (Hypodermic needles).



Hypodermic needle



23793
Hypodermic
safety needle

HYPODERMIC NEEDLES

Standard	Safety	Needle	Colour	Minimum order
23750	23790	18G - 1.2x38 mm	Pink	box of 100
23751	23791	19G - 1.1x38 mm	Cream	box of 100
23753	-	20G - 0.9x38 mm	Yellow	box of 100
23755	23793	21G - 0.8x38 mm	Green	box of 100
23757	-	22G - 0.7x32 mm	Black	box of 100
23758	23794	23G - 0.6x25 mm	Turquoise	box of 100
23760	-	25G - 0.5x16 mm	Orange	box of 100
23761	-	26G - 0.45x13 mm	Brown	box of 100
23762	-	27G - 0.4x13 mm	Grey	box of 100




STERILE EO

LATEX

CE

Medical Device

1370



SCALP - VEIN SETS

Standard	Safety	Needle	Colour	Minimum order
23766	-	19G - 1.1x19 mm	Cream	box of 100
23767	-	20G - 0.9x19 mm	Yellow	box of 100
23768	23774	21G - 0.8x19 mm	Green	box of 100
23769	-	22G - 0.7x19 mm	Black	box of 100
23770	23776	23G - 0.6x19 mm	Turquoise	box of 100
23771	-	25G - 0.5x19 mm	Orange	box of 100

TERUMO® AGANI™ HYPODERMIC NEEDLES - STERILE

Sterile hypodermic needles for single use, suitable for a wide range of applications.

Special geometry of the sharp lancet point facilitates the penetration of the skin and helps to minimize the trauma. Silicone-coated cannula provides good gliding properties for a smooth penetration and prevents any knocking effect.



Secure and easy handling increases safety and user comfort: the steeper point angle of the short bevel shortens the length of the needle bevel, reducing the risk of vascular perforation. Hub made of polypropylene, needle made of stainless steel. ISO colour-coded hub allows to recognize the diameter of the needle. Class IIa. MDR 2017/745. Shelf life at production: 5 years. Individually packed in blister. Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT, TR, RU.



AGANI™ HYPODERMIC NEEDLES				
GIMA Code	TERUMO® Code	Needle	Colour	Minimum order
37600	AN*1838R1	18Gx1 1/2" - 1.2x38 mm	pink	box of 100
37605	AN*2138R1	21Gx1 1/2" - 0.8x38 mm	green	box of 100
37607	AN*2232R1	22Gx1 1/4" - 0.7x32 mm	black	box of 100
37609	AN*2325R1	23Gx1" - 0.6x25 mm	blue	box of 100
37610	AN*2332R1	23Gx1 1/4" - 0.6 x 32 mm	blue	box of 100
37611	AN*2425R1	24Gx1" - 0.55 x 25 mm	purple	box of 100
37612	AN*2516R1	25Gx5/8" - 0.5 x 16 mm	orange	box of 100
37613	AN*2525R1	25Gx1" - 0.5 x 25 mm	orange	box of 100
37616	AN*2719R1	27Gx3/4" - 0.4 x 19 mm	grey	box of 100
37617	AN*3013R1	30Gx1/2" - 0.3 x 13 mm	yellow	box of 100



37616

TERUMO® 3 PART SYRINGES WITH AND WITHOUT NEEDLE

Sterile 3-part syringes for single use, with or without pre-connected sterile hypodermic needle. Syringes are intended for the aspiration of fluid or blood, or for the injection of fluids immediately after filling into the human body. Hub is made of polypropylene and master batch, while cannula is made of stainless steel. Latex free.

Transparent barrel provides easy recognition of air bubbles

and the clear printing allows for precise dosage management. Class IIa, Medical Device following MDR 2017/745.

Shelf life at production: 5 years. Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT, TR, RU.



TERUMO® SYRINGES WITHOUT NEEDLE - STERILE

GIMA code	TERUMO® Code	Capacity	Minimum order
CONCENTRIC LUER LOCK			
37630	MDSS03LE	3 ml	box of 100
37632	SS*05LE1	5 ml	box of 100
37634	SS*10LE1	10 ml	box of 100
37635	SS+20L1	20 ml	box of 50
37637	SS+30L1MP	30 ml	box of 25
37638	SS+50L1	50 ml - transparent Polypropylene barrel	box of 25
CONCENTRIC LUER SLIP			
37640	MDSS01SE	1 ml - Tuberculin	box of 100
37641	MDSS03SE	3 ml	box of 100
37643	MDSS05SE	5 ml	box of 100
ECCENTRIC LUER SLIP			
37645	MDSS10ESE	10 ml	box of 100
37646	MDSS20ESE	20 ml	box of 50
CATHETER TIP			
37648	SS+50C1	50 ml	box of 25

TERUMO® SYRINGES WITH NEEDLE - STERILE

GIMA code	TERUMO® Code	Capacity	Assembled needle	Minimum order
37652	SS*02SE21381	2.5 ml	21Gx1 1/2" - 0.8x40 mm	box of 100
37653	SS*02SE22381	2.5 ml	22Gx1 1/2" - 0.7x38mm	box of 100
37654	MDSS03S2138E	3 ml	21Gx1 1/2" - 0.8x38mm	box of 100
37655	MDSS05S2238E	5 ml	21Gx1 1/2" - 0.8x38mm	box of 100
37656	MDSS05S2332E	5 ml	22Gx1 1/2" - 0.7x38mm	box of 100
37659	MDSS10S2138E	10 ml	21Gx1 1/2" - 0.8x38mm	box of 100



37655

TUBERCULIN

37662	MDSS01S2613E	1 ml	26Gx1/2" - 0.45x13mm	box of 100
-------	--------------	------	----------------------	------------



37662

TUBERCULIN LUER

37664	MDSS01S2516E	1 ml	25Gx5/8" - 0.5x16mm	box of 100
-------	--------------	------	---------------------	------------



37648




INTRAMUSCULAR INJECTION TRAINING PAD

• 23797 INTRAMUSCULAR INJECTION TRAINING WEARABLE PAD - with strap
Intramuscular injection training pad, with strap for fixation around body limbs, to perform injection exercises. To increase realism of operation, the module is divided into skin, muscle layer and outer frame card structure. Surface is made of silicone with internal sponge, base is made of plastic.

23797

Size: 8.5x7.5x3 cm
Weight: 0.08 kg



18 cm strap
with velcro
closure

BD EMERALD™ SYRINGES

Minimizes Waste Disposal
Has up to 30% less material than other syringes

Flange Size and Shape
To provide stability and comfort during use

Textured Thumb Press
To reduce slippage during administration

CE 0318

Green Stopper and Bold Scale Markings
To aid in setting an accurate dose

Clear Barrel
To allow visualization of fluid

Robust Retaining Ring
To help prevent accidental plunger rod pullout

BD EMERALD™ SYRINGES

GIMA code	BD code	Capacity	NEEDLE	Minimum order
SYRINGES WITHOUT NEEDLES - LUER SLIP CENTRIC TIP				
23863	307727	2 ml		box of 100
23864	307731	5 ml		box of 100
23865	307736	10 ml		box of 100
SYRINGES WITH NEEDLES - LUER SLIP CENTRIC TIP				
23868	307728	2 ml	22G - 0.7x30 mm	box of 100
23871	307732	5 ml	21G - 0.8x40 mm	box of 100
23870	307737	10 ml	21G - 0.8x40 mm	box of 100

BD PLASTIPAK™ SYRINGES

CE 0318

23873



23878

BD PLASTIPAK™ SYRINGES

GIMA code	BD code	Capacity	MODEL	Minimum order
SYRINGES WITHOUT NEEDLE - LUER LOCK				
23873	301229	30 ml	Luer - lock	box of 60
23874	300865	50 ml	Luer - lock	box of 60
SYRINGES WITHOUT NEEDLE - CATHETER CONE				
23877	300867	50 ml	Catheter cone	box of 60
23878	300605	100 ml	Catheter cone	box of 25


BD MICROLANCE™ 3 NEEDLES

BD Microlance™ 3 Needles are designed for easy penetration, smooth gliding force during insertion/withdraw and reduced cannula friction while sliding through tissue.

GIMA code	BD code	BD MICROLANCE™ 3 NEEDLES	Minimum order
23739	303262	18G - 1.20x40 mm - pink	box of 100
23742	304432	21G - 0.80x40 mm - green	box of 100
23744	300800	23G - 0.60x25 mm - blue	box of 100
23745	300600	25G - 0.50x16 mm - orange	box of 100
23747	300635	27G - 0.40x13 mm - grey	box of 100
23748	304000	30G - 0.29x13 mm - yellow	box of 100



BD MICRO-FINE™ INSULIN NEEDLES AND SYRINGES




GIMA code	BD code	BD MICRO-FINE™ NEEDLES	Minimum order
23850	320140	BD Micro-Fine needles 4 mm - 32G	box of 100
23851	320594	BD Micro-Fine needles 5 mm - 31G	box of 100
23852	320592	BD Micro-Fine needles 8 mm - 31G	box of 100

BD MICRO-FINE™ SYRINGES			
23860	324826	BD Micro-Fine syringes 0.3 ml - 8 mm - 30G	box of 100
23862	324827	BD Micro-Fine syringes 0.5 ml - 8 mm - 30G	box of 300

Compatible with most of insuline pens available on the market (Sanofi-Aventis, Eli Lilly, Novo Nordisk, Owen Mumford)

INSULIN SYRINGE INSU/LIGHT



No-dead space syringe has double protector (needle cover and plunger cover) for better sterility

Available with assembled needle or fixed needle both with no-dead space. Insulin syringes with fixed needle have the needle directly fixed on the cone and have a double protection (needle cap and plunger cap), assuring the sterility of the device. Non-toxic and non-pyrogenic, made of transparent polypropylene with graduated insulin scale. Needle in stainless steel AISI 304 with triple bevelling.

GIMA code	INSULIN SYRINGE NO DEAD SPACE - ASSEMBLED NEEDLE	Minimum order
23802	Insulin syringe 1 ml - 25G - 0.5x16 mm	box of 100
23803	Insulin syringe 1 ml - 27G - 0.4x13 mm	box of 100
INSULIN SYRINGE NO-DEAD SPACE - FIXED NEEDLE		
23805	Insulin syringe 0.5 ml - 30G - 0.3x8 mm	box of 100

INSULIN PEN NEEDLES

Multilingual box: GB, FR, IT, DE, ES, PT, GR, PL, CZ, Arabic



STERILE EO CE 0197

Insulin Pen Needles are disposable, sterile needles used with insulin pens to control high blood sugar in adults and children with diabetes mellitus. With a fine design to minimize the discomfort of injection, ultra-thin needles provide painless and a non traumatic injection due to a silicone film. Pyrogen free.

GIMA code	GIMA INSULIN PEN NEEDLES - UNIVERSAL	Minimum order
23840	GIMA insulin pen needles 4 mm - 32G	box of 100
23841	GIMA insulin pen needles 5 mm - 31G	box of 100
23842	GIMA insulin pen needles 6 mm - 31G	box of 100
23843	GIMA insulin pen needles 8 mm - 31G	box of 100

Compatible with most of insulin pens available on the market (Sanofi-Aventis, Novo Nordisk, Eli Lilly, Owen Mumford)

MESO-RELLE MESOTHERAPY, SCLEROTHERAPY, FILLER NEEDLES



MESO-RELLE HYPODERMIC NEEDLES

Meso-relle disposable hypodermic needles are multiple sharpened to facilitate penetration. Stainless steel AISI 304 cannula and polypropylene connection for syringe. Sterilized with ethylene oxide. Sold in individual blister made of medical paper. **MADE IN JAPAN**



GIMA code	MESOTHERAPY LUER NEEDLES				Minimum order
	Gauge	Ø mm	Length mm	Colour	
23670	27G	0.40	4	grey	box of 100
23671	27G	0.40	6	grey	box of 100
23672	30G	0.30	4	yellow	box of 100
23673	30G	0.30	6	yellow	box of 100
23674	31G	0.26	4	light blue	box of 100
23675	31G	0.26	6	light blue	box of 100
23676	32G	0.23	4	pink	box of 100
23677	32G	0.23	6	pink	box of 100
23678	33G	0.20	4	green	box of 100
SCLEROTHERAPY/FILLER LUER NEEDLES					
23679	27G	0.40	12	grey	box of 100
23680	30G	0.30	12	yellow	box of 100
23682	31G	0.26	12	light blue	box of 100
23684	32G	0.23	12	pink	box of 100
FILLER LUER NEEDLES					
23686	27G	0.40	40	grey	box of 100
23687	30G	0.30	25	yellow	box of 100
23688	30G	0.30	40	yellow	box of 100

MICRO-MESOTHERAPY



MESO-RELLE MICRO-MESOTHERAPY

Microtherapy, an evolution of mesotherapy, is a recent technique for administration of drugs via superficial intraderma injection, improving intramuscular intra-venereal and subcutaneous injection technique.




Single use device with reversed cone shape and special ultrafine needle in its centre

GIMA code	MICRO-MESOTHERAPY NEEDLES - STERILE					Minimum order
	Type	Gauge	Ø mm	Length mm	Colour	
23690	Needle	27G		2.5	white	box of 25
23691	Needle	30G		2.5	yellow	box of 25
23694	S.I.T.*	31G	0.26	2.5	white	box of 25

*Skin Injection Therapy

MESO-RELLE MULTI-INJECTORS



MESO-RELLE MULTI-INJECTOR CONNECTORS

Disposable mesotherapy multi-injectors connectors sterilized with ethylene oxide. Packed in individual, transparent blister, colourless and transparent to enable full examination of the inner chamber. Perfectly resistant to high pressure.



GIMA code	MULTI-INJECTOR CONNECTORS				Minimum order
	Type	Way-out	Taper		
23700	Straight	3	Luer		box of 50
23701	Straight	5	Luer		box of 50
23703	Circular	5	Luer		box of 50

OPTIVA® I.V. CATHETERS WITH STRAIGHT HUB OR SIDE PORT - RADIOPAQUE

High flow rate
thin wall kink-resistant polyurethane and FEP radiopaque catheters

Secure grip and intuitive handling
ribbed housing and bevel-up indicator

Needle-in-vein aspiration
removable flash plug

Easy vein entry
ultrasharp V-point introducer needle reduces penetration force by 25 to 30%

Patient comfort and secure fixation
long hub for generous suturing or taping

Fast visualization
translucent flash chamber

Hygienic sideport access
with Snap Cap design to reduce disturbance to the catheter during injection valve usage

Stabilization

Intermittent therapies

OPTIVA® STRAIGHT HUB CATHETER - STERILE

OPTIVA® 2 SIDEPORT CATHETER - STERILE

MADE IN ITALY

High flow rate

Secure fixation

Easy vein entry

Hygienic sideport access

Stabilization

Intermittent therapies

GIMA code	Smiths code	I.V. CATHETER WITH STRAIGHT HUB			Minimum order
		Gauge	Gauge mm	Colour	
23720	5064-AI	18G	45 mm (1 3/4")	Green	box of 50
23721	5066-AI	20G	32 mm (1 1/4")	Pink	box of 50
23722	5060-AI	22G	25 mm (1")	Blue	box of 50

GIMA code	Smiths code	I.V. CATHETER WITH SIDE PORT			Minimum order
		Gauge	Gauge mm	Colour	
23730	1018-AI	18G	45 mm (1 3/4")	Green	box of 50
23731	1020-AI	20G	32 mm (1 1/4")	Pink	box of 50
23732	1022-AI	22G	25 mm (1")	Blue	box of 50

BD VENFLON™ IV CATHETERS - STERILE

BD VENFLON™ IV CATHETERS

I.V. Cannula Luer-Lock with a plug attached to the flow control plug. Designed to gain access to peripheral veins of the patient's blood system for rehydration, parenteral nutrition, medication delivery, blood transfusion, and monitoring purposes.

GIMA code	BD code	BD VENFLON™ I.V. CATHETERS					Minimum order
		Gauge	Gauge mm	Flow rate	Colour		
23715	391453	18G	45 mm (1 3/4")	80	Green		box of 50
23716	391452	20G	32 mm (1 1/4")	54	Pink		box of 50
23717	391451	22G	25 mm (1")	31	Blue		box of 50

Flash Back Function - During the insertion of the cannula, blood flow is visualized by blood flash back in the hub of the cannula



23716

Multilingual box and instructions



INFUSION SETS

INFUSION SETS

Infusion set includes:
ABS spike with air vent,
6 cm PVC chamber, 150 cm
PVC tube, ABS roller,
ABS luer lock terminal
without needle, latex free injection site.
Sterilized by Ethylene oxide.

- 26062 Y INFUSION SET - polybag - box of 500 (20 bags of 25)
- 26069 HOSE INFUSION SET - polybag - box of 500
- 26064 HOSE INFUSION SET - blister - box of 360 (18 bags of 20)



26062

NEEDLE-LESS VALVE

- 26059 NEEDLE-LESS VALVE - sterile - box of 500

Device for I.V. infusion and/or withdrawal,
needle-less closed circuit.

For protection of vascular access and of all
infusion systems equipped with male luer-lock
standard connections.

This closed and self-sealing protection
system allows to infuse or inject fluid,
make withdrawals eliminating the risk
of accidental puncture to operators
and reducing the risk of contamination
to patients.

Medical Device Class I.

Latex free
PVC free

26059



INFUSION SET FOR PUMP

- 26073 ARIES I.V. INFUSION SET - box of 150 (6 bags of 25)

Infusion set for the
administration of
solution into the
human body with
pump.

Tested with infusion pump
Ampall Argus 707 and Optima
PT.



BD SPINAL NEEDLES

Single use latex free sterile BD quincke
point spinal anesthesia needles.



BD

CE 0318



23707

GIMA code	BD code	BD QUINCKE POINT SPINAL ANESTHESIA NEEDLE - STERILE	Minimum order
23705	405248	18G 3.50 IN 1.2x90 mm - pink	box of 25
23707	405253	20G 3.50 IN 0.9x90 mm - yellow	box of 25
23709	405256	22G 3.50 IN 0.7x90 mm - black	box of 25
23710	405240	23G 3.50 IN 0.64x90 mm - turquoise	box of 25
23711	405257	25G 3.50 IN 0.5x90 mm - orange	box of 25

MULTI-VENOUS IV TRAINING ARMS



26048

Length 72.4 cm Weight 1.7 kg



Laerdal
helping save lives



26049

Length 73.7 cm Weight 1.79 kg

- 26048 LAERDAL MALE MULTI-VEIN IV TRAINING KIT

- 26049 LAERDAL FEMALE MULTI-VEIN IV TRAINING KIT

Lifelike arm reproduction with multi-vein system
designed for peripheral intravenous therapy.
Kit includes: arm, replacement skin and multi
vein system, simulated blood, blood bag with
tubing and connector, clamp and hook, 5 syringes,
manikin lubricant, carry case and directions for use.

POCKET SHARP CONTAINER

- 25780 POCKET SHARP CONTAINER 0.4 l - 12.3x6.2xh 16.5 cm - box of 70

Single use container for both sharp
and cutting waste disposal.

The permanent and inviolable
closure that prevents waste leakage
and the high container resistance
guarantee top quality and safety.

Equipped with a patented double
opening system, it can be used
by just one hand.

The double opening includes one side
with needles and insulin pen removal,
and the other side with a large
opening for entire syringe disposal.

Double opening system

Temporary closure

Definitive closure



25780

MADE IN ITALY

FEATURES 25780, 25781

Pocket size

Ergonomic shape (25780)

Practical opening system

Needles removal

Insulin pen removal device

Maximum filling level

Shock resistant

Puncture resistant

Sturdy handle

Possible side fastening

Autoclavable

at 134°C for 18'

No emission of any

harmful gas during the

incineration process

DISPO SHARP CONTAINER

- 25781 DISPO SHARP CONTAINER 0.5 l - 8.3x h 15.1 cm - box of 187

Single use conical shape puncture
resistant container for the disposal
of sharp and cutting waste.

The permanent and inviolable closure
that prevents waste leakage and the
high container resistance guarantee
top quality and safety.

Temporary and permanent closures with
a wide opening for the disposal
of larger materials, and three indents for
the needles and insulin pen removal.

Easy lid
opening

Temporary closure

Permanent closure



25781

COMPARTMENT TRAY

- 37725 COMPARTMENT TRAY 35x24,5xh 4,4 cm - box of 24

Compartment tray suitable for
carrying small containers and
minor hospital waste.

Tray fits PBS containers 25790-91-92
and CS container 25782, see page 74.



37725

+ 25790

MADE IN ITALY

PRO SPIKE FOR USED NEEDLES

- 27149 PRO SPIKE CONTAINER

27149

Small disposable plastic container for the separate
collection of spiky and cutting wastes (needles,
scalpels) on the operation field. It has a clever lid
with a safe closing system and a clip for belts and
backpacks attachment.

It can be disposed of with the contents.

Size: 90x54x185 mm

Weight: 88 g



GIMA SHARPS CONTAINERS

Full range of Sharps containers for the disposal of sharp or cutting waste. Equipped with temporary and definitive closure, certified hermetic seal, do not release any harmful emission during incineration.

The containers hold the UN mark for the transportation of dangerous materials.

- practical opening system
- temporary and permanent closure
- needles removal
- insulin pen removal device
- maximum filling level
- shock resistant
- puncture resistant
- easy to use

- easy to assemble
- excellent stability
- square or round shaped
- sturdy handle for easier transport
- customisable label for quantities
- autoclavable at 134°C for 18'
- no emission of any harmful gas during the incineration process



GIMA code	CS LINE SHARPS CONTAINERS	Minimum order	GIMA code	PBS SHARPS CONTAINERS	Minimum order
25782	CS 2 l - 14.8x14.8xh 19.5 cm	box of 60	25790	PBS 0.6 l - Ø 10xh 14.5 cm	box of 100
25783	CS 3 l - 16.7x16.7xh 23.2 cm	box of 50	25791	PBS 1.5 l - Ø 14xh 16.5 cm	box of 100
25784	CS 4 l - 16.7x16.7xh 28.2 cm	box of 40	25792	PBS 2 l - Ø 14xh 21 cm	box of 75
25785	CS 5 l - 21x21xh 20 cm	box of 30	25793	PBS 5 l - Ø 21.4xh 19.5 cm	box of 55
	CS PLUS LINE SHARPS CONTAINERS		25796	PBS 6 l - Ø 21.4xh 24 cm	box of 55
25787	CS Plus 10 l - 22x22xh 32.5 cm	box of 33	25794	PBS 7 l - Ø 25xh 19.5 cm	box of 40
25788	CS Plus 22 l - 37.5x27.2xh 34.5 cm	box of 10	25795	PBS 12 l - Ø 25xh 31.3 cm	box of 25
25799	Fixing system* - horizontal or vertical		TRAY FOR CONTAINERS		
Suitable for wall or trolley			37725	Compartment tray 32x24.5xh 4.4 cm	box of 24

CERTIFICATIONS



FDA K071517

ISO 23907-1:2019

*Suitable for 25790-2, 25782

37725 + 25790

WASTE CONTAINERS - 30-50-60 L

WASTE CONTAINERS - 3 sizes

Containers for disposal of special hospital waste in 3 sizes, with 2 different lids.

- extremely solid and shock resistant at high and low temperature
- hermetic seal by peripheral adhesive on the lid (single or double)
- lateral and central handles tested for twice the max weight
- stacking up to a maximum of 3 m with full load
- conformity: ISO 23907-1:2019

MADE IN ITALY



GIMA code	SC LINE WASTE CONTAINERS			
	Tide	Lid	Max load kg	Size
25774	Container 30 l	single	15	40x30xh 39 cm
25775	Container 30 l	double	15	40x30xh 39 cm
25776	Container 50 l	single	20	40x30xh 56 cm
25777	Container 50 l	double	20	40x30xh 56 cm
25778	Container 60 l	single	25	40x30xh 65 cm
25779	Container 60 l	double	25	40x30xh 65 cm

Single lid with adhesive


Double lid with adhesive

CONTAINERS UP TO 72 CM

LONG SHARPS CONTAINERS

Disposable containers for the elimination of long objects up to 72 cm, such as biopsy forceps and sharps. Large opening, for elimination of bulky objects, temporary and definitive closure, ergonomic handle. Autoclavable at 134°C for 18'.

MADE IN ITALY



GIMA code	LONG SHARPS CONTAINERS
25797	Magnum L 6.3 l - 17.2x12xh 61.7 cm
25798	Magnum XL 8.3 l - 17.2x12xh 84.5 cm
25799	Fixing system* - horizontal or vertical

*Suitable for wall or trolley

VEINSPY™ HAND-HELD VEIN FINDER

• 23450 VEINSPY™ HAND-HELD VEIN FINDER

VeinSpy™ is a portable vein finder for easy visualisation of veins in adult, pediatric, geriatric, obese and dark skin patients with difficult venous access. It works on Transverse Illumination Technology: light beams from a circular LED array penetrate the skin and act as a light source to achieve deep illumination of subcutaneous tissues. The veins in the illuminated area absorb the light and become prominent as dark lines.

The **VeinSpy™** device consists of two sets of coloured LEDs: one set emits orange light optimal for superficial veins, another set emits red light to find deeper veins.

Advantages of the VeinSpy™:

- reduces multiple needle prick and improves site selection for needle insertion
- no treatment delay due to difficult venous access
- reduces patient pain due to needle insertion
- increases patient comfort, avoids hematoma

• 23451 PLASTIC SLEEVES FOR 23450 - disposable - box of 50

TECHNICAL SPECIFICATIONS


Number of LED's:	22 Orange, 8 Red
Functional opening:	17 mm
Viewing diameter:	33 mm
Rechargeable battery:	Lithium Ion, 3.7V, 1150mAh
Battery charger:	110-240 V AC input, 5V DC output 750mA
Battery life:	1½ hour continuous use
Dimensions - weight:	125x65x22.7 mm - 90 g

How to use VeinSpy™

- 1- Wipe the area of the patient's skin with alcohol swab.
- 2- Insert VeinSpy™ into the disposable sleeve and switch ON.
- 3- Move the VeinSpy™ across and forth to get a clear view of the specific vein.



23450



- 4- After visualization of vein move the device so that needle can be inserted easily



- 5- Gently press the VeinSpy™ on the vein, insert needle into the vein through the open end of VeinSpy™

STANDARD ACCESSORIES

Rechargeable battery
50 disposable plastic sleeves

Battery charger
Multilingual manual: GB, FR, IT, PT, ES.

PROFESSIONAL VEIN FINDERS

• 23456 QV-500 PROFESSIONAL VEIN FINDER

• 23454 QV-600 PROFESSIONAL VEIN FINDER with display

Portable vein finding devices that can project blood vessels image on skin surface accurately and timely.


Ergonomic design with an advanced infrared imaging technology to provide real-time accurate vein digital image on patient's skin.

Both help doctors and nurses to easily find veins of obese or chemotherapy patients, hairy or dark skin people, people with varicose veins, fine blood vessels and vascular obstruction.

They can adapt 7 (23456) or 12 (23454) different colours to different ages, body shapes, skin colours, weights and various operating environments. Suitable for adults, children and infants. Both allow to observe the vascular thickness, length, curve and branch information to determine the best puncture location to avoid pain. Also available optional table stand or mobile trolley.

QV-600 additional features

- display with switchable languages (GB, FR, IT, ES, DE, PL, RU)
- possibility to photograph and store up to 6,000 images to track and analyze pathologies.



23456
+
23458



23456
+
23457


• 23457 TABLE STAND for 23454 and 23456

• 23458 MOBILE STAND for 23454 and 23456


STANDARD ACCESSORIES

Aluminium carrying case
Calibration card
Charging cable
Adapter
User manual GB, FR, IT, ES

23456



- 7 (23456) or 12 (23454) different colours to better adapt to skin and environment conditions



23454

TECHNICAL SPECIFICATIONS


	23456	23454
Light:	harmless infrared light	
Infrared wave:	760-940 nm	
Display:	-	2.4"
Image:	-	3D pseudo-colour, AI algorithm
Adjustable brightness:	5 levels	6 levels
Detection distance:	15-25 cm ± 0.3 mm	
Detection depth:	0-8 mm	0-10 mm
Noise:	≤20 dB	≤30 dB
Rechargeable Li-ion battery:	3400mAh	4800mAh
Power supply:	5V 2.0A, 100V-240V - 50/60Hz	
Size - weight	20x6.2xh 5.5 cm - 280 g	24x6xh 6.5 cm - 450 g

MILITARY TOURNIQUETS

• 25526 MILITARY TOURNIQUET - orange

• 25527 MILITARY TOURNIQUET - black

Mainly used for stopping massive haemorrhage of limbs on field conditions or in normal pre-hospital emergency. Disposable easy to use, one-handed application tourniquet improving the efficiency of hemostasis. It can be self-applied or applied by a first responder, allowing a quick and easy regulation of the pressure needed. Suitable for both upper and lower limbs. Made of woven synthetic fabric with thermoplastic clasp and secure windlass. Size: 95x4 cm. Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, Arabic.




25526



25527

25526

TOURNIQUETS IN LATEX, RUBBER AND LATEX-FREE



GIMA TOURNIQUETS

Efficient and simple solutions for professional medical use.

Latex: Ø 5 mm, for manual use

Velcro: for easy fastening

Velcro "D" buckle: self applicable made of a nylon band threads through plastic "D" buckle

Lock system: self applicable with a quick release system.

Multilingual box (excluding 25733): GB, FR, IT, DE, ES, PT, GR, Arabic.

Velcro system

VELCRO
25725

GIMA code	TOURNIQUETS	Material	Colour	Length cm	Width cm	Fastening systems	Origin
25733	TUBOLARE Tourniquet - box of 200	Latex	Yellow	33	-	Manual	Europe
25721	GIMA Tourniquet	Rubber band, terylene cover	Yellow	40	2.5	Velcro	Asia
25725	STRIP-VELCRO Tourniquet	Velcro	Green	45	2.5	Velcro+ring: use with one hand	Asia
25727	FAST Tourniquet	Rubber band, terylene cover	Red	40.5	3	Single lock system: quick release	Asia
25728	FAST Tourniquet	Rubber band, terylene cover	Green	40.5	3	Single lock system: quick release	Asia
25726	FAST Tourniquet - latex free	Spandex band, terylene cover	Blue	40.5	3	Single lock system: quick release	Asia
37500	FAST Tourniquet - daisy	Rubber band, terylene cover	Blue	40.5	3	Single lock system: quick release	Asia
37501	FAST Tourniquet - peace	Rubber band, terylene cover	Blue	40.5	3	Single lock system: quick release	Asia
37502	FAST Tourniquet - space	Rubber band, terylene cover	Blue	40.5	3	Single lock system: quick release	Asia

JETPULL 2 TOURNIQUETS WITH GRADUAL ADJUSTMENT




JETPULL 2 WITH GRADUAL ADJUSTMENT

Easy to use, allows single-hand operation.

Sophisticated 2 button system allows an easy opening of the fastening with 2 choices: smooth or fast, even when under tension. Fastening device is made using high class medical plastics (PVC free). Ribbon is in special material to avoid any skin bruises. Washable at 60/90° C and autoclavable.

Jetpull system



GIMA code	TOURNIQUETS	Material	Colour	Length cm	Width cm	Fastening systems	Origin
25722	JETPULL 2	Polyester/rubber	Blue	56	2.4	Double lock system: quick and slow release	Germany
25723	JETPULL 2 - latex free	Polyester/elastan	Green	56	2.4	Double lock system: quick and slow release	Germany
25724	JETPULL BABY	Polyester/rubber	Child design	56	2.4	Double lock system: quick and slow release	Germany

TAPE TOURNIQUET

• 37515 TAPE TOURNIQUET - latex free - disposable - box of 25

Latex free disposable tape tourniquet made of TPE with dot surface. Each box contains 25 pre-cut tourniquets 45x2.5 cm (18"x1") with 0.636 mm thickness.

Multilingual box:

GB, FR, IT, ES, PT, DE, NL, PL, RO, GR, EE, SE, HU, Arabic.



37515

PRE-CUT TOURNIQUETS

• 37510 PRE-CUT TOURNIQUET 45x2.5 cm - blue - box of 250
• 37511 PRE-CUT TOURNIQUET 45x2.5 cm - orange - box of 250
• 37512 PRE-CUT TOURNIQUET 45x2.5 cm - violet - box of 250

Latex free disposable tape tourniquet made of TPE with smooth surface. Thickness 0.636 mm

Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, RO, GR, EE, SE, HU, Arabic.



37510