



GIMA

MISSION® CHOLESTEROL METER

Code: 23932

Category: Urine and blood analyzers and test strips

Unit of sale: 1 pc.

Minimum order: 1

Type: IVD - Self Test



NSIS: 1462047

CND: W0201069099

EAN13: 0682607610907

Description: Provides high quality and reliable results for complete Lipid Panel testing in =2 min.

3 in 1 complete Lipid Panel device tests for Total Cholesterol (CHOL), High Density Lipoprotein (HDL) and Triglycerides (TRIG)

Also includes calculated tests:
Low Density Lipoprotein (LDL), CHOL/HDL ratio and Cardiac Risk Assessment

Reliable

- can test whole blood, serum or plasma
- wide hematocrit range of 0-50%
- small printer available for immediate printing of results

Easy to use

- large LCD display for easy reading
- colour coded test devices with engraved test name for easy identification

Supplied with 2 control devices, 4 batteries, carrying case and manual (GB, IT).
Available FR, ES, DE, PT manual in www.gimaitaly.com
Multilingual box: GB, FR, IT, ES, DE, PT, PL, SE, FI, GR, NL, RO, BG, RU, SA.

Technical Specifications:

Detection principle: Reflectance Photometry
Time to results: = 2 minutes
Memory with date/time: 200 results
Type of test: Quantitative
Tests: CHOL, HDL, TRIG
Calculated tests: LDL, CHOL/HDL ratio, Cardiac Risk Assessment
Specimen volume: 3-1 Lipid Panel Device: 35 µL, CHOL Individual Test Device: 10 µL
Measurement units: mg/dL or mmol/L
Measurement range:
- CHOL: 100-500 mg/dL (2.59-12.93 mmol/L)
- HDL: 15-100 mg/dL (0.39-2.59 mmol/L)
- TRIG: 45-650 mg/dL (0.51-7.34 mmol/L)
Pc interface for data transfer: Mini USB Port
Calibration: Code Chip
Languages on the screen: English abbreviations
Power source: 4 AAA or AC Adaptor (Mini USB, 5V dc, 50 mA)
Battery life: = 1,000 tests



GIMA

Automatic shut off: 5 minutes

Meter dimensions - weight: 137×79×26 mm (5.4" x 3.1" x 1.0") - 145 g (5.11 oz)

Display dimensions: 50×50 mm (2.0" x 2.0")