



GIMA

LUX MONITOR FOR LIPID, HEMOGLOBIN

Code: 23983

Category: Urine and blood analyzers and test strips

Unit of sale: 1 pc.

Minimum order: 1

Type: IVD - Self Test



NSIS: 1370069

CND: W0201069099

EAN13: 8032937231433

Description: The cardiovascular diseases are linked to the concentration of cholesterol in the blood. Many people don't know that their cholesterol is high or do not monitor it. Lux gives benefit to doctor to do a complete lipid screening plus hemoglobin in less than 3 minutes.

Easy
Only 45 µl for the complete lipid panel, 7 µl for the hemoglobin.

Effective
Doctor can test during consultation with patient.

Complete lipid panel in 3 minutes
Measurement of cholesterol, triglycerides, HDL, LDL at the same time in less than 3 minutes.

Additional hemoglobin

Software for download PC (GB, IT)
Optional data reading software (CD + cable) for easy download.

Lux monitor is supplied with carrying bag, pipette 15 µl, pipette 7 µl, 3xAAA batteries 1.5 V and manual (GB, IT).
Minimum software requirements: Windows 7.

Made in Italy.

Technical Specifications: Battery life: More than 500 measurements (lipid) or more than 1000 measurements (hemoglobine)
Method of measurements: optical method
Display: LCD (LED Backlight)
Operation buttons: 3 buttons (ON/OFF switch, up button and down button)
Memory: 1000 readings
Power source: 3XAAA 1,5 V alkaline manganese batteries