



GIMA

LACTATE SCOUT 4

Code: 23992

Category: Urine and blood analyzers and test strips

Unit of sale: 1 pc.

Minimum order: 1

Type: IVD - Self Test



NSIS: 1228660

CND: W01010218

EAN13: 8023279239928

Description: Ergonomic device for fast, accurate and reliable lactate measurement.

Easy to use:

- only 0.2 µl of capillary blood required
- pre-calibrated sensors
- device calibrated by simple coding
- lactate measurement, date/time, mode, temperature, heart rate and memory ID on a single display

Fast and accurate:

- enzymatic amperometric detection method
- results within 10 seconds
- measuring range: 0.5 - 25 mmol/L ($\pm 3\%$)

Performance measurement:

- single and step-test measurements (resting/exercise/recreation)
- stopwatch and countdown timer functions
- connection to heart rate monitors for linking to lactate value

Integrated bluetooth technology and PC software:

An USB key allows easy transfer to PC and a 4 language (GB, DE, IT, ES) software allows easy data management.

Supplied with nylon bag, empty spray bottle, Bluetooth module, USB key with software and manual, (GB, FR, IT, DE, ES, NO, SE, CZ, GR, NL, PL, TR).

Stores over 500 results

1000 tests using just 2 x CR2450 batteries

E-paper display for easier viewing on the move

Pocket size: 91 x 46 x 21 mm

Weight: 60 g

Multilingual box: GR, FR, IT, ES, DE

Integrated Bluetooth low energy connectivity