



# GIMA

## GIMA URINE ANALYZER with Bluetooth

**Code:** 24046  
**Category:** Blood and urine test  
**Unit of sale:** 1 pc.  
**Minimum order:** 1  
**Type:** IVD - Professional Use



**NSIS:** 2402832  
**CND:** W020101020101  
**EAN13:** 8023279240467

**Description:** GIMA URINE ANALYZER

Gima urine analyzer is a high-precision instrument with advanced technologies for clinical inspection of urine. It offers easy and quick operation, exact results and good repeatability.

**Features:**

- high-luminance, white LED, features long life and good stability
- large LCD display
- user-friendly interface
- monitoring the whole testing process, and audible-control prompt
- data transmission function.

**Application:**

- convenient point of care testing for doctor
- chronic diseases, long-term monitoring for nephrotic patients
- health examination
- use as telemedicine terminal.

**Tested parameters**

Glucose (GLU), Bilirubin (BIL), Specific gravity (SG), PH, Ketone body (KET), Occult blood (BLD), Protein (PRO), Urobilinogen (URO), Nitrite (NIT), Leukocytes (LEU), Ascorbic acid (VC).

Can be used with Urinalysis strips 24076.

**PHMS App**

Free App to save data, check report and results.

**Multilingual manual and box:** GB, FR, IT, ES, SE.

**Technical Specifications:**

**Working mode:** One-step  
**Internal language:** Italian and English  
**Display:** 2.4" LCD, resolution: 320x240  
**Data storage:** 500 sample data  
**Principle:** RGB tricolor test theory  
**Communication interface:** Micro USB interface, Bluetooth wireless communication  
**Repeatability:** CV =1%



# GIMA

**Stability:** CV =1%

**Record mode:** L CD display, flash ROM data storage

**Power supply:** Built-in rechargeable lithium battery 3.7 V, 1900 mAh

**Host computer:** DC 5 V, 1 A

**A dapter:** AC 100-240 V, 50/60 Hz

**Dimension - Weight:** 126x73.5xh 30 mm - 0.18 kg