



GIMA

GIMA URINE STRIPS 11 PARAMETERS - professional use

Code: 24074

Category: Urine and blood analyzers and test strips

Unit of sale: box of 100 pcs.

Minimum order: 1

Type: IVD - Professional Use



NSIS: 2420294

CND: W0101060204

EAN13: 4260033718525

Description: 11 parameter strips - tube of 100 strips

For professional use.

These strips can be used by professionals for visual tests

Global check/screening for general health condition.

For use in preliminary screening tests for diabetes, liver diseases, hemolytic diseases, urogenital and kidney disorders and metabolic abnormalities.

Multilingual box: GB, FR, IT, ES, PT, DE, NL, GR, NO, SE, FI, DK, RU, CZ, PL, TR, BG, MT, HU, HR, RO, SK, SL, LV, LT, GA, ET and instructions: GB, FR, IT, ES, PT, DE, NL.

Made in Germany.

Technical Specifications:

Parameters:

- Glucose
- Ascorbic Acid
- Ketones
- Protein
- Ph-Value
- Blood
- Nitrite
- Leucocytes
- Spec. Gravity
- Bilirubin
- Urobilinogen