



GIMA

VET URINE STRIPS 11 parameters

Code: 24100

Category: Vet test strips and devices

Unit of sale: box of 100 pcs.

Minimum order: 1

Type: No medical device



EAN13: 4260033718181

Description: VETERINARY URINE STRIPS - 11 parameters - bottle of 100 strips

For examination of glucose, ascorbic acid, ketones, protein, ph-value, blood, nitrite, leukocytes, specific gravity, bilirubin, urobilinogen.

Each box includes a manual for interpretation of results in each of the following animal: horse, pig, cattle, sheep, goat, dog, cat, bunny and guinea pig.

Multilingual box: GB, FR, IT , ES, PT, NL, DE.

Multilingual instructions: GB, FR, DE, NL.

Ascorbic acid protected.

Performs highest quality standard.

Made in Germany.