



GIMA

ROUND SAFETY NEEDLE 21G - automatic lancets

Code: 24451
Category: Lancets
Unit of sale: box of 100 pcs.
Minimum order: 1
Type: Medical device
Class: II A
NSIS: 2213718
CND: V010401
EAN13: 6970934250274



Description: Easy to use, comfortable, safety lancets to meet various capillary blood requirements.

- **ROUND SAFETY LANCETS**
- **Colour:** yellow
- **Depth of penetration:** 2.2 mm
- **Gauge:** 21G
- **Blood Flow:** 40-60 µl
- **Finger stick:** YES
- **Glucose:** YES
- **Hematocrit:** YES
- **Bleeding:** High

- **single-use:** the build-in-self-destructing design can prevent reuse of the device

- **no accidental injury risk:** needle is hidden before and after use to prevent accidental injuries

- **sterilization guaranteed:** the lancet is Gamma sterilized and the needle tip is protected with a sterility tab until its use allowing a long shelf life

- **no blood contamination**

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, CZ, SE, RO, SK, NL, DK, FI, EE, HR, HU, LT, LV, NO, SI, BG, RU, GR, Arabic.