



GIMA

STREP-A TEST - dual colour - strip

Code: 24519

Category: Professional diagnostic test strips and devices

Unit of sale: box of 50 pcs.

Minimum order: 1

Type: IVD - Professional Use



NSIS: 1702284

CND: W0105090103

EAN13: 8023279245196

Description: STREP-A DUAL COLOUR TEST - strip - box of 50 tests

Dual colour rapid test for Strep A uses the antibody labelled colour particles coated at two separate locations in the test device.

Optical immunoassay technology allows greater sensitivity and less false negative results without affecting specificity.

Standard procedure and easy to view results. If Group A Streptococcus is present in the sample, a visible red Test Line will appear to indicate a positive result and a blue Control Line will also appear to indicate the test is valid.

Kit includes 50 test strips, 50 tubes, 50 swabs, 2 reagents and 2 control solutions.

WARNING: SHORT EXPIRY DATE FOR TESTS. Most tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.

Technical Specifications:

- Relative sensitivity: 95.1% (95%CI: 89.6-98.2%)
- Relative specificity: 97.8% (95%CI: 95.8-99%)
- Overall accuracy: 97.1% (95%CI: 95.3-98.4%)