



# GIMA

## INFLUENZA A/B TEST PANEL

**Code:** 24537

**Category:** Professional diagnostic test strips and devices

**Unit of sale:** box of 30 pcs.

**Minimum order:** 1

**Type:** IVD - Professional Use

**NSIS:** 1704450

**CND:** W0105099004

**EAN13:** 8023279245370



**Description:** INFLUENZA A/B TEST - box of 30 tests (3 packs of 10 pcs.)

**For professional use only.**

**Influenza A/B is a rapid in vitro test to detect type A and B influenza antigens (nucleoprotein) extracted from the nasal specimens.**

**The test is designed to detect nucleoproteins of Flu-A core subtypes, independent of H and N surface proteins (including H1N1, H3N2, H5N1, H5N3, H7N3, H9N2).**

**A sterile nasal swab allows an easy collection of specimen.**

**Results in 10-15 minutes.**

**Includes rapid influenza A/B test panel, swab, test tube, extraction buffer.**

**WARNING: SHORT EXPIRY DATE FOR TESTS.**

**Tests have 24 months expiry at production time. Our stock may vary from 10 to 22 months.**