



GIMA

MULTI DRUG "CUP" TEST for 7 drugs - urine

Code: 24546

Category: Drugs of abuse test strips and devices

Unit of sale: box of 25 pcs.

Minimum order: 1

Type: IVD - Professional Use



NSIS: 1374555

CND: W0102160511

EAN13: 8023279245462

Description: MULTI DRUG "CUP" TEST - urine - box of 25 tests for 7 drugs

For professional use.

Detects following parameters: (MTD/MDMA)
+<COC>+<AMP500>+<MET500>+<THC>+<MOP>.

- easy to use: no urine handling, just insert the key into the active test and read results in 5 min.
- tightly sealed urine cup: no leaking or spilling of urine specimen
- unique design: ensures that urine used in the test, is separated from residual urine, that will remain uncontaminated, to be used for a further test
- key packed separately to ensure that only technician can take the test

Technical Specifications:

- MOP
 - CUT OFF: 300 ng/mL
 - Detection time*: 5-7 days
- COC
 - CUT OFF: 300 ng/mL
 - Detection time*: 2-5 days
- THC
 - CUT OFF: 50 ng/mL
- AMP
 - CUT OFF: 500 ng/mL
 - Detection time*: 2-5 days
- MDMA
 - CUT OFF: 500 ng/mL
 - Detection time*: 5-7 days
- MTD
 - CUT OFF: 300 ng/mL
 - Detection time*: 5-7 days
- MET
 - CUT OFF: 500 ng/mL
 - Detection time*: 3-5 days

* Drug use over detection time not measured