



# GIMA



## MULTI DRUG PANEL TESTS - 8 drugs - urine

**Code:** 24549

**Category:** Drugs of abuse test strips and devices

**Unit of sale:** box of 25 pcs.

**Minimum order:** 1

**Type:** IVD - Professional Use

**NSIS:** 1737723

**CND:** W0102160511

**EAN13:** 8023279245493

**Description:** DRUGS OF ABUSE, URINE TESTS FOR PROFESSIONAL USE

• MULTI DRUG PANEL TESTS - urine - box of 25 tests for 8 drugs

Accurate results in 5 min. on the following drugs: amphetamine, cocaine, marijuana cannabis, methamphetamine, opiates, phencyclidine, MDMA and 6-AM.

Cut-off levels meet SAMHSA guidelines.

Manual: GB, FR, IT, ES, DE, PT, GR, Arabic.

**Technical Specifications:**

**DRUG**

• COC

- Cut off: 150 ng/ml

- Detection time: 2-5 days

• THC

- CUT OFF: 50 ng/mL

- Detection time: 5-8 days

• AMP

- CUT OFF: 500 ng/mL

- Detection time: 2-5 days

• MDMA

- CUT OFF: 500 ng/mL

- Detection time: 5-7 days

• MET

- CUT OFF: 500 ng/mL

- Detection time: 3-5 days

• PCP

- CUT OFF: 25 ng/mL

- Detection time: -

• OPI

- CUT OFF: 2000 ng/mL

- Detection time: 5-7 days

• 6-AM

- CUT OFF: 10 ng/mL

- Detection time: -