



GIMA

MULTI DRUG PANEL TESTS - 6 drugs - urine

Code: 24550

Category: Drugs of abuse test strips and devices

Unit of sale: box of 25 pcs.

Minimum order: 1

Type: IVD - Professional Use

NSIS: 1737755

CND: W0102160511

EAN13: 8023279245509



Description: DRUGS OF ABUSE, URINE TESTS FOR PROFESSIONAL USE

• MULTI DRUG PANEL TESTS - urine - box of 25 tests for 6 drugs

Accurate results in 5 min. on the following drugs: amphetamine, cocaine, marijuana cannabis, methamphetamine, opiates, phencyclidine.

Cut-off levels meet SAMHSA guidelines.

Manual: GB, FR, IT, ES, DE, PT, GR, Arabic.

Technical Specifications:

DRUG

- **COC**
 - Cut off: 150 ng/ml
 - Detection time: 2-5 days
- **THC**
 - CUT OFF: 50 ng/mL
 - Detection time: 5-8 days
- **AMP**
 - CUT OFF: 500 ng/mL
 - Detection time: 2-5 days
- **MET**
 - CUT OFF: 500 ng/mL
 - Detection time: 3-5 days
- **PCP**
 - CUT OFF: 25 ng/mL
 - Detection time: -
- **OPI**
 - CUT OFF: 2000 ng/mL
 - Detection time: 5-7 days