



GIMA

MULTI DRUG PANEL TESTS - 12 drugs - urine

Code: 24551

Category: Drugs of abuse test strips and devices

Unit of sale: box of 25 pcs.

Minimum order: 1

Type: IVD - Professional Use



NSIS: 1737760

CND: W0102160511

EAN13: 8023279245516

Description: DRUGS OF ABUSE, URINE TESTS FOR PROFESSIONAL USE

• MULTI DRUG PANEL TESTS - urine - box of 25 tests for 12 drugs

Accurate results in 5 min. on the following drugs: amphetamine, cocaine, marijuana cannabis, methamphetamine, opiates, phencyclidine, barbiturates, benzodiazepines, ecstasy, morphine, methadone, tricyclic antidepressants.

Cut-off levels meet SAMHSA guidelines.

Manual: GB, FR, IT, ES, DE, PT, GR, Arabic.

Technical Specifications:

DRUG

- MOP
 - CUT OFF: 300 ng/mL
 - Detection time: 5-7 days
- COC
 - CUT OFF: 150 ng/mL
 - Detection time: 2-5 days
- THC
 - CUT OFF: 50 ng/mL
 - Detection time: 5-8 days
- AMP
 - CUT OFF: 500 ng/mL
 - Detection time: 2-5 days
- MDMA
 - CUT OFF: 500 ng/mL
 - Detection time: 5-7 days
- MTD
 - CUT OFF: 300 ng/mL
 - Detection time: 5-7 days
- MET
 - CUT OFF: 500 ng/mL
 - Detection time: 3-5 days
- BZO
 - CUT OFF: 300 ng/mL
 - Detection time: 7-10 days
- BAR



GIMA

- **CUT OFF: 300 ng/mL**
- **Detection time: 2 days - 4 weeks**
- **PCP**
- **CUT OFF: 25 ng/mL**
- **Detection time: -**
- **OPI**
- **CUT OFF: 2000 ng/mL**
- **Detection time: 5-7 days**
- **TCA**
- **CUT OFF: 1000 ng/mL**
- **Detection time: -**