



GIMA

MULTI DRUG SALIVA MIDSTREAM TEST - 6 drugs, 8 parameters - professional

Code: 24566

Category: Drugs of abuse test strips and devices

Unit of sale: box of 20 pcs.

Minimum order: 1

Type: IVD - Professional Use



NSIS: 1702286

CND: W0102160511

EAN13: 8023279245660

Description: MULTI DRUG SALIVA MIDSTREAM TEST - BOX OF 20 TESTS

A multi parameters, all-in-one test for simultaneous detection of 6 different drugs in 10 minutes.

Cost effective system for the screening of workers for illicit substance consumption.

Professional use.

Multilingual box and manual: GB, FR, IT, ES.

Technical Specifications:

Test: Opiates (OPI/MOP)

- Morphine: 40 Cut-off (ng/ml)

- 6 Monoacetylmorphine: 4 Cut-off (ng/ml)

Test: Cocaine (COC)

- Benzoyllecgonine: 30 Cut-off (ng/ml)

Test: Amphetamine (AMP)

- d-Amphetamine: 40 Cut-off (ng/ml)

Test: Oxycodone (OXY)

- Oxycodone: 40 Cut-off (ng/ml)

Test: Methamphetamine (MET)

- d-Methamphetamine: 40 Cut-off (ng/ml)

- 3,4-Methylenedioxymethamphetamine (MDMA): 50 Cut-off (ng/ml)

Test: Marijuana (THC)

- 11-nor-delta-8-TCH-9 COOH: 10 Cut-off (ng/ml)