



GIMA

MULTI DRUG CUP - 7 drugs + S.V.T. for 24570

Code: 24576

Category: Drugs of abuse test strips and devices

Unit of sale: box of 25 pcs.

Minimum order: 1

Type: IVD - Professional Use



NSIS: 1691117

CND: W0102160511

EAN13: 8023279245769

Description: MULTI DRUG CUP - 7 DRUGS + S.V.T.

Simply scan one of the two multi drug tests on drug reader and get a complete analysis on your PC.

- automatic analysis of “visual read” test panels
- automatic test menu recognition (2D barcode)
- in-process system quality control check
- image capture, data management/export functionality, archived results

Main benefits:

- automation: eliminates results interpretation variability/subjectivity
- accurate: performance advantages vs. visual interpretation
- speed: scanning process requires 15-30 sec.
- flexible: semi-portable, simple installation, Windows 7, 8, 10 compatible.

Professional use.

Technical Specifications:

- DRUGS**
- MTD
 - CUT OFF NG/M:L 300
 - Detection time: 5-7 days
 - MDMA
 - CUT OFF NG/M:L 500
 - Detection time: 5-7 days
 - MET
 - CUT OFF NG/M:L 500
 - Detection time: 3-5 days
 - COC
 - CUT OFF NG/M:L 300
 - Detection time: 2-5 days
 - THC
 - CUT OFF NG/M:L 50
 - AMP
 - CUT OFF NG/M:L 500
 - Detection time: 2-5 days
 - MOP
 - CUT OFF NG/M:L 300



GIMA

- Detection time: 5-7 days

S.V.T.

- **OX**
- **S.G.**
- **PH**
- **NiT**
- **Glut**
- **CRE**