



GIMA

MULTI DRUG CUP - 8 drugs + S.V.T. for 24570

Code: 24577

Category: Drugs of abuse test strips and devices

Unit of sale: box of 25 pcs.

Minimum order: 1

Type: IVD - Professional Use



NSIS: 1691125

CND: W0102160511

EAN13: 8023279245776

Description: MULTI DRUG CUP - 8 DRUGS + S.V.T.

Simply scan one of the two multi drug tests on drug reader and get a complete analysis on your PC.

- automatic analysis of "visual read" test panels
- automatic test menu recognition (2D barcode)
- in-process system quality control check
- image capture, data management/export functionality, archived results

Main benefits:

- automation: eliminates results interpretation variability/subjectivity
- accurate: performance advantages vs. visual interpretation
- speed: scanning process requires 15-30 sec.
- flexible: semi-portable, simple installation, Windows 7, 8, 10 compatible.

Professional use.

Technical Specifications:

DRUGS

- MTD
 - CUT OFF: 300 ng/mL
 - Detection time: 5-7 days
- MDMA
 - CUT OFF: 500 ng/mL
 - Detection time: 5-7 days
- MET
 - CUT OFF: 500 ng/mL
 - Detection time: 3-5 days
- COC
 - CUT OFF: 300 ng/mL
 - Detection time: 2-5 days
- THC
 - CUT OFF: 50 ng/mL
- AMP
 - CUT OFF: 500 ng/mL
 - Detection time: 2-5 days
- MOP
 - CUT OFF: 300 ng/mL



GIMA

- Detection time: 5-7 days
- BUP
- CUT OFF NG/ML 10
- Detection time: 5-7 days

S.V.T.

- OX
- S.G.
- PH
- NiT
- Glut
- CRE