



GIMA

PH VAGINAL TEST - strip

Code: 29110

Category: Gyn tests

Unit of sale: box of 10 pcs.

Minimum order: 1

Type: IVD - Professional Use

NSIS: 1691103

CND: W0101060299

EAN13: 8023279291100

Description:

Quick test to determine vaginal pH. An easy routine test to check vaginal pH, results in just 1 minute, to make sure It is within the physiological limits (< 4.5).
An increase in the pH level may indicate a change in the vaginal environment. This could mean there is a possible vaginitis.
The test is ideal for pregnant women and for women with particular diseases such as diabetes.
For professional use.

Standard accessories: The kit contains: 10 cotton swabs + 10 reagent strips in an individual box with user manual (English + Italian)

