



GIMA

CHLAMYDIA RAPID TEST - professional

Code: 29117

Category: Gyn tests

Unit of sale: box of 20 pcs.

Minimum order: 1

Type: IVD - Professional Use



NSIS: 1519227

CND: W0105090101

EAN13: 6970277511032

Description: Provide accurate results in 10 minutes, giving doctors immediate diagnosis. High sensitivity and specificity.

Each box includes 20 tests, reagents, female cervical swabs, tubes, droppers, pipette and user manual.

For professional use.

Technical Specifications:

- Male/Female swab specimen
 - Add reagent A and B to specimen
 - Add 3 drops of extracted specimen
 - Read results at 10 minutes
- Male urine specimen
 - Add reagent B and A to specimen
 - Add 3 drops of extracted specimen
 - Read results at 10 minutes