



GIMA

PAP TEST KIT - sterile

Code: 29702
Category: Kit pap test
Unit of sale: box of 100 pcs.
Minimum order: 1
Type: Medical device
Class: I S
NSIS: 21841
CND: U0899
EAN13: 8023279297027



Description: Pap-Test is a painless examination to detect the presence of any anomalous cells on the cervix. It allows an early diagnosis of uterus neck cancer.
The functional configuration of the kit allows a ready to go operation.

The kit includes: two latex gloves, one perno speculum (medium size), one brush, one Ayre spatula (wooden), two slides and one box for them.

Individually packed in sterile pouch.

Instructions: GB, IT, PT.