



GIMA

PARKER OPERATIVE HALOGEN OTOSCOPE

Code: 31474
Category: GIMA Oscopes
Unit of sale: 1 pc.
Minimum order: 1
Type: Medical device
Class: I
NSIS: 151040
CND: Z12149006
EAN13: 8023279314748



Description: Black model with halogen 2.5 V bulb and bayonet connection system with operative head.

- Reliable and indestructible high quality instrument
- Invaluable for everyday ear nose and throat use
- Easily adjustable Neutral 2.5V vacuum light (made in Germany)
- Quality 3 times magnifying lens
- Ideal for pneumatic tympanic membrane tests, can also be used for rhinoscopic examinations

Supplied with 3 autoclavable ear funnels in a soft bag. Requires 2xC batteries (not included) or rechargeable battery.

Multilingual box and manual: GB, FR, IT, ES, PT, GR, DE, Arabic.