



GIMA

EDAN SE-1 ECG - 1 channel with monitor

Code: 33330
Category: 1-3 channel ECG
Unit of sale: 1 pc.
Minimum order: 1
Type: Medical device
Class: II A
NSIS: 213455
CND: Z12050302
EAN13: 8023279333305



Description: Edan SE-1 Single channel portable ECG.

- foldable LCD screen vivid waveform display
- built-in high resolution thermal printer
- built-in rechargeable Li-ion Battery, AC/DC supply
- complete digital filter, avoid baseline drift, AC, EMG interfere
- automatic baseline adjustment
- auto/manual operation selectable
- lead-off detection and alarm

Technical Specifications:

- Lead: Standard 12 leads
- Acquisition mode: Single
- Sensitivity: 2.5, 5, 10, 20, mm/mV
- Input impedance: > 50 MOhm
- Filter:
 - EMG Filter 35 Hz (-3 dB)
 - AC/DFT Filter 50 Hz / 60 Hz (-20 dB)
- Recording mode: Auto/Manual
- Rhythm lead: Any lead selectable
- Safety standard: IEC I/CF
- Power supply:
 - AC: 100 V - 115 V / 220 V - 240 V, 50/60 Hz
 - Rechargeable Li-ion battery: voltage 14.8 V
 - More than 3 hours continuous operation with battery supply
- Thermal paper: 50 mm x 30 m
- Paper speed: 25 mm/s, 50 mm/s ($\pm 3\%$)
- Dimension: 300 x 260 x 75 mm
- Weight: 2,3 kg with battery included
- Display: 192 x 64 pixels LCD
- Selectable display language: GB, IT, FR, ES, PL, RU, TR, UA
- Interface: RS 232

Standard accessories:

- Patient Cable
- 1 roll of paper



GIMA

1 roll of paper

- 4 limb electrodes
- 1 bottle of gel
- Rechargeable battery
- 6 chest limb electrodes
- Manual: GB, IT (on request: FR, ES, PT, TR)