



GIMA

SPIROBANK II BASIC + SOFTWARE

Code: 33510
Category: Spirometers
Unit of sale: 1 pc.
Minimum order: 1
Type: Medical device
Class: II A
NSIS: 2494320
CND: Z12150102
EAN13: 8023279335101



Description: SPIROBANK II BASIC + MIR SPIRO BASIC PC SOFTWARE

Complete and rapid pocket spirometer for the General Practitioner to carry out spirometry testing anywhere.

No calibration required, use disposable turbine FlowMIR with integrated mouthpiece (33507) or reusable turbine (33526).

- **measures the most important spirometric parameters including: FVC, VC, IVC, IC, ERV, FEV1, FEV1%, PEF, FEF 25-75, FET, EVOL, TEST POST BD**
- **traffic lights for immediate test interpretation**
- **memory up to 10,000 tests**
- **rechargeable long life battery (40 hours continuous operation)**
- **real time tests via USB**
- **selectable display language (GB, IT, ES, FR, DE, PT)**

Supplied with MIR Spiro Basic PC software, USB cable, one disposable turbine and manual GB, IT (FR, ES, PT, DE on www.gimaitaly.com) in a carrying bag.

Made in Italy.

Technical Specifications: Power supply: 3.7 V, 110 mA lithium battery , USB rechargeable
Keypad: 6-key membrane keyboard
Display: LCD monochrome 160x80 pixel
Flow range: +16 l/s
Flow accuracy: $\pm 5\%$ or 200 ml/s whichever is greater
Max volume: 10 l
Volume accuracy: $\pm 2.5\%$ or 50 ml whichever is greater
Data transmission: USB 2.0
Size - weight: 160x55x25 mm - 140 g