


REINFORCED SURGICAL GOWNS 35 g/m2 135 cm - size XL - sterile



Product Code	21016
Unit of sale	box of 30 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	346340
CND	T020402
EAN/UPC	8023279210163

Description

SURGICAL SPECIALISTIC GOWNS - STERILE

Reinforced gowns 35 g
Size XL - length 135 cm

Specialistic, water repellent, breathable gowns for a barrier effect to fluids and bacteria.

- cotton cuffs 8 cm
- ultrasonic welding
- neck velcro closure
- 4 waist straps for internal and external closure

Pre folded
For easy programmed dressing

Welded reinforcement
Water-proof, breathable reinforcement on front and on sleeves ultrasonic welded.

All gowns include two dry-paper towels.
Conforms to UNI EN ISO 13795.

Made in Italy.