

SURGICAL GOWNS 50 g/m2 120 cm - size M - sterile



Product Code	21020
Unit of sale	box of 30 pcs
Minimum order	1
Type	Medical device
Class	I S

RDM (NSIS)	346336
CND	T020401
EAN/UPC	8023279210200

Description

SURGICAL SPECIALISTIC GOWNS - STERILE

Surgical gowns 50 g
Size M - length 120 cm

Specialistic, water repellent, breathable gowns for a barrier effect to fluids and bacteria.

- cotton cuffs 8 cm
- ultrasonic welding
- neck velcro closure
- 4 waist straps for internal and external closure

Pre folded
For easy programmed dressing

All gowns include two dry-paper towels.
Conforms to UNI EN ISO 13795.

Made in Italy.