

SURGICAL GOWNS 50 g/m2 140x160 cm - size XXL - sterile



Product Code	21049
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	No

RDM (NSIS)	2727250
CND	T020401
EAN/UPC	6973866550153
UMDNS	11897

Description

Surgical gowns 50g - sterile

- Size: XXL
- Length: 140 cm
- Width: 160 cm
- Length of sleeve: 61 cm
- Packing size: 66 x 37 x 38 cm

Disposable surgical gowns that guarantee comfort and protection.

- fluid repellent
- abrasion resistant
- ultrasonic sewing
- book folding
- knitted cuff/velcro at neck
- anti microorganism and antibody fluids
- antistatic treatment
- great air permeability
- high strength and reliability
- set-in sleeve

Material: SMMS (spunbond + 2 meltblown + spunbond non wovens).

Weight per piece: 120 ± 2g

Sterile gowns are individually packed.

Sterilized by oxide Ethylene.

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, SE, RO, PL, Arabic.