

KNEE-HIGHS SOCKS - L/XL - medium compression - black



| | |
|---------------|----------------|
| Product Code | 41101 |
| Unit of sale | 1 pair |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 2636023 |
| CND | M030405 |
| EAN/UPC | 8033201400333 |
| PARAF | 925563975 |

Description

Knee-Highs designed for dynamic women who need the preventive action of elastic compression but don't want to compromise on style.

- Medium compression Knee-Highs
- Colour: Black
- Size: Large/X-Large

This elegant product line features a true heel, an anti-shrink cuff and extremely soft fabrics. This line offers two different compression levels (medium and strong) and is available in beige or black.

- Medium compression level: 70 den - 10-14 mm/Hg at the ankle. Ideal for those who would like to initiate a prevention programme, but are mainly interested in the comfort and well-being offered by elastic compression stockings.

Composition: 85% polyamide, 15% elastane.

Multilanguage box and instructions: GB, FR, IT, ES, PT, DE, GR, HU

Made in Italy.