

KNEE-HIGHS SOCKS - L/XL - medium compression - beige



Product Code	41103
Unit of sale	1 pair
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2636023
CND	M030405
EAN/UPC	8033201400302
PARAF	925563999

Description

Knee-Highs designed for dynamic women who need the preventive action of elastic compression but don't want to compromise on style.

- Medium compression Knee-Highs
- Colour: Beige
- Size: Large/X-Large

This elegant product line features a true heel, an anti-shrink cuff and extremely soft fabrics. This line offers two different compression levels (medium and strong) and is available in beige or black.

- Medium compression level: 70 den - 10-14 mm/Hg at the ankle. Ideal for those who would like to initiate a prevention programme, but are mainly interested in the comfort and well-being offered by elastic compression stockings.

Composition: 85% polyamide, 15% elastane.

Multilanguage box and instructions: GB, FR, IT, ES, PT, DE, GR, HU

Made in Italy.