

LUX 100 CLAR TYPE HEAD MIRROR with battery case



Product Code	31109
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2509389
CND	L160299
EAN/UPC	8023279311099

Description

LUX-100 battery head mirror - Ø 100 mm - 6 V
(Battery mounted on head mirror band)

Quality head mirrors Ø 100 mm to perfectly illuminate the examination area in every situation.
Long life 6 V lamp provides bright and natural light, quality silver coated mirror allows optimum light reflection.
Mirror and light beam can be easily and precisely adjusted.

Supplied with multilingual manual: GB, FR, IT, ES, PT, DE, GR, Arabic.

Standard accessories

- Mirror with head band
- Long life bulb
- Head band mounted battery case for 2 x 3V CR-123A (Battery not included)
- Cordura bag