



Product Code	41127
Unit of sale	1 pair
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2636026
CND	M030405
EAN/UPC	8033201400463
PARAF	925564229

Description

Stay ups with gradual compression, silicone lace top, honeycomb mesh

- Medium compression Stay Ups
- Colour: Beige
- Size: X-Large

Stay Ups designed for dynamic women who need the preventive action of elastic compression but don't want to compromise on style.

This elegant product line features a true heel, an anti-shrink cuff and extremely soft fabrics.

This line offers two different compression levels (medium and strong) and is available in beige or black.

- Medium compression level: 70 den - 10-14 mm/Hg at the ankle Ideal for those who would like to initiate a prevention programme, but are mainly interested in the comfort and well-being offered by elastic compression stockings.

Composition: 84% polyamide, 16% elastane.

Multilanguage box and instructions: GB, FR, IT, ES, PT, DE, GR, HU

Made in Italy.