

HEINE DELTAONE DERMATOSCOPE - black - K-210.28.305



Product Code	31140
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2724742
CND	Z12040108
EAN/UPC	4053755196685

Description

HEINE DELTAONE DERMATOSCOPE - black - K-210.28.305

Pocket dermatoscope for non-invasive transient contact and non-contact examination of the skin. It produces extraordinary colour rendering, due to the quality of LEDHQ.
The polarisation allows doctors to work glarefree and without disrupting reflections. It combines diagnostic precision and mobility.

- Features:
- 10x magnification
 - extremely sharp image with the new achromatic HEINE optical system
 - true field of view of Ø 21mm (lens 22mm)
 - 3 levels of light intensity
 - connection to free HEINE DERM App
 - fast change from contact to non-contact examination
 - autoclavable contact plate with scale
 - dioptre compensation
 - toggle function: switch between polarised and non-polarised light at the touch of a button
 - rechargeable Li-ion battery.

Supplied with contact plate with scale, rechargeable Li-ion battery, USB power supply and zipper case.

Technical Specifications

Diopter: -4 to +4 dpt
Illuminance: >17.000 lux in 15 mm distance without contact plate
Colour temp: 5.700 K
Colour rendering index (CRI): = 90
Working distance: 15 mm (non contact mode) or contact to skin
LED life: > 50.000 h
Power supply: Li-ion built-in rechargeable battery, 3.25 W
Input: USB 2.0 Micro B, 5V DC, 650 mA
Charging time: 75 min
Size: 128x57xh 40 mm
Protection: IP20