



|               |                |
|---------------|----------------|
| Product Code  | 31146          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2790695       |
| CND        | Z12040108     |
| EAN/UPC    | 4053755197163 |

Description

HEINE DELTA 30 DERMATOSCOPE - K-230.28.305

Delta 30 dermatoscope stands out due to an excellent colour rendering in LEDHQ, a large field of view and an extremely sharp image with new achromatic optical system.  
Polarisation provides glare-free and nonreflective working conditions. It is comfortable to operate with its angled, ergonomic design and intuitive control.

- 10X magnification
- large field of view with real Ø 30mm (lens 32mm)
- crystal clear image, no matter if you prefer a closer or greater working distance
- 3 levels of light intensity
- fast change from contact to noncontact examination
- no cross-contamination: autoclavable contact plate with scale
- toggle function: switch between polarised and non-polarised illumination made at the touch of a button
- connection to free HEINE DERM App

Supplied with contact plate with scale, rechargeable Li-ion battery, USB power supply and case.

Technical Specifications

Diopter: -4 to +4 dpt  
Illuminance: >25.000 lux in 15 mm distance without contact plate  
Colour temp: 5.700 K  
Colour rendering index (CRI): = 90  
Working distance: 15 mm (non contact mode) or contact to skin  
LED life: > 50.000 h  
Power supply: Li-ion built-in rechargeable battery , 6 W  
Input: USB 2.0 Micro B, 5V DC, 1.2 A or charging station  
Charging time: 120 min  
Size: 195x55xh 70 mm  
Protection: IP20