

# SUCTION CATHETER ch/fr 16 - 50 cm - sterile



|               |                |
|---------------|----------------|
| Product Code  | 22116          |
| Unit of sale  | box of 100 pcs |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II A           |

|        |     |
|--------|-----|
| UK-REP | No  |
| CH-REP | Yes |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2834777       |
| CND        | A06010103     |
| EAN/UPC    | 8023279221169 |
| PARAF      | 977810047     |

## Description

PVC DEHP FREE\* STERILE SUCTION CATHETERS with conical adapter  
\*DEHP is an additive used in the plastic industry for PVC which could be harmful to human reproductive system. Latex is totally DEHP free

Suction catheter ch/fr 16 - orange  
Length: 50 cm

Special technology: inner Ø does not shrink at edge

Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, PL, CZ, SK, HU, EE, LT, LV, DK, SE, FI, RO, BG, GR, HR, SI