

HEINE DELTA 20 T USB+TR DERMATOSCOPE 3.5V Li-Ion immersion + polariz.



Product Code	31161
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	2790708
CND	Z12040108
EAN/UPC	4053755185146
UMDNS	18-021

Description

HEINE DELTA 20® T LEDHQ BETA 4 USB DERMATOSCOPE 3.5V (Li-ion) for polarization and immersion with rechargeable handle for mains socket and USB cable - K-262.28.388

Precision HEINE Optics
- high quality, achromatic optical system provides unsurpassed sharpness and resolution. Focusing system provides from 10X to 16X magnification for detection of even the smallest details and structures.

LED in HQ
- the new standard in LED illumination defining optimal light intensity, homogeneity and colour rendering for each individual examination
- Colour temperature: 5.000 K, Colour Rendering Index (CRI) > 80 on a scale of 100. This ensures true tissue colour for the most accurate diagnosis.

Toggle between polarised and non-polarised illumination modes
- alternation between two modes with the press of a button, no exchange of contact plates necessary. Allows recognition of the “Blink Sign” when observing crystalline structures and presence of milia cysts, increasing diagnostic function.

Optional small contact plate for difficult to access lesions.

Digital documentation capabilities
- SLR photo adapters connect any Delta 20 dermatoscope to major brands of digital cameras.

LED Thermal Management
- ensure long-term performance and product stability. Carefully selected and applied heatdiffusing materials guarantee up to 50.000 hours.

Supplied with hard case, dermatoscopy compendium.