

2 LED HEADBAND ILLUMINATED MAGNIFIER with 5 lens



Product Code	31168
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279311686

Description

2 LED TWO-WAY REGULATION HEADWEARING MAGNIFIER with 5 lenses

Two-way regulation to adjust magnifier.
Provided with 5 lenses 1.0x, 1.5x, 2.0x, 2.5x, 3.5x, which combined in group of two, ensure 11 different bigger magnifications and wider applications. For example, with lenses 1.0X and 2.5X magnification is 3.5X.
Headband is made of PVC leather.
It works with 3xAAA batteries, not included.

Multilingual manual: GB, FR, IT, ES, PT, DE, SE, GR.