

PANTYHOSES - S - strong compression - beige



|               |                |
|---------------|----------------|
| Product Code  | 41170          |
| Unit of sale  | 1 pair         |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2636030       |
| CND        | M030405       |
| EAN/UPC    | 8033201400845 |
| PARAF      | 925565362     |

Description

Pantyhose with gradual compression, reinforced body, honeycomb mesh

- Strong compression Pantyhoses
- Colour: Beige
- Size: Small

Pantyhoses designed for dynamic women who need the preventive action of elastic compression but don't want to compromise on style.

This elegant product line features a true heel, an anti-shrink cuff and extremely soft fabrics.  
This line offers two different compression levels (medium and strong) and is available in beige or black.

- Strong compression level: 140 den - 15-21 mm/Hg at the ankle  
Ideal for those suffering from the initial signs of ankle swelling and discomfort. This level of compression is recommended for initial onset of varicose veins. The product's compression level supports and massages the leg, providing effective relief throughout the day. While this product line boasts an elegant design, particular attention has also been paid to maximum elasticity, which is ensured through the use of a honeycomb fabric (net fabric).

Doctor consult is recommended before use.

Composition: 79% polyamide, 21% elastane.

Multilanguage box and instructions: GB, FR, IT, ES, PT, DE, GR, HU

Made in Italy.