



Product Code	31186
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2807141
CND	Z12040108
EAN/UPC	8961102789779
PARAF	950833780

Description

GIMA POCKET LED DERMATOSCOPE 10X

Modern, compact and lightweight LED pocket dermatoscope with 10X magnification for evaluation of cutaneous lesions.
Precision 10X optical system with achromatic lenses provides distortion-free sharp image.
Graduated skin contact plate to measure pigmented skin lesions offers accurate diagnoses.
Focal range: $\pm 5D$.
Focal length: 25 mm.
It can be powered by two 1.5V AA alkaline batteries or rechargeable battery (both not included).

Supplied in a rigid case with contact plate and multilingual manual: GB, FR, IT, ES, PT, DE, NL, PL, LV, HR, BG, GR.